Neuroform™ Microdelivery Stent System

Instructions For Use

QUESTIONS AND ANSWERS

Q: The Neuroform Microdelivery Stent System seems to be binding with the guidewire, making it difficult to advance the System. What should I do?

A: Use soft guidewires rather than support guidewires because soft guidewires facilitate maneuverability of the Neuroform Microdelivery Stent System and deployment of the Stent. Excess tension can build up in the guidewire resulting in increased friction in the System. Alleviate the friction by slightly retracting the guidewire and 3F Microdelivery Catheter to remove any accumulated tension. If excessive friction continues, confirm that the 3F Microdelivery Catheter saline flush is functioning. With use, guidewires can become kinked and lose their lubricious coatings. If excessive friction persists, consider removing and discarding the guidewire and Neuroform Microdelivery Stent System and replacing it with a new one.

Generally, once the Neuroform Microdelivery Stent System is tracking forward over the guidewire, take advantage of the momentum and continue tracking to a target site that is distal to the aneurysm neck. It is easier to move the Neuroform Microdelivery Stent System from a distal to proximal location across the aneurysm neck instead of trying to reposition it by advancing the Neuroform Microdelivery Stent System.

Q: What is the optimal position of the Stent with respect to the aneurysm?

A: Generally, try to position the Stent so that each end of the Stent is secured in relatively normal areas of the parent vessel. The Stent will be more stable if each end of the Stent is anchored in at least 4-6mm of normal vessel. For example, for aneurysms located in the supraclinoid carotid, it may be better to secure the Stent by deploying the distal end in the M1 (middle cerebral artery, first segment) than trying to deploy it in the few millimeters between the aneurysm and the ICA (internal carotid artery) bifurcation.

When deploying the Stent, care should be taken to use a view that best shows the parent vessel distal to the aneurysm, so that the distal end of the Stent can be accurately deployed with respect to the aneurysm. This view may be different from the view used to advance the Neuroform Microdelivery Stent System, or the view used as a working position for aneurysm embolization.

Deploy the distal end of the Stent as precisely as possible with respect to the aneurysm neck to assure at least 4mm of each end of the Stent lies along the parent vessel. Ensure accurate measurement of the aneurysm neck so that, when the Stent is properly sized per Table 7, the proximal end of the Stent will deploy at the correct location, even if it is difficult to see it because of curves in the vessel.

Q: Which Stent size should I choose if I intend to place the Stent in a vessel that has a different diameter between the proximal and distal ends of the Stent? Example: Vessel increases from 2mm PCA (posterior communicating artery) to a 3.4mm basilar.

A: Choose the Stent sized for the larger vessel. In this example, choose the 3.5mm Stent. This Stent can be deployed safely in the smaller PCA and will be well anchored in the basilar artery.

Q: Is there any problem with deploying the Stent across a branch vessel? Can the Stent be safely deployed across the anterior choroidal artery? What about lenticulostriate arteries or perforators arising from the basilar?

A: No adverse events resulting from branch vessel occlusion or emboli to "jailed" vessels have been observed in the limited clinical study conducted on this Stent (26 patients followed through 6 months). Stents have been placed extending from the M1 (middle cerebral artery, first segment) to the ICA (internal carotid artery) without problems.

Q: A coil was placed only partially into the aneurysm and won't retract back into the embolization microcatheter. Part of the coil is outside the aneurysm. What should I do?

A: Do not attempt to remove the partially deployed coil in this situation. The best solution is to place a second Neuroform Microdelivery Stent System in the parent vessel and parallel to the embolization microcatheter, gently pull back the microcatheter containing the coil to unroll, but not stretch, the coil. Then deploy the second Stent to hold the unrolled coil against the vessel wall, and finally detach the coil. This may be safer than trying to forcibly remove the coil and risk dislodging the coil mass and Stent into the parent vessel. Strict attention to heparinization and antiplatelet medication is important.

Q: A loop, or several loops, of a coil (especially a small diameter coil such as 2mm) are protruding through the interstices of the Stent, and I am unable to reposition it. What should I do?

A: If the risk of leaving part of the coil in the parent vessel is unacceptable, place a second Stent inside the first Stent to pin the herniated coil portion against the wall of the vessel. Three-dimensional angiography using an orthogonal view (i.e., "down the

SMART Therapeutics

SHAPING A MINIMALLY INVASIVE FUTURE

Neuroform™ Microdelivery Stent System

Instructions For Use

barrel") may be helpful to assess whether or not a coil loop is inside the lumen of the Stent, in the parent vessel, or between the wall of the vessel and the Stent. Strict attention to heparinization and antiplatelet medication is important.

Q: I have accidentally started to deploy the Stent, but it is not in the location that I wanted. What should I do?

A: The safest course of action generally is not to try repositioning the Stent, but to continue to deploy the Stent where it is, and then deploy a second Stent at the desired location. Safely deploying a Stent, even in an undesired location will minimize vascular injury. Animal studies have demonstrated that the Stent endothelializes in less than 30 days.

Q: I misjudged the positioning of the Stent and have deployed it with one end adjacent to the aneurysm rather than in the normal part of the parent vessel? What should I do?

A: Leave the guidewire in place, remove the 3F Microdelivery Catheter and 2F Stabilizer Catheter, and insert and deploy a second Stent starting from inside from the first Stent to the normal portion of the parent vessel (telescoping Stents). The second Stent should be of the same diameter or larger than the first.

WARRANTY

SMART Therapeutics warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond SMART Therapeutics control directly affect the device and the results obtained from its use. SMART Therapeutics obligation under this warranty is limited to the repair or replacement of this device and SMART Therapeutics shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. SMART Therapeutics neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. SMART Therapeutics assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such device.

SMART Therapeutics

SHAPING A MINIMALLY INVASIVE FUTURE

Manufactured By:
SMART Therapeutics, Inc.
2551 Merced Street
San Leandro, CA 94577

Customer Service (US and International):
1-510-614-4410

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Humanitarian Use Devices

- [General Information](#)
- [Listing of CDRH Humanitarian Device Exemption Summaries of Safety and Possible Benefit](#)
- [Other Resources](#)

General Information

On June 28, 1996, FDA issued a final rule to carry out provisions of the Safe Medical Devices Act of 1990 regarding humanitarian use devices (HUDs). This regulation became effective on October 24, 1996. An HUD is a device that is intended to benefit patients by treating or diagnosing a disease or condition that affects fewer than 4,000 individuals in the United States per year. A device manufacturer's research and development costs could exceed its market returns for diseases or conditions affecting small patient populations. FDA, therefore, developed and published this regulation to provide an incentive for the development of devices for use in the treatment or diagnosis of diseases affecting these populations.

The regulation provides for the submission of an humanitarian device exemption (HDE) application, which is similar in both form and content to a premarket approval (PMA) application, but is exempt from the effectiveness requirements of a PMA. An HDE application is not required to contain the results of scientifically valid clinical investigations demonstrating that the device is effective for its intended purpose. The application, however, must contain sufficient information for FDA to determine that the device does not pose an unreasonable or significant risk of illness or injury, and that the probable benefit to health outweighs the risk of injury or illness from its use, taking into account the probable risks and benefits of currently available devices or alternative forms of treatment. Additionally, the applicant must demonstrate that no comparable devices are available to treat or diagnose the disease or condition, and that they could not otherwise bring the device to market.

An approved HDE authorizes marketing of the HUD. However, an HUD may only be used after IRB approval has been obtained for the use of the device for the FDA approved indication. The labeling for an HUD must state that the device is an humanitarian use device and that, although the device is authorized by Federal Law, the effectiveness of the device for the specific indication has not been demonstrated.

Listing of CDRH Humanitarian Device Exemption Summaries of Safety and Possible Benefit

HDE Number Approval Date, and Docket Number	Device Name	Company Name and Address	Device Description/Device Indications
H020002 09-SEPT-02 02M-0409	Neuroform Microdelivery Stent System	SMART Therapeutics, Inc. 2551 Merced St. San Leandro, CA 94577 USA	The Neuroform Microdelivery Stent System is intended for use with embolic coils for the treatment of wide neck, intracranial, saccular aneurysms arising from a parent vessel with a diameter of greater than or equal to 2mm and less than or equal to 4.5mm that are not amenable to treatment with surgical clipping. Wide neck aneurysms are defined as having a neck of 4mm or a dome-to-neck ratio of <2.
H010004 09-AUG-02 02M-0361	NEUROLINK® System, including NEUROLINK® Stent & Delivery Catheter and NEUROLINK® Balloon Dilatation Catheter	Guidant Corporation Santa Clara, CA	The NEUROLINK® System is indicated for the treatment of patients with recurrent intracranial stroke attributable to atherosclerotic disease refractory to medical therapy in intracranial vessels ranging from 2.5 to 4.5 mm in diameter with ≥ 50% stenosis and that are accessible to the stent system

H000002 19-DEC-01 02M-0311	VISX Excimer Laser System and Custom Contoured Ablation Pattern (C-CAP) Method™	VISX, Inc.	For the treatment of asymmetrical ablation patterns from previous laser refractive surgery caused by decentration of the treatment as viewed on the Zeiss Humphrey® topography unit and treated with the STAR S3 ActiveTrak™ Excimer Laser System in patients: <ul style="list-style-type: none"> • who exhibit symptomatology supportive of visual defect: reduced best spectacle-corrected visual acuity, debilitating glare, monocular diplopia (double vision), and/or debilitating halos; and, • who pre-operatively have at least a 6 µm difference on the elevation topography, from the lowest point to the highest point, over a 6.5 mm diameter or over the patient's pupil diameter as measured by the Zeiss Humphrey topographer, whichever is larger
H000007 05-APR-02 02M-0167	Amplatzer® PFO Occluder	AGA Medical Corporation Golden Valley, MN	For the non-surgical closure of a patent foramen ovale (PFO) in patients with recurrent cryptogenic stroke due to presumed paradoxical embolism through a PFO and who have failed conventional drug therapy.
H010005 22-MAR-02 02M-0121	Ascension® PIP	Ascension Orthopedics, Inc. Austin, TX	For use in arthroplasty of the proximal interphalangeal (PIP) joint when the patients has soft tissue and bone that can provide adequate stabilization and fixation under high demand loading conditions after reconstruction; and needs a revision of a failed PIP prosthesis, or has pain, limited motion, or joint subluxation/dislocation secondary to damage or destruction of the articular cartilage.
H010002 17-OCT-01 01M-0482	OP-1™ Implant	Stryker Biotech Hopkinton, MA	The device is indicated for use as an alternative to autograph in recalcitrant long bone nonunions where use of autograph is unfeasible and alternative treatments have failed.
H010001 28-AUG-01 01M-0392	Avanta Metacarpophalangeal (MCP) Joint Implant Finger Prosthesis	Avanta Orthopaedics, Inc. San Diego, California	The device is indicated for use in arthroplasty of the MCP joint when either the: <ol style="list-style-type: none"> 1. patient is in need of a revision of failed MCP prosthesis(es); or 2. patient expects to place his/her hands under loading situations which preclude the use of an alternative implant in the painful osteo-arthritic and post traumatic MCP joint.
H000004 19-MAR-01 01M-0485	PROSTALAC Hip Temporary Prosthesis	DePuy Orthopaedics, Inc., a Johnson & Johnson Company	This device is indicated for use as a short-term total hip replacement (THR) in patients who need a two-stage procedure to treat a confirmed infection of their THR and where vancomycin and tobramycin are the most appropriate antibiotics for treatment of the infection based on the susceptibility pattern of the infecting microorganism(s).
H990013 21-FEB-01 01M-0201	Composite Cultured Skin (CCS)	Ortec International, Inc. New York, NY 10032	For use in patients with mitten hand deformities due to Recessive Dystrophic Epidermolysis Bullosa (RDEB) as an adjunct to standard autograft procedures (i.e., skin grafts and flaps) for covering wounds and donor sites created after the surgical release of hand contractures (i.e., mitten hand deformities).
H000001 10-Jan-01 010M-0062	JOMED JOSTENT® Coronary Stent Graft	JOMED AB	For Use in the treatment of free perforations, defined as free contrast extravasation into the pericardium, in native coronary vessels or saphenous vein bypass grafts ≥ 2.75 mm in diameter.
H990012 11-May-00 00M-1354	TAS Ecarin Clotting Time Test	Cardiovascular Diagnostics, Inc.	To be used to determine the anticoagulant effect of recombinant hirudin (r-hirudin) during cardiopulmonary bypass in patients who have heparin induced

			thrombocytopenia (HIT).
H990014 31-March-00 00M-1451	Enterra™ Therapy System (formerly named Gastric Electrical Stimulation (GES) System	Medtronic, Inc.	For the treatment of chronic, intractable (drug refractory) nausea and vomiting secondary to gastroparesis of diabetic or idiopathic etiology.
H990008 09-March-00 00M-1228	Telescopic Plate Spacer (TPS) Spinal System	Interpore Cross International	To replace normal body structures following a vertebrectomy/corpectomy of the spine for metastatic disease in the cervical and/or cervico-thoracic spine (C ₃ -T ₂). The TPS Spinal System implants are intended to correct spinal alignment and stabilize the spinal operative site during fusion. TPS Spinal System implants attach to the spine anteriorly by means of their trapezoidal shape and by screws joined with a plate and spacer component.
H990011 01-Feb-00 00M-0599	CardioSEAL® Septal Occlusion System	Nitinol Medical Technologies, Inc.	For closure of a patent foramen ovale (PFO) in patients with recurrent cryptogenic stroke due to presumed paradoxical embolism through a patent foramen ovale and who have failed conventional drug therapy.
H980006 10-Dec-99 99M-5539	TheraSphere®	MDS Nordion, Inc., Kanata, Ontario, Canada	For radiation treatment or as a neoadjuvant to surgery or transplantation in patients with unresectable hepatocellular carcinoma (HCC) who can have placement of appropriately positioned hepatic arterial catheters.
H990007 07-Dec-99 99M-5327	BioGlue® Surgical Adhesive	CryoLife, Inc.	For use as an adjunct in the surgical repair of acute thoracic aortic dissections.
H980007 30-Sept-99 99M-4810	Shelhigh Pulmonic Valve Conduit Model NR-4000 with "No-React®" Treatment	Shelhigh, Inc.	For replacement of the diseased, damaged, or absent pulmonic artery in small children or infants up to age 4 years, with Transposition of the Great Arteries, Truncus Arteriosus, Tetralogy of Fallot with associated cardiac anomalies or with Pulmonary Atresia, or replacement of failed conduits in young patients with accelerated conduit failure.
H990005 28-Sept-99 99M-4763	CardioSEAL® Septal Occlusion System	Nitinol Medical Technologies, Inc.	For the treatment of patients with complex ventricular septal defects (VSD) of a significant size to warrant closure, but that, based on location, cannot be closed with standard surgical transatrial or transarterial approaches.
H990003 20-Sept-99 99M-4619	Acticon™ Neosphincter	American Medical Systems, Inc.	For the treatment of severe fecal incontinence in post-pubescent males and females who have failed, or are not candidates for, less invasive forms of restorative therapy.
H990004 8-Sept-99 99M-4134	CardioSEAL® Septal Occlusion System	Nitinol Medical Technologies, Inc.	For the treatment of patients with complex single ventricle physiology who have undergone a fenestrated Fontan palliation procedure and require closure of the fenestration.
H980008 19-Feb-99 99M-0255	VOCARE® Bladder System	NeuroControl Corporation Valley View, Ohio 44125	For the treatment of patients who have clinically complete spinal cord lesions (ASIA Classification) with intact parasympathetic innervation of the bladder and are skeletally mature and neurologically stable, to provide urination on demand and to reduce post-void residual volumes of urine. Secondary intended use is to aid in bowel evacuation
H980005 28-Dec-98 99M-0150	VOCARE® Bladder System	NeuroControl Corp. Valley View, OH 44125	For the treatment of patients who have clinically complete spinal cord lesions (ASIA Classification) with intact parasympathetic innervation of the bladder and are skeletally mature and neurologically stable, to provide urination on demand and to reduce post-void residual volumes of urine.
H980002 28-Sep-98 98M-0895	Avanta Proximal Interphalangeal (PIP) Finger Prosthesis	Avanta Orthopaedics, Inc	For use in arthroplasty of the PIP joint when either the: (1) patient is in need of a revision of failed PIP prosthesis (es); or (2) patient expects to place his/her hands under loading situations which preclude the use of an alternative implant in the painful osteo-arthritis and post traumatic arthritis PIP joint.
H970005 30-Apr-98 98M-0452	Perma-Flow® Coronary Graft, Model 2C10	Possis Medical, Inc. Minneapolis, MN 55433	For single or multiple vessel coronary artery bypass in patients who are receiving coronary bypass grafting but who have inadequate autologous conduit to complete the

			required revascularization.
H970004 06-Apr-98 98M-0435	Excorim3 Immunoadsorption System	Cobe BCT, Inc. Lakewood, CO 80215	For use in the treatment of patients with hemophhilia A and B who have Factor VIII or Factor IX inhibitor titers above 10 Bethesda Units/ml (BU/ml). The purpose of the system is to lower the inhibitor levels so that routine clotting factor replacement therapy can be considered. It may be used in an acute setting (to control bleeding during an acute hemmorrhage or for emergency surgery) or as a preventive measure to prepare patients for elective surgery.
H970003 16-Dec-97 98M-0163	Urostim	William E. Kaplan, M.D. and Ingrid Richards, R.N., MSN Chicago, IL 80614	For use in children for the treatment of neurogenic bladder disease secondary to spina bifida
H970001 30-Sep-97 98M-0164	King's College Hospital (KCH) Fetal Bladder Drainage Catheter	Rocket Medical PLC England.	For urinary tract decompression following the diagnosis of post-vesicular obstructive uropathy in fetuses 18 to 32 weeks gestational age
H960001 14-Feb-97 97M-0124	Harrison Fetal Bladder Stent Set (Lowery Modification)	Cook OB/GYN3 Spencer, IN 47460	For fetal urinary tract decompression following the diagnosis of fetal post-vesicular obstructive uropathy in fetuses of 18 to 32 weeks gestational age

Other HDE Information

- Humanitarian Device Exemptions (HDE) Regulation: Questions and Answers; Final Guidance for Industry Text or PDF
- HDE Checklist for Filing Decisions Text or PDF

[Text](#) [PDF](#)

[Text](#) [PDF](#)

[Back to Device Evaluation Page](#)

Updated September 17, 2002

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(2) PALMAZ Balloon-Expandable Stent for Renal Arteries - P890017/S10 [原文]

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New Device Approval

PALMAZ® Balloon-Expandable Stent for Renal Arteries - P890017/S10

This is a brief overview of information related to FDA's approval to market this product. See the links below to the Summary of Safety and Effectiveness and product labeling for more complete information on this product, its indications for use, and the basis for FDA's approval.

Product Name: PALMAZ® Balloon-Expandable Stent for Renal Arteries
Manufacturer: Cordis Corporation
Address: 14201 N.W. 60th Avenue, Miami Lakes, Florida 22014
Approval Date: July 10, 2002
Approval Letter: <http://www.fda.gov/cdrh/pdf/P890017/S10a.pdf>

What is it and when is it used? The PALMAZ® Balloon-Expandable Stent for Renal Arteries (PALMAZ® Stent) is a stainless steel, expandable, wire-mesh tube, mounted over a balloon on the end of a long thin flexible tube called a delivery catheter (Powerflex® Plus delivery catheter). The PALMAZ® Stent is used in patients who have a narrowing in an artery caused by the deposition of a fatty material (atherosclerotic plaque) in a section of the renal artery (blood vessel leading to the kidney), which was not successfully opened in an earlier procedure called balloon angioplasty. This narrowing could restrict blood flow to the kidneys; complete blockage of blood flow can lead to the loss of kidney function.

How does it work? The PALMAZ® Stent, mounted on the delivery catheter, is inserted into a blood vessel in the groin (the femoral artery) and is advanced to the narrowed section of the renal artery. When the PALMAZ Stent is correctly positioned within the renal artery, the balloon is inflated, causing the stent to expand. Expansion of the stent pushes the plaque aside, opening the narrowed section of the artery, restoring normal blood flow to the kidney. The balloon on the delivery catheter is deflated and the delivery catheter is removed from the patient. The PALMAZ® Stent remains permanently implanted within the renal artery, acting as a scaffold to support the newly opened section.

What will it accomplish? Expansion of the PALMAZ® Stent within the narrowed section of the renal artery opens the narrowing, allowing more blood flow to the kidney. If the narrowing is not treated, a person can develop high blood pressure (the kidney's response to decreased blood flow) and if the narrowing is severe, loss of kidney function may occur over time.

When should it not be used? The PALMAZ® Stent should not be used in patients with bleeding disorders or in patients who have transplanted kidneys or kidneys that have had bypass surgery.

Additional information: Summary of Safety and Effectiveness and labeling are available at:
<http://www.fda.gov/cdrh/pdf/P890017/S10.html>

Other:

National Institutes of Health, National Library of Medicine:

- Shows placement of a stent in a coronary artery—similar procedure to placing a stent in a renal artery :
<http://www.nlm.nih.gov/medlineplus/ency/presentations/100160.6.htm>
- Renal artery stenosis with drawings of kidney: <http://www.nlm.nih.gov/medlineplus/ency/article/001273.htm#visualFile>
- Acute renal occlusion: <http://www.nlm.nih.gov/medlineplus/ency/article/001272.htm#visualFile>
- Description of X-ray with contrast: <http://www.nlm.nih.gov/medlineplus/ency/article/003781.htm>

Updated 12/26/2002