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
FemExam[®]

pH and Amines TestCard™

For the assessment of elevated vaginal
pH and the presence of amines.

Instructions For Use

PRODUCT NO.: 52-7464-5
52-7464-25
52-7464-50

DISPOSABLE: Discard after single use.

CAUTION: For professional use only.

PACKAGED: 5; 25; and 50 individual devices per box.

STORE AT ROOM TEMPERATURE:

59°-77°F (15°-25°C).

USE WITH VAGINAL FLUID ONLY.

FOR IN VITRO DIAGNOSTIC USE

INDICATIONS FOR USE:

The FemExam[®] Testcard™ contains two qualitative, colorimetric tests for use in the characterization of a vaginal fluid sample:

1. A pH test that differentiates vaginal fluid pH less than pH 4.7 from vaginal fluid pH equal to or greater than pH 4.7.
2. A test that detects alkali volatilized amines in vaginal fluid.

For *In Vitro* Diagnostic Use Only. For Professional Use Only.

SUMMARY AND EXPLANATION OF THE TEST:

Vaginal discharge and odor are among the most frequent gynecological complaints and two of the most common reasons for an adult woman to visit a physician (1-3). Because vulvovaginal infections can result in serious clinical sequelae (4-9), vulvovaginal symptoms and signs warrant careful evaluation and appropriate therapy. The diagnostic approach to a patient with vaginal discharge should include a detailed sexual history and ascertainment of the source of the discharge (i.e., cervical or vaginal) (1,10,11). The clinical evaluation of vaginal discharge should include examination of the following parameters: appearance (color, amount and consistency); pH; release of a fishy amine odor upon the addition of 10% KOH (the whiff test); the presence of clue cells; and wet mount microscopy (1,10-12).

Elevated vaginal fluid pH (pH 4.7 or higher) and the detection of a fishy odor upon adding 10% KOH to vaginal fluid represent two of the four Amsel criteria (13). The other two Amsel criteria include examination of the vaginal discharge for color, amount and consistency, and microscopic evaluation of vaginal fluid for clue cells.

The FemExam[®] TestCard™ is designed to provide indications of elevated vaginal fluid pH (pH $\geq$ pH 4.7) and the presence of volatile vaginal fluid amines. It is intended to be used as part of the clinical evaluation of vaginal discharge. The formation of a blue plus sign (+) within each test circle indicates a positive test result; a blue minus sign within each test circle indicates a negative test result.

CooperSurgical
Shelton, CT 06484
Phone: (203) 929-6321
Toll Free: (800) 243-2974
Fax: (203) 926-9435

PRINCIPLES OF THE FEMEXAM® TESTCARD™:

Principle of the pH Test:

The FemExam® TestCard pH test system contains a traditional colorimetric pH indicator, nitrazine yellow, which produces a visual, geometric color change within two minutes of specimen application. When contacted with a vaginal fluid specimen at or above pH 4.7, the pH test produces a blue plus sign against a greenish-yellow background. When contacted with a vaginal fluid specimen below pH 4.7, the pH test produces a blue minus sign against a greenish-yellow background.

A narrow minus sign across the center of the pH test serves as the positive procedural control and is produced when the pH test is wetted with a vaginal fluid specimen. The positive control contains nitrazine yellow specifically formulated to change color (from greenish-yellow to blue) regardless of whether the pH of the vaginal fluid specimen is normal or elevated. Although the area surrounding the undeveloped plus sign within the pH test circle also contains nitrazine yellow, it is specifically designed not to change color, regardless of whether the pH of the vaginal fluid specimen is normal or elevated, and serves as the negative background control.

Principle of the Amine Test:

The FemExam® TestCard™ amine test system contains a traditional colorimetric pH indicator, bromocresol green (BCG), which produces a visual, geometric color change within two minutes of specimen application. When contacted with a vaginal fluid specimen that contains alkali volatilizable amines, the amine test produces a blue plus sign against a yellow background. When contacted with a vaginal fluid specimen that does not contain alkali volatilizable amines, the amine test produces a blue minus sign against a yellow background.

The TestCard™ amine test employs a very thin dried film of BCG in the yellow test circle which is surrounded by a thick black ring covered with a dried alkali, sodium aluminate. The BCG on the yellow amine test circle produces a visual color change (blue plus sign) in response to the pH elevation produced when the BCG is contacted with a vaginal fluid specimen containing amines volatilized by contact with the sodium aluminate on the thick black ring. The amine test detects volatile amines at concentrations above 0.5mM.

Sodium aluminate is a strong alkalization reagent analogous to the potassium hydroxide (KOH) used to perform the whiff test capable of volatilizing vaginal fluid amines. The sodium aluminate on the thick black ring is physically separated from the BCG in the yellow circle. Additionally, the BCG is specifically formulated to prevent the pH indicator from turning color when contacted by an alkalized vaginal fluid sample that does not contain volatile amines. However, the BCG is permeable to alkaline amine vapor released by contact of a vaginal fluid sample containing alkali volatilizable amines with the sodium aluminate on the thick black ring and produces a blue color change. A narrow blue minus sign across the center of the amine test serves as the positive procedural control and is produced when the amine test is wetted with a vaginal fluid specimen. The positive control contains BCG specifically formulated to change color (from yellow to blue) regardless of whether the vaginal fluid specimen does or does not contain alkali volatilizable amines. Although the area surrounding the undeveloped plus sign within the test circle also contains BCG, it is specifically designed not to change color upon the addition of the vaginal fluid specimen, regardless of whether it contains alkali volatilized amines, and serves as the negative background control.

REAGENTS:

- pH Test - Nitrazine Yellow (1.3 µg / test)
- Amine Test - Bromocresol Green (1.8 µg / test)
- Amine Test Alkali Ring - Sodium Aluminate (equivalent to 29.2 µg HCL / test)

WARNINGS AND PRECAUTIONS:

- For *In Vitro* Diagnostic Use Only.
- Do not use a FemExam® TestCard™ after the expiration date printed on the foil pouch or if the foil pouch is not sealed.
- Do not remove the FemExam® TestCard™ from the foil pouch until ready to use. Once the foil pouch has been opened, the TestCard™ must be used within 30 minutes.
- Cards with a blue plus sign (+) or blue minus sign (–) in either test circle prior to the application of a clinical specimen should not be used and should be discarded.
- This product is intended only for vaginal fluid use.
- Do not touch or collect vaginal fluid near the cervix.
- Do not use vaginal fluid specimens that contain blood. The FemExam® TestCard™ should not be used in women who are menstruating or otherwise have vaginal blood.

WARNINGS AND PRECAUTIONS CON'T:

- The thick black ring surrounding the yellow amine test circle contains an alkali, sodium aluminate, and should be handled with care. Do not expose unprotected skin to the alkaline sodium aluminate since it may cause irritation.
- Used FemExam® TestCard™, sample swabs and other associated materials should be considered potentially infectious and discarded in appropriate receptacles. Products should never be reused.
- Patient vaginal swabs are not appropriate for any other purpose, including bacterial culture, after performing the FemExam® TestCard™

STORAGE AND STABILITY:

Foil pouches should be stored at room temperature 59° - 77°F (15° - 25°C) to ensure that the FemExam® TestCard™ remains stable until the expiration date printed on the foil pouch.

INDICATIONS OF INSTABILITY:

Indications of possible product deterioration are a blue plus sign (+) or blue minus sign (-) in either test circle prior to the application of a clinical specimen.

SPECIMEN COLLECTION:

Collect an undiluted vaginal fluid sample with a cotton swab. The swab can be moistened with vaginal fluid by gently stroking the walls of the vagina. Moisten the swab thoroughly with the sample. The sample should be tested immediately.

NOTE: Do not use vaginal fluid specimens that contain blood. Do not touch or collect vaginal fluid near the cervix.

MATERIALS:

Materials Provided:

- Individual foil pouched FemExam® pH and Amines TestCard™
- Instructions for Use

Material Not Provided:

- Timer
- Cotton Swabs

DIRECTIONS FOR USE:

CAUTION: APPLY THE SWAB FIRST TO THE pH TEST AND THEN TO THE AMINE TEST. The alkali reagent in the thick black ring surrounding the yellow amine test circle may interfere with pH test results if the amine test is done first.

DIRECTIONS FOR USE CON'T:

The FemExam® TestCard™ contains two specimen application areas one for the **pH test** (Test Number 1) and one for the **amine test** (Test Number 2).

1. Remove a TestCard™ from the protective foil pouch prior to use.

The TestCard™ should be used within 30 minutes once the foil pouch has been opened.

2. Perform the pH test first.

Using a circular motion, gently rub the moistened swab over the entire surface of the pH test several times until the entire surface is completely wetted. The same swab may then be used to perform the amine test.

3. Apply the swab to the amine test.

Using a circular motion, gently rub the moistened swab over the entire surface of the amine test starting from the outside black ring and proceeding towards the center of the yellow amine test circle. The alkali is picked up on the swab which is then rubbed over the surface of the yellow amine test circle. The entire surface of the amine test should be completely wetted with the vaginal fluid specimen.

CAUTION: The black ring surrounding the yellow amine test circle contains a dry alkali, sodium aluminate, which may cause skin irritation. Do not contact with exposed skin.

4. The results of the pH test and the amine test are interpreted within 2 minutes after addition of the vaginal fluid specimen to both tests.

NOTE: Volatilized amines are transient and the amine test is correctly interpreted as positive even if the blue plus sign forms and subsequently fades within 2 minutes.

RESULTS:

pH Test:

The formation of a blue PLUS sign (+) within the pH test circle indicates a positive pH test ($\text{pH} \geq 4.7$); a blue MINUS sign (-) within the pH test circle indicates a negative pH test ($\text{pH} < 4.7$). A positive pH test is indicative of elevated vaginal fluid pH ($\text{pH} \geq \text{pH } 4.7$). A negative pH test is indicative of normal vaginal fluid pH ($\text{pH} < \text{pH } 4.7$). Elevated vaginal fluid pH is one of the four Amsel criteria.

RESULTS CON'T:

Amine Test:

The formation of a blue PLUS sign (+) within the amine test circle indicates a positive amine test (presence of amines volatilized by an alkali); a blue MINUS sign (-) within the amine test circle indicates a negative amine test (absence of amines volatilized by an alkali). A positive amine test is indicative of the presence of alkali volatilized amines. The amine test detects volatile amine concentrations above 0.5mM. A negative amine test is indicative of the absence of alkali volatilized amines. The presence of volatile vaginal fluid amines is one of the four Amsel criteria.

Test Result Combinations:

There are four possible FemExam® TestCard™ results: (1) positive pH test and negative amine test; (2) positive pH test and positive amine test; (3) negative pH test and positive amine test; and (4) negative pH test and negative amine test.

QUALITY CONTROL:

Each test on the Fem Exam® TestCard™ contains both a built-in positive procedural control and a negative background control. The controls are performed simultaneously with the tests by the addition of an undiluted vaginal fluid specimen to each test circle.

Positive Procedural Control:

The formation of a blue minus sign (-) in each test circle within 2 minutes is the positive procedural control. If no minus sign forms within 2 minutes, this indicates a possible failure of that test due to either an inadequate vaginal fluid specimen or test reagent deterioration. Results for that test should be disregarded and/or that test repeated.

Negative Background Control:

The background of each test should not change color although there may be some minor color variations or darkening of the background color as a result of wetting the test with the vaginal fluid specimen. A substantial change in background color of a test indicates a potential failure of that test due to the presence of an interferent. Results for that test should be disregarded and/or that test repeated.

External Controls:

Solutions containing the following materials may also be used as separate controls. If the controls do not perform as expected repeat with another FemExam® TestCard™ or call Customer Service at CooperSurgical (1-800-234-2974). Follow state, local and federal regulations for quality control.

pH Test Negative Control:

Citric Acid, Monohydrate	1.68 grams
Sodium Citrate, Dihydrate	1.18 grams
Bovine Serum Albumin	8.00 grams
Deionized Water	Dissolve to final volume of 40ml.

pH Test Positive Control:

Citric Acid, Monohydrate	0.12 grams
Sodium Citrate, Dihydrate	0.42 grams
Bovine Serum Albumin	8.00 grams
Deionized Water	Dissolve to final volume of 40ml.

Amine Test Negative Control:

N-Propanol	0.10 grams
Methylcellulose	2.50 grams
Citric Acid, Monohydrate	3.93 grams
Sodium Citrate, Dihydrate	3.07 grams
Deionized Water	Dissolve to final volume of 100ml.

Amine Test Positive Control:

N-Propanol	4.00 grams
Methylcellulose	2.50 grams
Citric Acid, Monohydrate	0.09 grams
Sodium Citrate, Dihydrate	0.32 grams
Trimethylamine, HCL	0.05 grams
Deionized Water	Dissolve to final volume of 100ml.

LIMITATIONS OF THE PROCEDURE:

- Test results may be affected by improper specimen collection and/or handling.
- Elevated vaginal fluid pH and/or the presence of volatile vaginal fluid amines is not diagnostic for a vaginal disorder.
- The FemExam® TestCard™ should not be used in women who are menstruating or otherwise have vaginal blood.
- Failure to follow instructions may adversely affect performance of the FemExam® TestCard™ as described in the package insert.

EXPECTED VALUES:

The FemExam® TestCard™ pH test can differentiate vaginal fluid pH < 4.7 from vaginal fluid pH ≥ 4.7. A positive FemExam® TestCard™ pH test means that the pH of the vaginal fluid is ≥ 4.7. A negative pH test result means that the pH of the vaginal fluid is < pH 4.7.

The FemExam® TestCard™ amine test can detect volatile vaginal fluid amines at a concentration above 0.50mM. A positive FemExam® TestCard™ amine test means that the vaginal fluid contains volatile vaginal amines at a concentration above 0.50mM. A negative amine test result means that the vaginal fluid contains volatile amines at a concentration below 0.50mM.

Paired Subpopulations:

A total of 604 women were tested, among these women, 416 had vaginal complaints; 40 were pregnant; 98 had candidiasis; 43 had trichomoniasis; and 94 had engaged in sexual intercourse or used vaginal products within 24 hours of testing (some women had more than one of these conditions). There were no significant differences in performance of the FemExam® TestCard™ among the subpopulations of women with these conditions.

PERFORMANCE CHARACTERISTICS:

Method Comparison:

The FemExam® TestCard™ pH and amine tests were evaluated in a statistically designed study conducted by 12 site personnel at 5 geographically separate clinical sites. While a total of 629 women were enrolled and consecutively tested, individual tables and analysis utilize fewer than 629 women. Incomplete data was available on individual women, thereby reducing the number of specimens used in a particular analysis. The only protocol exclusion criteria was the presence of blood in the vagina. Three different production lots of the FemExam® TestCard™ were used in the clinical study.

FemExam® TestCard™ pH Test:

In the clinical study, independent investigators evaluated the performance of the FemExam® TestCard™ pH test to differentiate vaginal fluid pH less than pH 4.7 from vaginal fluid pH equal to or greater than pH 4.7. The performance of the FemExam® TestCard™ pH test was established by comparison to a commercial pH test strip (ColorpHast® pH test strip) in 607 women (TABLE 1). Any result equal to or greater than pH 4.7 on the pH test strip was considered a positive test; any result less than pH 4.7 was considered a negative test.

TABLE 1:

FEMEXAM® TESTCARD™ pH TEST COMPARED TO A pH TEST STRIP

COMPARATIVE METHOD	RESULT	pH TEST STRIP		N
		POSITIVE	NEGATIVE	
FEMEXAM® TESTCARD™ pH TEST	POSITIVE	263	40	607
	NEGATIVE	46	258	

Positive agreement between the FemExam® TestCard™ pH test and the pH test strip was 85.1%; negative agreement was 86.6%; and overall agreement was 85.8%.

FemExam® TestCard™ Amine Test:

In the clinical study, independent investigators also evaluated the performance of the FemExam® TestCard™ amine test to detect alkali volatilizable amines in vaginal fluid. The performance of the FemExam® TestCard™ amine test was established by comparison to the whiff test in 617 women (TABLE 2). The FemExam® TestCard™ amine test detects volatile vaginal fluid amines at concentrations above 0.50mM.

TABLE 2:

FEMEXAM® TESTCARD™ AMINE TEST COMPARED TO THE WHIFF TEST

COMPARATIVE METHOD	RESULT	REFERENCE METHOD: WHIFF TEST		N
		POSITIVE	NEGATIVE	
FEMEXAM® TESTCARD™ AMINE TEST	POSITIVE	166	36	617
	NEGATIVE	30	385	

Positive agreement between the FemExam® TestCard™ amine test and the whiff test was 84.7%; negative agreement was 91.4%; and overall agreement was 89.3%.

CUTOFF PRECISION:

FEMEXAM® TESTCARD™ pH TEST:

Analytical Samples:

The FemExam® TestCard™ pH test produces a positive result when the pH of the vaginal fluid is ≥ pH 4.7. This performance was defined in laboratory studies employing blinded, coded buffer solutions at eight different pH values (pH 4.0, 4.3, 4.4, 4.5, 4.6, 4.7, 4.8, and 5.0) (TABLE 3). A total of 104 observations were made.

TABLE 3:
FEMEXAM® TESTCARD™ pH TEST
CUTOFF PRECISION

COMPARATIVE METHOD	RESULT	CODED BUFFER SOLUTIONS		N
		POSITIVE (pH ≥ 4.7)	NEGATIVE (pH < 4.7)	
FEMEXAM® TESTCARD™ pH TEST	POSITIVE	35	0	104
	NEGATIVE	4	65	

100% of the observations correctly interpreted the FemExam® TestCard™ pH test as negative using coded buffer solutions with a pH < 4.7. 90% of the observations correctly interpreted the FemExam® TestCard™ pH test as positive using coded buffer solutions with a pH ≥ 4.7. When compared to the coded buffer solutions, the FemExam® TestCard™ pH test produced an overall agreement of 96%.

Clinical Samples:

The performance of the FemExam® TestCard™ pH test was also verified in the clinical study using vaginal fluid specimens from 607 women (TABLE 1). Clinical study participants determined the pH of the vaginal fluid specimens using the ColorpHast® pH test strip and recorded the numerical pH value on the Case Report Form. These pH test strips have 14 individual pH color keys covering the pH range of 4.0-7.0 which are defined by a color chart provided by the manufacturer. Applying the Mann-Whitney Rank Sum Test, no statistically significant differences (positive, negative or overall) were found between the results of the pH test strip and the FemExam® TestCard™ pH test.

FEMEXAM® TESTCARD™ AMINE TEST:

Analytical Samples:

The FemExam® TestCard™ amine test produces a positive result for volatile vaginal fluid amines at concentrations above 0.50mM. This sensitivity was defined in laboratory studies employing serially diluted calibrated test solutions of individual amines and control solutions containing no amines (TABLE 4). Amines were tested at 11 concentrations: 0, 0.05mM, 0.10mM, 0.25mM, 0.5mM, 1.0mM, 2.0mM, 3.0mM, 4.0mM, 5.0mM and 10.0mM. Coded solutions were tested in triplicate by investigators blinded to the contents of the test solutions. A total of 165 observations were made.

TABLE 4:
FEMEXAM® TESTCARD™ AMINE TEST
CUTOFF PRECISION

COMPARATIVE METHOD	RESULT	STANDARDIZED AMINE SOLUTION ^(a)		N
		POSITIVE (≥ 0.50mM)	NEGATIVE (< 0.50mM)	
FEMEXAM® TESTCARD™ AMINE TEST	POSITIVE	105	0	165
	NEGATIVE	0	60	

(a) Solutions of five amines (trimethylamine, methylamine, isobutylamine, amylamine and phenylethylamine) were tested.

The minimum concentration of each volatile amine which produced a positive FemExam® TestCard™ amine test result in 100% of the tests was 0.50mM. All concentrations of each amine at a concentration of 0.50mM or higher produced a positive FemExam® TestCard™ amine test result. When compared to the standardized amine solutions, the FemExam® TestCard™ amine test produced an overall agreement of 100%.

Clinical Samples - Quantitative Amine Analysis:

The performance of the FemExam® TestCard™ amine test was verified by quantitating volatile amines in vaginal fluid specimens. Vaginal fluid specimens from 123 of the 629 women who participated in the clinical study were quantitatively analyzed for volatile amines. Of these 123 vaginal fluid specimens, 109 contained amines at concentrations above 0.50mM. The results of the FemExam® TestCard™ amine test were compared to the results of the whiff test in these 109 vaginal fluid specimens with quantitatively determined volatile amine concentrations above 0.50mM to verify the cutoff precision of the FemExam® TestCard™ amine test in clinical samples.

**TABLE 5:
FEMEXAM® TESTCARD™ AMINE TEST CUTOFF
PRECISION COMPARED TO THE WHIFF TEST IN
VAGINAL FLUID SPECIMENS QUANTITATED FOR
AMINE CONTENT**

COMPARATIVE METHOD	RESULT	REFERENCE METHOD: WHIFF TEST		N
		POSITIVE	NEGATIVE	
FEMEXAM® TESTCARD™ AMINE TEST	POSITIVE	89	0	109
	NEGATIVE	2	18	

Positive agreement between the FemExam® TestCard™ amine test and the whiff test was 98%; negative agreement was 100%; and overall agreement was 98% in 109 vaginal fluid specimens with volatile amine concentrations above 0.50mM.

Clinical Samples:

The performance of the FemExam® TestCard™ amine test was also verified in the clinical study using vaginal fluid specimens from 617 women (TABLE 2). Clinical study participants determined the presence of detectable volatile vaginal fluid amines using the whiff test and recorded the results on the Case Report Form. Applying the Mann-Whitney Rank Sum Test, no statistically significant differences (positive, negative or overall) were found between the results of the whiff test and the FemExam® TestCard™ amine test.

REPRODUCIBILITY STUDY:

Analytical Samples:

The total reproducibility of the FemExam® TestCard™ was evaluated in three separate studies performed at three different sites by different types of users: (1) at 3 physician offices by 4 first time users; (2) at 2 clinical sites by 2 experienced, trained health care providers; and (3) in-house by 3 trained users. The study was conducted following the general guidelines outlined in the National Committee for Laboratory Standards (NCCLS) EP5-T2 (Vol. 12 No. 4). The samples tested contained different levels of volatile amines and pH buffers corresponding to a negative test result and a positive test result.

**TABLE 6:
REPRODUCIBILITY OF THE FEMEXAM®
TESTCARD™ pH AND AMINE TESTS USING
ANALYTICAL SAMPLES**

EXPECTED RESULTS	pH	AMINE CONC.	OBSERVATIONS	RESULTS	
				CORRECT	%
POSITIVE	4.7		111	111	100
POSITIVE	5.5		111	111	100
NEGATIVE	4.0		111	111	100
NEGATIVE		0.0	111	109	98.2
POSITIVE		1.0	111	110	99.1
POSITIVE		5.0	111	111	100

The results of this reproducibility study confirms the first time users, trained clinicians and trained in-house users are able to perform and interpret FemExam® TestCard™ pH and amine tests accurately.

Clinical Samples:

The reproducibility of the FemExam® TestCard™ was also tested using duplicate vaginal fluid specimens obtained from 94 women at two separate clinical sites. Each of the two swabs was used to perform a FemExam® TestCard™ pH and amine test. The results from each FemExam® TestCard™ were compared (TABLE 7).

**TABLE 7:
REPRODUCIBILITY OF THE FEMEXAM®
TESTCARD™ pH AND AMINE TESTS USING
VAGINAL FLUID SAMPLES**

FEMEXAM® TESTCARD™	OBSERVATIONS	RESULTS	
		AGREEMENT	%
pH TEST	196	188	96
AMINE TEST	198	186	94

The two FemExam® TestCard™ pH tests produced an overall agreement of 96% using duplicate vaginal fluid specimens. The two FemExam® TestCard™ amine tests produced an overall agreement of 94% using duplicate vaginal fluid specimens.

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