LAP Mentor™

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Gaumard®
Simulators for Health Care Education


Emergency care training starts with **HAL®**

- Easy to use
- Fully-functional in transit
- Wireless control and documentation

HAL® is mobile, wireless, and fully responsive...even in transit


Airway

- Oral or nasal intubation
- Programmable airway to control tongue edema, pharyngeal swelling, and laryngospasm
- Perform tracheostomy or needle cricothyrotomy
- Sensors detect depth of intubation
- Automatic unilateral chest rise with right mainstem intubation
- Multiple upper airway sounds synchronized with breathing


Breathing

- Control rate and depth of respiration and observe chest rise
- Ventilation is measured and logged
- Gastric distention with excess bag-valve-mask ventilation
- Select independent left, right, upper, and lower lung sounds
- Multiple lung sounds are synchronized with selectable breathing patterns
- Accommodates assisted ventilation, including bag-valve-mask and mechanical support
- Four needle decompression sites
- Left and right unilateral chest rise simulate pneumothoraces



Circulation



- Multiple heart sounds
- Chest compressions are measured and logged
- Blood pressure can be taken bilaterally using a cuff, palpation, or auscultation
- Radial pulse stops when cuff pressure exceeds systolic pressure
- Korotkoff sounds audible between systolic and diastolic pressures


HAL's electrically conductive skin regions allow the use of *real* equipment for obtaining an ECG, performance of temporary pacing, and defibrillation

Attach real electrodes and view HAL's ECG generated in real time. HAL's ECG features physiologic variations in rhythm, never repeating textbook patterns. Here pacing therapy converts HAL's profound bradycardia into paced ventricular rhythm.


Attach AED pads directly to HAL's conductive skin. Your AED will display HAL's ECG, analyze his cardiac rhythm, and advise action. Program HAL's response to defibrillation. Stack shocks as needed. AED shown converting HAL's ventricular fibrillation.


Use monophasic or biphasic defibrillators on HAL's skin. Shock HAL using your defibrillator just like a real patient. HAL even distinguishes between defibrillation and synchronized cardioversion. Shock resolves pulseless ventricular tachycardia.

Control HAL® using a wireless tablet PC to quickly generate multiple life threatening situations ... and track the actions of care providers

- Use our preprogrammed scenarios. Modify them or create your own
- Specify HAL's condition on the Details page
- Save each condition as a physiologic state on the Palette page

- Combine physiologic states to build a scenario
- Select the time HAL is in each state
- Build in delays between states or specify smooth transitions
- When activated HAL will progress from one state to the next. Use stylus control to jump between states
- Actions are logged and time stamped for evaluation and debriefing

HAL® Mobile Emergency Care Team Trainer

Rush HAL from the accident scene, while care providers diagnose and treat his condition using real monitoring and resuscitation equipment

Control HAL at distances up to 300 meters and between rooms and floors of conventional buildings

HAL smoothly transitions between physiologic states in response to commands from a wireless PC

HAL is fully functional where you need him:

- Back-board or EMS cart
- EMS vehicles
- Life flights
- Simulated disaster scenes
- Your simulation center

Use our scenarios designed by emergency medical professionals or quickly create your own

HAL has conductive skin regions so you can:

- Apply real electrodes and AED pads
- Use real EMS equipment
- See HAL's ECG on your AED

HAL is perfect for competency based programs:

- Sensors track student actions
- Changes to HAL's condition and care provided are time stamped and logged for evaluation and debriefing
- Instructors evaluate caregiver interventions with a single click and insert notes on a real time performance log
- Save caregiver performance, forward to any standard printer, or send anywhere via the internet

Palpation


Bilateral carotid, radial, brachial, femoral, and pedal pulses operate continuously. Pulse strengths vary with HAL's blood pressure and pulses are synchronized with the ECG

Venous Access

Bilateral IV training arms; adjust HAL's response to boluses and/or IV infusion

Trauma

Four needle decompression sites. Left and right unilateral chest rises simulate pneumothoraces

HAL® Mobile Emergency Care Team Trainer

Simulator

- Full-size adult
- Self-contained respiratory and circulatory functions
- Two-way wireless communications

Control

- Wireless tablet PC with stylus control
- Preprogrammed scenarios
- Build your favorite scenarios

Accessories



- Two (2) interchangeable batteries
- 100-240 VAC charger
- Blood pressure cuff
- Six (6) neck collars
- Six (6) replaceable pneumothorax decompression sites
- Instructions for use
- CD-ROM tutorial
- Carrying case with rollers

HAL S3000

FOB Factory, Miami, Florida, USA

Other

- One year limited warranty; extended warranty to three years
- User manual and demonstration video make it easy to become a HAL expert
- Tutorial workshops available
- Stay abreast of the latest in patient simulation by joining the HAL online community at www.gaumard.com/HAL
- Patents pending; HAL is a registered trademark of Gaumard Scientific
- HAL support personnel available weekdays 8:00 a.m. to 4:30 p.m. ET

Gaumard® Scientific
14700 SW 136th Street
Miami, Florida 33196-5691

Toll Free USA
8:00 a.m. - 4:30 p.m. ET
Monday - Friday
800.882.6655

Worldwide
305.971.3790 / 305.666.8548

Fax
305.667.6085
Twenty-four hours

Web Site
www.gaumard.com

E-mail
sima@gaumard.com

©Gaumard Scientific, 2004
All rights reserved

Newborn HAL® Neonatal Simulator with total mobility

GAUMARD SCIENTIFIC

www.gaumard.com