

- 589 8. Transfer the contents of both flasks to a 50 mL conical tube. Repeat step 7 with an
590 additional 5 mL 1X PBS and transfer to the 50 mL conical tube.
- 591 9. Immediately add 20 mL estrogen-free DMEM to the 50 mL conical tube to inhibit
592 further cellular digestion by residual trypsin.
- 593 10. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
594 for an additional 5 minutes.
- 595 11. Aspirate media from pellet and re-suspend it in 4 mL estrogen-free DMEM,
596 drawing the pellet repeatedly through a 25 mL serological pipette to break up
597 clumps of cells.

598 At this point, cells are ready to be divided into the ongoing tissue culture and estrogen-free
599 conditioning groups.

600 9.2.1 Ongoing Tissue Culture Maintenance

- 601 1. Add 20 mL RPMI to two T150 flasks.
- 602 2. Add 220 µl G418 to the RPMI in the T150 flasks
- 603 3. Add 1 mL of cell suspension from 9.2 step 11 to each flask.
- 604 4. Place T150 flasks in tissue culture incubator (see conditions in Section 9.0) and
605 grow to 80% to 90% confluence (approximately 48 to 72 hrs).
- 606 5. Tissue culture medium may need to be changed 24 hours after addition of G418 to
607 remove cells that have died because they do not express reporter plasmid.
- 608 6. G418 does not need to be added to the flasks a second time.
- 609 7. Repeat Section 9.2 steps 1-11 for ongoing tissue culture maintenance.

610 9.2.2 Conditioning in Estrogen-free Medium

- 611 1. Add 20 mL estrogen-free DMEM to two T150 flasks.
- 612 2. Add 150 µL G418 to the estrogen-free DMEM in the T150 flasks.
- 613 3. Add 1 mL of cell suspension from Section 9.2 step 11 to each flask.


- 614 4. Tissue culture medium may need to be changed 24 hours after addition of G418 to
615 remove cells that have died because they do not express reporter plasmid.
- 616 5. G418 does not need to be added to the flasks a second time.
- 617 6. Place the T150 flasks in the incubator (see conditions in **Section 9.0**) and grow to
618 80% to 90% confluence (approximately 48 to 72 hrs).

619 9.2.3 Plating Cells Grown in Estrogen-free DMEM for Experimentation

- 620 1. Remove the T150 flasks that have been conditioned in estrogen-free DMEM for
621 48 to 72 hours from the incubator.
- 622 2. Aspirate the medium, then rinse the cells with 5 mL 1X PBS.
- 623 3. Aspirate 1X PBS, then add 3 mL 1X trypsin to the flasks, gently swirling the flask
624 to coat all cells with the trypsin.
- 625 4. Place the flasks in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 626 5. Detach cells by hitting the side of the flask sharply against the palm or the heel of
627 the hand.
- 628 6. Confirm cell detachment by examination under an inverted microscope. If cells
629 have not detached, return the flask to the incubator for 2 additional minutes, then
630 hit the flask again.
- 631 7. After cells have detached, add 5 mL 1X PBS and transfer the suspended cells
632 from the T150 flask to a 50 mL conical tube. Add an additional 5 mL PBS to the
633 flask, gently swirl around the flask, and then transfer to the 50 mL conical tube.
- 634 8. Immediately add 20 mL estrogen-free DMEM to each conical tube to inhibit
635 further cellular digestion by residual trypsin.
- 636 9. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
637 for an additional 5 minutes.
- 638 10. Aspirate the media from the pellet and re-suspend it in 20 mL DMEM, drawing
639 the pellet repeatedly through a 25 mL serological pipette to break up any clumps
640 of cells.

- 641 11. Pipette 15 μ L of the cell suspension into the “v” shaped slot on the
642 hemocytometer. Ensure that the solution covers the entire surface area of the
643 hemocytometer grid, and allow cells to settle before counting.
- 644 12. Using 100x magnification, view the counting grid.
- 645 13. The counting grid on the hemocytometer consists of nine sections, four of which
646 are counted (upper left, upper right, lower left, and lower right, see **Figure 9-1**
647 Each section counted consists of four by four grids. Starting at the top left and
648 moving clockwise, count all cells in each of the four by four grids. Some cells
649 will be touching the outside borders of the square, but only count those that touch
650 the top and right borders of the square. This value is then used in the calculation
651 below to get to the desired concentration of 200,000 cells/mL.

Figure 9-1 Hemocytometer Counting Grid.



652

The volume of each square is 10^{-4} mL, therefore:

653

$\text{Cells/mL} = (\text{average number per grid}) \times 10^{-4} \text{ mL} \times 1/(\text{starting dilution})$.

654

655.

Starting dilution: 20 mL (for T150 flasks)

656

657 Harvested cells for a T150 flask are suspended in 20 mL of estrogen-free DMEM and sampled

658 for determination of concentration of cells/mL.

659 Example Calculation:

660 • Grids 1, 2, 3, and 4 are counted and provide the following data:

661 ○ 50, 51, 49, and 50: average number of cells per grid is equal to 50.

662 Cells/mL = 50 cells per grid ÷ 10^{-4} mL volume of grid = 50×10^4 cells/mL (or 500,000
663 cells/mL)

664 Total # of Cells Harvested = 500,000 cells/mL x 20 mL

665 Desired Concentration (or Concentration_{Final}) = 200,000 cells/mL

666 Formula: (Concentration_{Final} x Volume_{Final} = Concentration_{Initial} x Volume_{Initial})

667

668 Concentration_{Final} = 200,000 cells/mL

669 Concentration_{Initial} = 500,000 cells/mL

670 Volume_{Initial} = 20 mL

671 Volume_{Final} – to be solved for.

672

673 Therefore: 200,000 cells/mL x Volume_{Final} = 500,000 cells/mL x 20 mL

674 Solving for Volume_{Final} we find = 50 mL

675

676 Therefore, add 30 mL of DMEM Growth media to the cell suspension for a total volume of 50
677 mL, which will yield the desired concentration of 200,000 cells/mL for plating.

678 14. This dilution scheme will give a concentration of 200,000 cells/mL. 200 µL of
679 this cell suspension is used for each well of a 96-well plate (i.e., 40,000 cells per
680 well).

681 15. Remove a 96-well plate from its sterile packaging. Use a repeater pipetter to
682 pipette 200 µL of cell suspension into each well except the outside ring of wells.

683 16. Use a repeater pipetter to pipette 200 µL of estrogen-free DMEM to the outer
684 wells of the 96-well plate.

685 17. Incubate plate(s) in an incubator (see conditions in **Section 9.0**) for a minimum of
686 24 hours, but no longer than 48 hours before dosing.

687 Two T150 flasks containing cells at 80% to 90% confluence will typically yield sufficient cells
688 to fill four 96-well plates (not including the perimeter wells).

689 **10.0 PREPARATION OF TEST SUBSTANCES**

690 The solvent used for dissolution of test substances is 100% DMSO. All test substances should be
691 allowed to equilibrate to room temperature before being dissolved and diluted. Test substance
692 solutions (except for reference standards and controls) should not be prepared in bulk for use in
693 subsequent tests. Test substances are to be used within 24 hours of preparation. Solutions should
694 not have noticeable precipitate or cloudiness.

695 All information on weighing, solubility testing, and calculation of final concentrations for test
696 substances, reference standards and controls is to be recorded in the study notebook.

697 **10.1 Determination of Test Substance Solubility**

- 698 1. Prepare a 100 mg/mL solution of the test substance in 100% DMSO in a 15 mL
699 conical tube.
- 700 2. Vortex to mix.
- 701 3. If the substance does not go into solution after vigorous vortexing, sonicate for 10
702 min.
- 703 4. If the test substance does not dissolve at 100 mg/mL, add an additional 9 mL of
704 DMSO to the conical tube (10 mg/mL), then vortex and/or sonicate for 10 min. as
705 before.
- 706 5. If the test substance does not dissolve at 10 mg/mL solution, prepare a 1 mg/mL
707 solution in a 15 mL conical tube. Vortex and/or sonicate for 10 min. as before.

708 6. If the test substance does not dissolve at 1 mg/mL, add an additional 9 mL of
709 DMSO to the conical tube (100 µg/mL), then vortex and/or sonicate for 10 min.
710 as before.

711 7. Continue testing, using 1/10 less substance in each subsequent attempt until test
712 substance is solubilized in DMSO.

713 Once the test substance has fully dissolved in 100% DMSO, the solubility of the test substance
714 must be determined in the 1% DMSO/99% estrogen-free DMEM mixture used for LUMI-
715 CELL® ER testing.

716 8. Add 4 µL of the highest concentration of the test substance/DMSO solution to a
717 13 mm test tube.

718 9. Add 400 µL estrogen-free DMEM to the test tube and vortex gently,

719 10. If cloudiness or precipitate develop, vortex for up to 10 minutes.

720 11. If vortexing does not dissolve test substance, sonicate test substance for up to 10
721 minutes.

722 12. If test substance has visible precipitate or is cloudy return to **10.1 step 7** to try the
723 next lower concentration for the test substance.

724 The Testing Facility shall forward the results from the solubility tests assay to the SMT through
725 the designated contacts in electronic format and hard copy upon completion of testing.

726 **10.2 Preparation of Reference Standards, Control and Test Substances**

727 All “dosing solutions” of test substance concentrations are to be expressed as µg/mL in the study
728 notebook and in all laboratory reports.

729 All information on preparation of test substances, reference standards and controls is to be
730 recorded in the study notebook.

731 10.2.1 Preparation of Reference Standard and Positive Control Stock Solutions

732 Stock solutions of E2 and methoxychlor are prepared in 100% DMSO and stored at room
733 temperature for up to three years or until the expiration date listed in the certificate of analysis
734 for that substance.

735 10.2.1.1 E2 Stock Solution

736 The final concentration of the E2 stock solution is 1.0×10^{-2} µg/mL. Prepare the E2 stock as
737 shown in **Table 10-1**.

738 **Table 10-1 Preparation of E2 Stock Solution**

Step #	Action	DMSO	E2 Concentration
1	Make a 10 mg/mL stock solution in 100% DMSO in a 4mL vial.		10 mg/mL
2	Transfer 10 µL E2 solution from Step #1 to a new 4 mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	100 µg/mL
3	Transfer 10 µL E2 solution from Step #2 to a new 4mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	1 µg/mL
4	Transfer 10 µL E2 solution from Step #3 to a 13 mm test tube to create the working solution.	Add 990 µL of 100% DMSO. Vortex to mix.	0.1 µg/mL

739 10.2.1.2 Methoxychlor Stock Solution

740 The final concentration of the methoxychlor stock solution is 313 µg/mL.

741 To prepare the methoxychlor stock solution, proceed as follows:

- 742 1. Make a 10 mg/mL stock solution of Methoxychlor in 100% DMSO in a 4 mL
743 vial.
- 744 2. Remove 94 µL of the methoxychlor solution and place it in a new 4 mL vial.
- 745 3. Add 2.906 mL of 100% DMSO to the 4mL vial and gently vortex to mix.

746 10.2.2 Preparation of Reference Standard and Positive Control Dosing Solutions for Range
747 Finder Testing

748 Range finder testing is conducted on 96-well plates using three concentrations of E2 in duplicate
749 as the reference standard. Six replicate wells are used for the DMSO control. All wells on the 96
750 well plate are used during range finder testing.

751 Store dosing solutions at room temperature. Use within 24 hours of preparation.

752

752 10.2.2.1 *Preparation of E2 Reference Standard Dosing Solutions for Range Finder Testing*
753 In preparation for making E2 dosing solutions, label two sets of three glass 13 mm test tubes
754 with the numbers one through three and place them in a test tube rack. Tube number 1 will
755 contain the highest concentration of E2 (Table 10-2).

756 **Table 10-2 Preparation of E2 Reference Standard Dosing Solution for Range Finder**
757 **Testing**

Tube Number	E2	Estrogen-free DMEM ¹	Final Volume	Final E2 Concentration
1	4 µl of 1.0×10^{-2} µg/mL working solution	400 µL	404 µL	1.00×10^{-3} µg/mL
2	8 µL of 1.0×10^{-3} µg/mL from Tube #1	400 µL	408 µL	2.00×10^{-5} µg/mL
3	1 µL of 1.0×10^{-3} µg/mL from Tube#1	4000 µL	4001 µL	1.00×10^{-6} µg/mL

758 ¹Vortex all tubes to mix media and E2.

759

760 10.2.2.2 *Preparation of DMSO Control Dosing Solution for Range Finder Testing*

761 1. Add 4 µL of 100% DMSO to six 13 mm tubes (solvent/negative controls).

762 2. Add 400 µL of estrogen-free DMEM to each tube and vortex vigorously.

763 10.2.3 Preparation of Test Substance Dosing Solutions for Range Finder Testing

764 Range finder experiments are used to determine the concentrations of test substance to be used
765 during comprehensive testing. Agonist range finding for coded substances consists of six point,
766 logarithmic serial dilutions run in duplicate.

767 Label two sets of six glass 13 mm test tubes with the numbers 1 through 6 and place them in a
768 test tube rack. Tube number 1 will contain the highest concentration of test substance (Table 10-
769 3).

770

770 **Table 10-3 Preparation of Test Substance Dosing Solutions for Range Finder Testing**

Tube Number	100% DMSO	Test Substance ¹	Transfer ²	Estrogen-free DMEM ³	Final Volume
1	-	4 µL of test substance solution from Section 10.1 step 14	4 µL	400 µL	404 µL
2	90 µL	10 µL of test substance solution from Section 10.1 step 14	4 µL	400 µL	404 µL
3	90 µL	10 µL from Tube #2	4 µL	400 µL	404 µL
4	90 µL	10 µL from Tube #3	4 µL	400 µL	404 µL
5	90 µL	10 µL from Tube #4	4 µL	400 µL	404 µL
6	90 µL	10 µL from Tube #5	4 µL	400 µL	404 µL

771 ¹Vortex tubes #2 through 5 before removing test substance solution to place in the next tube in the series.

772 ²Transfer test substance/DMSO solutions to a new set of 13 mm test tubes.

773 ³Vortex all tubes to mix media and test substance solution.

774

775 Determination of whether a substance is positive in range finder testing and selection of starting
776 concentrations for comprehensive testing will be discussed in **Section 12.0**.

777 10.2.4 Preparation of Reference Standard and Positive Control Dosing Solutions for
778 Comprehensive Testing

779 Comprehensive testing is conducted on 96-well plates using 10 concentrations of E2 in duplicate
780 as the reference standard. Four replicate wells for the DMSO control and three replicate wells for
781 the methoxychlor control are included on each plate.

782 Store dosing solutions at room temperature. Use within 24 hours of preparation.

783 10.2.4.1 *Preparation of E2 Reference Standard Dosing Solutions for Comprehensive Testing*

784 In preparation for making E2 double serial dilutions, label two sets of 11 glass 13 mm test tubes
785 with the numbers 1 through 11 and place them in a test tube rack. Tube number 1 will contain the
786 highest concentration of E2 (**Table 10-4**).

787

787 **Table 10-4 Preparation of E2 Reference Standard Dosing Solution for Comprehensive**
788 **Testing**

Tube Number	100% DMSO	E2 ¹	Discard	Estrogen-free DMEM ²	Final Volume
1	-	4 µL of 1.0 x 10 ⁻² µg/mL working solution	-	400 µL	404 µL
2	4 µL	4 µL of 1.0 x 10 ⁻² µg/mL stock solution	-	400 µL	404 µL
3	4 µL	4 µL from Tube #2	-	400 µL	404 µL
4	4 µL	4 µL from Tube #3	-	400 µL	404 µL
5	4 µL	4 µL from Tube #4	-	400 µL	404 µL
6	4 µL	4 µL from Tube #5	-	400 µL	404 µL
7	4 µL	4 µL from Tube #6	-	400 µL	404 µL
8	4 µL	4 µL from Tube #7	-	400 µL	404 µL
9	4 µL	4 µL from Tube #8	Discard Tube #9	-	-
10	4 µL	4 µL from Tube #9	-	400 µL	404 µL
11	4 µL	4 µL from Tube #10	Remove and discard 4 µL from Tube #11	400 µL	404 µL

789 ¹Vortex tubes #2 through 11 before removing E2 solution to place in the next tube in the series.

790 ²Vortex all tubes to mix media and E2.

791

792 10.2.4.2 *Preparation of Methoxychlor Control Dosing Solution for Comprehensive Testing*

793 1. Add 4 µL of the 313 µg/mL methoxychlor to three separate 13 mm tubes.

794 2. Add 400 µL of estrogen-free DMEM to each tube and vortex vigorously.

795 10.2.4.3 *Preparation of DMSO Control Dosing Solution for Comprehensive Testing*

796 1. Add 4 µL of 100% DMSO to four 13 mm tubes (solvent/negative controls).

797 2. Add 400 µL of estrogen-free DMEM to each tube and vortex vigorously.

798 10.2.5 Preparation of Test Substance Dosing Solutions for Comprehensive Testing

799 Comprehensive testing experiments are used to determine whether a substance possesses ER

800 agonist activity in the LUMI-CELL® ER test method. Agonist comprehensive testing for coded

801 substances consists of 11 point, double serial dilutions, with each concentration tested in
802 triplicate wells of the 96-well plate.

803 Start the 11-point serial dilution series at a single log dilution higher than the concentration
804 giving the highest adjusted RLU value during the range finder (e.g., if the highest adjusted RLU
805 value occurred at a concentration of 0.01 mg/mL, start the serial dilution at 0.1 mg/mL).

806 Label three sets of 11 glass 13 mm test tubes with the numbers 1 through 11 and place them in a
807 test tube rack. Tube number 1 will contain the highest concentration of test substance (**Table 10-**
808 **4**).

809 **Table 10-4 Preparation of Test Substance Dosing Solutions for Comprehensive Testing**

Tube Number	100% DMSO	Test Substance ¹	Discard	Estrogen-free DMEM ²	Final Volume
1	-	4 µL of highest concentration of test substance solution	-	400 µL	404 µL
2	4 µL	4 µL of highest concentration of test substance solution	-	400 µL	404 µL
3	4 µL	4 µL from Tube #2	-	400 µL	404 µL
4	4 µL	4 µL from Tube #3	-	400 µL	404 µL
5	4 µL	4 µL from Tube #4	-	400 µL	404 µL
6	4 µL	4 µL from Tube #5	-	400 µL	404 µL
7	4 µL	4 µL from Tube #6	-	400 µL	404 µL
8	4 µL	4 µL from Tube #7	-	400 µL	404 µL
9	4 µL	4 µL from Tube #8	-	400 µL	404 µL
10	4 µL	4 µL from Tube #9	-	400 µL	404 µL
11	4 µL	4 µL from Tube #10	4 µL	400 µL	404 µL

810 ¹Vortex tubes #2 through 11 before removing test substance solution to place in the next tube in the
811 series.

812 ²Vortex all tubes to mix media and test substance solution.

813

814 11.0 GENERAL PROCEDURES FOR THE TESTING OF CODED SUBSTANCES

815 Range finder experiments are used to determine the concentrations of test substance to be used
816 during comprehensive testing. Comprehensive testing experiments are used to determine whether
817 a substance possesses ER agonist activity in the LUMI-CELL® ER test method.

818 General procedures for range finder and comprehensive testing are nearly identical. For specific
819 details (such as plate layout) of range finder testing see **Section 12.0**. For specific details of
820 comprehensive testing, see **Section 13.0**.

821 **11.1 Application of Reference Standard, Controls, and Test Substances**

- 822 1. Remove the 96-well plates from the incubator, inspect them using an inverted
823 microscope. Only use plates in which the cells in all wells receive a score of 1
824 according to **Table 11-1**.
- 825 2. Remove medium by inverting the plate onto blotter paper. Gently tap plate against
826 the bench surface to remove residual liquid trapped in the wells.
- 827 3. Add 200 µL of medium, reference standard, control, or test substance to each well
828 (see **Sections 12.0** and **13.0** for specific plate layouts).
- 829 4. Return plates to incubator and incubate (see **Section 9.0** for details) for 19 to 24
830 hours to allow maximal induction of luciferase activity in the cells.

831 **11.1.1 Preparation of Excel® Data Analysis Template**

- 832 1. In Excel®, open a new “AgICCVAMTemplate” and save it with the appropriate
833 project name as indicated in the NICEATM Style Guide.
- 834 2. Add appropriate information regarding the assay to the “Compound Tracking”
835 tab.
- 836 3. Enter substance testing information to the “List” page (i.e., Project /Sample ID,
837 Concentration, and Comments (or compound name). This should populate the
838 “Template”, “Compound Mixing” and “Visual Inspection” tabs with the
839 appropriate information for the experiment.
- 840 4. Save the newly named project file.
- 841 5. Print out either the “List” or “Template” page for help with dosing the 96-well
842 plate. Sign and date the print out and store in study notebook.

843

843 **11.2 Visual Evaluation of Cell Viability**

- 844 1. 19 to 24 hours after dosing the plate, remove the plate from the incubator and
845 remove the media from the wells by inverting the plate onto blotter paper. Gently
846 tap plate against the bench surface to remove residual liquid trapped in the wells.
- 847 2. Use a repeat pipetter to add 50 µL 1X PBS to all wells. Immediately remove PBS
848 by inversion.
- 849 3. Using an inverted microscope, inspect all of the wells used in the 96-well plate
850 and record the visual observations using the scores in **Table 11-1**.

851 **Table 11-1 Visual Observation Scoring**

Viability Score	Brief Description ¹
1	Normal Cell Morphology and Cell Density
2	Altered Cell Morphology and/or Small Gaps between Cells
3	Altered Cell Morphology and/or Large Gaps between Cells
4	Few (or no) Visible Cells
1P	Score of 1 with Precipitate
2P	Score of 2 with Precipitate
3P	Score of 3 with Precipitate
4P	Score of 4 with Precipitate
5P	Unable to View Cells Due to Precipitate

852 ¹Reference micrographs will be provided by NICEATM.
853

854 **11.3 Lysis of Cells for LUMI-CELL® ER**

- 855 1. Apply the reflective white backing tape to the bottom of the 96-well plate (this
856 will increase the effectiveness of the luminometer).
- 857 2. Add 30 µL 1X lysis reagent to the assay wells and place the 96-well plate on an
858 orbital shaker for one minute.
- 859 3. Remove plate from shaker and measure luminescence (as described in **Section**
860 **11.5**).

861 **11.4 CellTiter-Glo® Assessment of Cell Viability**

862 When considered necessary, a quantitative evaluation of cell viability will be performed with the
863 Promega CellTiter-Glo® assay system. CellTiter-Glo® uses luminescence as an indicator of the

864 number of cells per plate and therefore must be conducted in parallel with the LUMI-CELL® ER
865 test method (i.e., both test methods cannot be conducted on the same plate).

- 866 1. Dose and incubate cells under the same conditions as for LUMI-CELL® ER.
- 867 2. Remove plates from incubator and discard the medium by inverting the plate onto
868 blotter paper. Gently tap plate against the bench surface to remove residual liquid
869 trapped in the wells.
- 870 3. Use a repeat pipettor to add 50 µL 1X PBS to all assay wells. Immediately
871 remove PBS by inversion.
- 872 4. Examine all wells used under an inverted microscope. Make notes of any well
873 with codes described in **Table 11-1**.
- 874 5. Place white backing tape on the bottom of the 96-well plate.
- 875 6. Add 100 µL estrogen-free DMEM to each well containing cells.
- 876 7. Add 100 µL CellTiter-Glo® reagent to each well containing cells.
- 877 8. Place plate on an orbital shaker for 1 minute to induce cell lysis.
- 878 9. Incubate (see **Section 9.0** for details) for 10 min.
- 879 10. Measure luminescence promptly. Do *not* add luciferase reagent to CellTiter-Glo®
880 plates.

881 **11.5 Measurement of Luminescence**

882 Luminescence is measured in the range of 300 to 650 nm, using an injecting luminometer and
883 with software that controls the injection volume and measurement interval. Light emission from
884 each well is expressed as RLU per well. The luminometer output is saved as raw data in an
885 Excel® spread sheet. A hard copy of the luminometer raw data should be signed, dated and stored
886 in the study notebook.

887

888

888 **11.6 Data Analysis**

889 LUMI-CELL® ER uses an Excel® spreadsheet to collect and adjust the RLU values obtained
890 from the luminometer and a GraphPad Prism® template to analyze and graph data. The Excel®
891 spreadsheet subtracts background luminescence (average DMSO solvent control RLU value)
892 from test substance, reference standard and control RLU values. Plate induction is calculated
893 using these corrected RLU values. Test substance, reference standard, and control RLU values
894 are then adjusted relative to the highest E2 reference standard RLU value, which is set to 10,000.
895 After adjustment, values are transferred to GraphPad Prism® for data analysis and graphing.

896 11.6.1 Correction and Adjustment of Luminometer Data

897 The following steps describe the procedures required to populate the Excel® spreadsheet that has
898 been configured to collect and adjust the RLU values obtained from the luminometer.

- 899 1. Open the raw data file and the corresponding experimental Excel® spreadsheet
900 from **Section 11.1.1**.
- 901 2. Copy the raw data using the Excel® copy function, then paste the copied data into
902 cell C22 of the “RAW DATA” tab in the experimental Excel® spreadsheet using
903 the Paste Special – Values command. This position corresponds to position B2 in
904 the table labeled Table 1 in this tab.
- 905 3. Examine the DMSO data in Table 1 of the Excel® spreadsheet to determine
906 whether there are any potential outliers. See **Section 11.6.2** for further explanation
907 of outlier determinations.
- 908 4. If an outlier is identified, perform the following steps to remove the outlier from
909 calculations:
- 910 ▪ correct the equation used to calculate DMSO background in Table 1
911 (e.g., if outlier is located in cell G24, adjust the calculation in cell H42 to
912 read =AVERAGE(F24:H24,I24))
 - 913 ▪ then correct the equation used to calculate the average DMSO value in
914 Table 2 (e.g., following the above example, adjust cell M44 to read
915 =AVERAGE(F36:H36,I36))

- 916 ▪ then correct the equation used to calculate the standard deviation of the
917 DMSO value in Table 2 (e.g., following the above example, adjust cell
918 M45 to read =STDEV(F36,H36,I36))
- 919 5. Excel® will automatically subtract the background (the average DMSO control
920 value) from all of the RLU values in Table 1 and populate Table 2 with these
921 adjusted values.
- 922 6. To calculate plate induction, identify the cell containing the E2 replicate that has
923 the highest RLU value and the cell containing the RLU values for the same
924 concentration in the corresponding E2 replicate (e.g., if the highest RLU value for
925 E2 is located in cell E23, the corresponding cell would be E22).
- 926 7. Click into cell D16 and enter the cell number from the previous step into the
927 numerator.
- 928 8. Click on the “ER Agonist Report” tab.
- 929 9. The data for the E2 reference standard, methoxychlor, and DMSO replicates
930 populate the left portion (columns A – F) of the spreadsheet. The data is
931 automatically placed in an Excel® graph.
- 932 10. To set the highest RLU value for the reference standard to 10,000 RLU, go to cell
933 D2 of “ER Agonist Report” tab and check the formula contained within that cell.
934 The divisor should be the cell number of the cell containing the highest averaged
935 E2 RLU value (column E).
- 936 11. Copy the data into GraphPad Prism® for determination of outliers, graphing, and
937 analysis as indicated in the NICEATM Prism® user’s guide.

938 11.6.2 Determination of Outliers

939 The Study Director will use good statistical judgment for determining “unusable” wells that will
940 be excluded from the data analysis and will provide an explanation in the study notebook for any
941 excluded data. This judgment for data acceptance will include Q-test analysis.

942 The formula for the Q test is:

943 Outlier – Nearest Neighbor

944 Range (Highest – Lowest)

944 where the outlier is the value proposed for exclusion, the nearest neighbor is the value closest to
945 the outlier, and the range is the range of the three values.

946 If the value of this ratio is greater than 0.94 (the Q value for the 90% confidence interval for a
947 sample size of three) or 0.76 (the Q value for the 90% confidence interval for a sample size of
948 four), the outlier may be excluded from data analysis. For E2 reference standard replicates
949 (sample size of two), any adjusted RLU value for a replicate at a given concentration of E2 is
950 considered and outlier if its value is more than 20% above or below the adjusted RLU value for
951 that concentration in the historical database.

952 11.6.3 Acceptance Criteria

953 Acceptance or rejection of a test is based on evaluation of reference standard and control results
954 from each experiment conducted on a 96-well plate. Results are compared to quality controls
955 (QC) for these parameters derived from the historical database, which are summarized below.

- 956 • Induction: Plate induction, as measured by dividing the averaged highest E2
957 reference standard RLU value by the averaged DMSO control RLU value, must
958 be greater than three-fold.
- 959 • Reference standard results: Calculated E2 reference standard EC₅₀ values must be
960 within 2.5 times the standard deviation of the historical database EC₅₀ mean
961 value.
- 962 • Solvent control results: Solvent control RLU values must be within 2.5 times the
963 standard deviation of the historical solvent control mean RLU value.
- 964 • Positive control results: Methoxychlor control RLU values must be within 2.5
965 times the standard deviation of the historical methoxychlor control mean RLU
966 value.


967 An experiment that fails any single acceptance criterion will be discarded and repeated.

968

968 11.6.4 Calculation of Relative EC₅₀ Values

969 Relative EC₅₀ values are calculated as directed in the NICEATM Prism® users guide. The
970 concentration curve must have a sigmoidal shape and reach saturation at the highest and lowest
971 concentrations tested (**Figure 11-1**) for Prism® to calculate a valid relative EC₅₀ value.

972 **Figure 11-1 Example Concentration Curve for Calculation of Relative EC₅₀ Value**



973

974 The mathematical model used by Prism® to calculate a relative EC₅₀ value uses the Hill function
975 as described in **Section 6.0**.

976 11.6.5 Calculation of Absolute EC₅₀ Values


977 Calculate an absolute EC₅₀ value for all substances that have a positive response that reaches
978 50% of the E2 reference standard response (Eli Lilly 2005).

979 To calculate the absolute EC₅₀:

- 980 1. Convert the test substance data to the percentage of the maximum agonist
981 response for E2, the reference substance, excluding any negative values.
- 982 2. Find the concentration of the test substance that is 50% of the maximum E2
983 response.

984

984 **Figure 11-2 Example Concentration Curve for Calculation of Absolute EC₅₀ Value**




985

986 **12.0 RANGE FINDER TESTING**

987 Agonist range finding for coded substances consists of six point, logarithmic serial dilutions,
988 with each concentration tested in duplicate wells of the 96-well plate. **Figure 12-1** contains a
989 template for the plate layout to be used in agonist range finder testing.

990

990 **Figure 12-1 Agonist Range Finder Test Plate Layout**



- - Three Point E2 Reference Standard – Replicate 1
- - Three Point E2 Reference Standard – Replicate 2
- - DMSO Control (1% v/v)
- - Range Finder for Sample #2
- - Range Finder for Sample #1
- - Range Finder for Sample #3
- - Range Finder for Sample #4
- - Range Finder for Sample #5
- - Range Finder for Sample #6

991 Evaluate whether range finder experiments have met the acceptance criteria (see **Section 11.6.3**)
 992 and graph the data as described in the NICEATM Prism® users guide.
 993 To determine starting concentrations for comprehensive testing use the following criteria: