

191 **1.0 PURPOSE**

192 This protocol is designed to evaluate coded test substances for potential estrogen receptor (ER)
193 antagonist activity using the LUMI-CELL® ER assay.

194 **2.0 SPONSOR**

195 The National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative
196 Toxicological Methods (NICEATM), P.O. Box 12233 Research Triangle Park, NC 27709

197 William S. Stokes, DVM, DACLAM

198 Rear Admiral, U.S. Public Health Service

199 Chief Veterinary Officer, USPHS

200 Director, NICEATM

201 National Institute of Environmental Health Sciences, NIH, DHHS

202 Bldg. 4401, Room 3129, MD EC-14

203 79 T.W. Alexander Drive

204 Research Triangle Park, NC 27709

205 Phone: 919-541-7997

206 Fax: 919-541-0947

207 Email: stokes@niehs.nih.gov

208

209 Raymond Tice, Ph.D.

210 Deputy Director, NICEATM

211 National Institute of Environmental Health Sciences

212 MD EC-17, P.O. Box 12233

213 Research Triangle Park, NC 27709

214 Phone: 919-541-4482

215 FAX: 919-541-0947
216 Email: tice@niehs.nih.gov
217
218 David Allen, Ph.D.
219 Principal Investigator
220 ILS, Inc./Contractor supporting NICEATM
221 National Institute of Environmental Health Sciences
222 MD EC-17, P.O. Box 12233
223 Research Triangle Park, NC 27709
224 Phone: 919-316-4587
225 FAX: 919-541-0947
226 Email: allen7@niehs.nih.gov
227
228 Frank Deal, M.S.
229 Staff Toxicologist
230 ILS, Inc./Contractor supporting NICEATM
231 National Institute of Environmental Health Sciences
232 MD EC-17, P.O. Box 12233
233 Research Triangle Park, NC 27709
234 Phone: 919-316-4587
235 FAX: 919-541-0947
236 Email: dealf@niehs.nih.gov
237
238

238 Patricia Ceger, M.S.
239 Project Coordinator/Technical Writer
240 ILS, Inc./Contractor supporting NICEATM
241 National Institute of Environmental Health Sciences
242 MD EC-17, P.O. Box 12233
243 Research Triangle Park, NC 27709
244 Phone: 919-316-4556
245 Fax: 919-541-0947
246 E-Mail: cegerp@niehs.nih.gov

247

248 **2.1 Substance Inventory and Distribution Management**

249 Cynthia Smith, Ph.D.
250 Chemistry Resources Group Leader
251 National Institute of Environmental Health Sciences
252 MD EC-06, P.O. Box 12233
253 Research Triangle Park, NC 27709
254 Phone: 919-541-3473
255

256 **3.0 DEFINITIONS**

- 257 • **Dosing Solution:** The test substance, control substance, or reference standard
258 solution which is to be placed into the tissue culture wells for experimentation.
- 259 • **Raw Data:** Raw data includes information that has been collected but not
260 formatted or analyzed, and consists of the following:
- 261 ○ Data recorded in the Study Notebook

- 262 ○ Computer printout of initial luminometer data
- 263 ○ Other data collected as part of GLP compliance, e.g.:
- 264 ▪ Equipment logs and calibration records
- 265 ▪ Test substance and tissue culture media preparation logs
- 266 ▪ Cryogenic freezer inventory logs
- 267 • **Soluble:** Test substance exists in a clear solution without visible cloudiness or
- 268 precipitate.
- 269 • **Study Notebook:** The study notebook contains recordings of all activities related
- 270 to the conduct of the LUMI-CELL® ER TA antagonist assay.
- 271 • **Test Substances:** Substances supplied to the testing laboratories that are coded
- 272 and distributed such that only the Project Officer, Study Management Team
- 273 (SMT), and the Substance Inventory and Distribution Management have
- 274 knowledge of their true identity. The test substances will be purchased, aliquoted,
- 275 coded, and distributed by the Supplier under the guidance of the NIEHS/NTP
- 276 Project Officer and the SMT.

277 4.0 **TESTING FACILITY AND KEY PERSONNEL**

278 4.1 **Testing Facility**

279 Xenobiotic Detection Systems, Inc. (XDS), 1601 E. Geer St., Durham, NC 27704

280 4.2 **Key Personnel**

- 281 • Study Director: John Gordon, Ph.D.
- 282 • Laboratory Technician(s): Cynthia Matherly
- 283 • Scientific Advisor: Mike Denison, Ph.D.
- 284 • Quality Assurance Director: Andrew Chu
- 285 • Safety Manager: Andrew Chu
- 286 • Facility Management: George Clark, Ph.D., Robert Clark, Tina Ginter

287 • Laboratory Director: Andrew Chu

288 5.0 IDENTIFICATION OF TEST AND CONTROL SUBSTANCES

289 5.1 Test Substances

290 Test substances are coded and will be provided to participating laboratories by the Substance
291 Inventory and Distribution Management team.

292 5.2 Controls

293 Controls for the ER antagonist protocol are as follows:

294 *Vehicle control (dimethyl sulfoxide [DMSO]):* 1% v/v dilution of DMSO (CASRN 67-68-5)
295 diluted in tissue culture media.

296 *Ral/E2 reference standard for range finder testing:* Three concentrations (1.25×10^{-1} , 1.25×10^{-3} ,
297 and 5.00×10^{-5} $\mu\text{g/mL}$) of raloxifene HCl (raloxifene, CASRN 84449-90-1) plus a fixed
298 concentration (2.5×10^{-5} $\mu\text{g/mL}$) of 17β -estradiol (E2, CASRN: 50-28-2).

299 *Ral/E2 reference standard for comprehensive testing:* A serial dilution of raloxifene HCl
300 (raloxifene, CASRN 84449-90-1) plus a fixed concentration (2.5×10^{-5} $\mu\text{g/mL}$) of 17β -estradiol
301 (E2, CASRN: 50-28-2), consisting of nine concentrations of Ral/E2 in duplicate wells.

302 *E2 control:* 17β -estradiol, CASRN: 50-28-2, 2.5×10^{-5} $\mu\text{g/mL}$ in tissue culture media used as a
303 base line negative control.

304 *Flavone/E2 Control:* Flavone, CASRN 525-82-6, 25 $\mu\text{g/mL}$, with E2 2.5×10^{-5} $\mu\text{g/mL}$ in tissue
305 culture media used as a weak positive control.

306 6.0 OVERVIEW OF GENERAL PROCEDURES FOR ANTAGONIST TESTING

307 All experimental procedures are to be carried out under aseptic conditions and all solutions,
308 glassware, plastic ware, pipettes, etc., shall be sterile. All methods and procedures shall be
309 documented in the study notebook.

310 Antagonist range finder testing is conducted on 96-well plates using three concentrations of
311 Ral/E2 (1.25×10^{-1} , 1.25×10^{-3} , and 5.00×10^{-5} $\mu\text{g/mL}$ raloxifene with 2.50×10^{-5} $\mu\text{g/mL}$ E2) in

312 duplicate as the reference standard, with three replicate wells for the E2 and DMSO controls.
313 Range finder testing uses all wells of the 96-well plate.
314 Comprehensive testing is conducted on 96-well plates using nine concentrations of Ral/E2 in
315 duplicate as the reference standard (**Table 6-1**). Three replicate wells for the DMSO control,
316 Flavone/E2 and E2 controls are included on each plate. In order to avoid edging effects¹, wells
317 on the perimeter of the plate are not used for experiments. These wells should contain tissue
318 culture media only (i.e., no cells).

319 **Table 6-1 Concentrations of Ral/E2 Reference Standard**
320 **Used in Range Finder and Comprehensive Testing**

Raloxifene Concentrations ¹	E2 Concentrations
1.25×10^{-2}	2.5×10^{-5}
6.25×10^{-3}	2.5×10^{-5}
3.13×10^{-3}	2.5×10^{-5}
1.56×10^{-3}	2.5×10^{-5}
7.81×10^{-4}	2.5×10^{-5}
3.91×10^{-4}	2.5×10^{-5}
1.95×10^{-4}	2.5×10^{-5}
9.77×10^{-5}	2.5×10^{-5}
4.88×10^{-5}	2.5×10^{-5}

321 ¹Concentrations are presented in µg/mL.
322

323 Visual observations for cell viability are conducted for all experimental plates just prior to
324 LUMI-CELL® ER evaluation. CellTiter-Glo® based cell viability testing (when used) is
325 conducted concurrently in parallel plates on the same day, as outlined in **Section 11.4**.

326 Luminescence data, measured in relative light units (RLUs), is corrected for background
327 luminescence by subtracting the mean RLU value of the vehicle control (DMSO) wells from the
328 RLU measurements for each of the other wells of the 96-well plate. Data is then transferred into

¹Edging effects are variations in response seen in the outermost wells in a cell culture plate. These variations are believed to be due to variations in temperature, evaporation, etc., that may occur in these wells that would ultimately affect cellular growth and health (ICCVAM, 2001).

329 Excel® data management spreadsheets and GraphPad PRISM® 4.0 statistical software, graphed,
330 and evaluated for a positive or negative response as follows:

- 331 • A response is considered positive for antagonist activity when the average
332 adjusted RLU for a given concentration is less than the mean RLU value minus
333 three times the standard deviation for the E2 control.
- 334 • Any luminescence at or above this threshold is considered a negative response.

335 For substances that are positive at one or more concentrations, the concentration of test substance
336 that causes a half-maximal response (the relative IC₅₀) is calculated using a Hill function
337 analysis. The Hill function is a four-parameter logistic mathematical model relating the
338 substance concentration to the response (typically following a sigmoidal curve) using the
339 equation below

$$340 \quad Y = \text{Bottom} + \frac{\text{Top} - \text{Bottom}}{1 + 10^{(\log \text{IC}_{50} - X) \text{HillSlope}}}$$

341 where Y = response (i.e., relative light units); X = the logarithm of concentration; Bottom = the
342 minimum response; Top = the maximum response; log IC₅₀ = the logarithm of X as the response
343 midway between Top and Bottom; and HillSlope describes the steepness of the curve. The model
344 calculates the best fit for the Top, Bottom, HillSlope, and IC₅₀ parameters. See **Section 13.6.5** for
345 more details.

346 Acceptance or rejection of a test is based on evaluation of reference standard and control results
347 from each experiment conducted on a 96-well plate. Results for these controls are compared to
348 historical results compiled in the historical database, as seen in **Section 16.0**.

349 **6.1 Range Finder Testing**

350 Antagonist range finding for coded substances consists of a six point, logarithmic serial dilution
351 using duplicate wells per concentration. Concentrations for comprehensive testing are selected
352 based on the response observed in range finder testing. If necessary, a second range finder test
353 can be conducted to clarify the optimal concentration range to test (see **Section 14.0**).

354

354 **6.2 Comprehensive Testing**

355 Comprehensive antagonist testing for coded substances consists of 11 point, double serial
356 dilutions, with each concentration tested in triplicate wells of the 96-well plate. Three separate
357 experiments are conducted for comprehensive testing on three separate days, except during
358 Phases III and IV of the validation effort, in which comprehensive testing experiments are
359 conducted once (see **Section 15.0**).


360 **7.0 MATERIALS FOR LUMI-CELL® ER ANTAGONIST TESTING**

361 This section provides the materials needed to conduct LUMI-CELL® testing, with associated
362 brand names/vendors² in brackets.

363 **7.1 BG1Luc4E2 Cells:**

364 Human ovarian cancer cell line stably transfected with a plasmid containing an estrogen response
365 element (**Figure 7-1**) [XDS].

366 **Figure 7-1 pGudLuc7.ERE Plasmid.**



367

368

²Brand names and vendors should not be considered an endorsement by the U.S. Government or any member of the U.S. Government; such information is provided as examples.

368 **7.2 Technical Equipment:**

369 All technical equipment may be obtained from Fisher Scientific International, Inc. (Liberty Lane
370 Hampton, NH, USA 03842). Equivalent technical equipment from another commercial source
371 can be used.

- 372 • Analytical balance (Cat. No. 01-910-320)
- 373 • Berthold Orion 1 Microplate Luminometer [Berthold CatNo.: Orion 1 MPL3] or
374 equivalent and dedicated computer
- 375 • Biological safety hood, class II, and stand (Cat. No. 16-108-99)
- 376 • Centrifuge (low speed, tabletop with swinging bucket rotor) (Cat. No. 04-978-50
377 centrifuge, and 05-103B rotor)
- 378 • Combustion test kit (CO₂ monitoring) (Cat. No. 10-884-1)
- 379 • Drummond diaphragm pipetter (Cat. No. 13-681-15)
- 380 • Freezers, -20°C (Cat. No. 13-986-150), and -70°C (Cat. No. 13-990-86)
- 381 • Hand tally counter (Cat. No. 07905-6)
- 382 • Hemocytometer, cell counter (Cat. No. 02-671-5)
- 383 • Light microscope, inverted (Cat. No. 12-561-INV)
- 384 • Light microscope, upright (Cat. No. 12-561-3M)
- 385 • Liquid nitrogen flask (Cat. No. 11-675-92)
- 386 • Micropipetter, repeating (Cat. No. 21-380-9)
- 387 • Pipetters, air displacement, single channel (0.5 –10µl (Cat. No. 21-377-191), 2 –
388 20 µl (Cat. No. 21-377-287), 20 – 200 µl (Cat. No. 21-377-298), 200 - 1000 µl
389 (Cat. No. 21-377-195))
- 390 • Refrigerator/freezer (Cat. No. 13-986-106A)
- 391 • Shaker for 96-well plates (Cat. No. 14-271-9)
- 392 • Sodium hydroxide (Cat. No. 5318-500)

- 393 • Sonicating water bath (Cat. No. 15-335-30)
 - 394 • Tissue culture incubator with CO₂ and temperature control (Cat. No. 11-689-4)
 - 395 • Vacuum pump with liquid trap (side arm Erlenmeyer) (Cat. No. 01-092-29)
 - 396 • Vortex mixer (Cat. No. 12-814)
- 397 Equipment should be maintained and calibrated as per GLP guidelines and individual laboratory
398 SOPs.

399 7.3 Reference Standard, Controls, and Tissue Culture Supplies

400 All tissue culture reagents must be labeled to indicate source, identity, storage conditions and
401 expiration dates. Tissue culture solutions must be labeled to indicate concentration, stability
402 (where known), and preparation and expiration dates.

403 Equivalent tissue culture media and sera from another commercial source can be used, but must
404 first be tested as described in **Section 17.0** to determine suitability for use in this test method.

405 The following are the necessary tissue culture reagents and possible sources based on their use in
406 the pre-validation studies:

- 407 • 17 β -estradiol (CAS RN: 50-28-2) [Sigma-Aldrich, Cat. No. E8875]
- 408 • CellTiter-Glo® Luminescent Cell Viability Assay [Promega Cat. No. G7572]
- 409 • Cryovial, 2 mL (Corning Costar) [Fisher Scientific Cat. No. 03-374-21]
- 410 • Culture tube 13 x 100mm (case) [Thomas Scientific Cat. No.: 10009186R38]
- 411 • Culture tube, 15 mL conical (Corning Costar) [Fisher Scientific Cat. No. 05-538-
412 51]
- 413 • Culture tube, 50 mL conical (Corning Costar) [Fisher Scientific Cat. No. 05-
414 526C]
- 415 • DMSO, U.S.P. analytical grade. [Sigma-Aldrich, Cat. No. 34869-100ML]
- 416 • Dulbecco's Modification of Eagle's Medium (DMEM), containing 4.5 g/L
417 glucose, with sodium pyruvate, without phenol red or L-glutamine
418 [Mediatech/Cellgro, Cat. No. 17-205-CV]

- 419 • Fetal Bovine Serum [Mediatech/Cellgro Cat. No. MT 35-010-CV]
- 420 • Fetal Bovine Serum, charcoal/dextran treated, triple 0.1 µm sterile filtered
- 421 [Hyclone, Cat. No. SH30068.03]
- 422 • Flavone (CASRN: 525-82-6) [Sigma-Aldrich, Cat. No. F2003]
- 423 • Gentamycin Sulfate (G418), 50 mg/mL [Mediatech/Cellgro Cat. No. 30-234-CR]
- 424 • L-glutamine, 29.2 mg/mL [Cellgro; Cat. No. 25005-CI]
- 425 • Luciferase Assay System (10-Pack) [Promega Cat. No. E1501]
- 426 • Lysis Solution 5X [Promega, Cat. No. E1531]
- 427 • Penicillin/streptomycin solution, 5000 I.U. penicillin, 5000 µg/mL streptomycin
- 428 [Cellgro, Cat. No. 30-001-CI].
- 429 • Phosphate buffered saline (PBS, 1X) without calcium and magnesium [Cellgro,
- 430 Cat. No. 21-040-CV]
- 431 • Pipettes, serological: 2.0 mL [Sigma-Aldrich, Cat. No. P1736], 5.0 mL [Sigma-
- 432 Aldrich, Cat. No. P1986], 25 mL [Sigma-Aldrich, Cat. No. P2486]
- 433 • Raloxifene (CASRN 84449-90-1) [Sigma-Aldrich Cat. No. R1402]
- 434 • RPMI 1640 medium, containing L-glutamine [Mediatech, Cat. No. 10-040-CV]
- 435 • Tissue culture flasks (Corning-Costar): 25 cm² (T25) [Fisher Cat. No. 10-126-28];
- 436 75 cm² (T75) [Fisher Cat. No. 10-126-37]; and 150 cm² (T150) [Fisher Cat. No.
- 437 10-126-34]
- 438 • Tissue culture plates (Corning-Costar): 96-well [Thomas Scientific Cat. No.
- 439 6916A05]
- 440 • Trypsin (10X), 2.5% in Hank's balanced salt solution (HBSS), without calcium
- 441 and magnesium, without phenol red [Cellgro, Cat. No. 25-054-CI].

442 All reagent lot numbers and expiration dates must be recorded in the study notebook.

443

444 **8.0 PREPARATION OF TISSUE CULTURE MEDIA AND SOLUTIONS**

445 All tissue culture media and media supplements must be quality tested before use in experiments
446 (see **Section 15.0**).

447 **8.1 RPMI 1640 Growth Medium (RPMI)**

448 RPMI 1640 is supplemented with 0.9% Pen-Strep and 8.0% FBS to make RPMI growth medium
449 (RPMI).

450 Procedure for one 549 mL bottle:

- 451 1. Remove FBS from -70°C freezer, and Pen-Strep from -20°C freezer and allow to
452 equilibrate to room temperature.
- 453 2. Add 44 mL of FBS and 5 mL Pen-Strep to the bottle of RPMI 1640.
- 454 3. Label RPMI bottle as indicated in **Section 7.3**

455 *Store at 2-8°C for no longer than six months or until the shortest expiration date of any media*
456 *component.*

457 **8.2 Estrogen-Free DMEM Medium**

458 DMEM is supplemented to contain 4.5% charcoal/dextran treated FBS, 1.9% L-glutamine, 0.9%
459 Pen-Strep.

460 Procedure for one 539 mL bottle:

- 461 1. Remove charcoal/dextran treated FBS from -70°C freezer, and L-glutamine and
462 Pen-Strep from -20°C freezer and allow to equilibrate to room temperature.
- 463 2. Add 24 mL of charcoal/dextran treated FBS, 10 mL L-glutamine, and 5 mL Pen-
464 Strep to one 500 mL bottle of DMEM.
- 465 3. Label estrogen-free DMEM bottle as indicated in **Section 7.3**

466 *Store at 2-8°C for no longer than six months or until the shortest expiration date of any media*
467 *component.*

468

468 **8.3 1X Trypsin Solution**

469 1X Trypsin solution is prepared by dilution from a 10X premixed stock solution. The 10X stock
470 solution should be stored in 10 mL aliquots in a -20°C freezer.

471 Procedure for making 100 mL of 1X trypsin:

- 472 1. Remove a 10mL aliquot of 10X trypsin from -20°C freezer and allow to
473 equilibrate to room temperature.
- 474 2. Aliquot 1 mL Trypsin (10X) along with 9 mL of 1X PBS into ten 15 mL
475 centrifuge tubes.
- 476 3. Label 1X trypsin aliquots as indicated in **Section 7.3**

477 *1X Trypsin should be stored at -20°C.*

478 **8.4 1X Lysis Solution**

479 Lysis solution is prepared by dilution from a 5X premixed stock solution. Both the 5X and 1X
480 solutions can be repeatedly freeze-thawed.

481 The procedure for making 10 mL of 1X lysis solution:

- 482 1. Thaw the 5X Promega Lysis solution and allow it to reach room temperature.
- 483 2. Remove 2 mL of 5X solution and place it in a 15 mL conical centrifuge tube.
- 484 3. Add 8 mL of distilled, de-ionized water to the conical tube.
- 485 4. Cap and shake gently until solutions are mixed.

486 *Store at -20°C for no longer than 1 year from receipt.*

487 **8.5 Reconstituted Luciferase Reagent**

488 Luciferase reagent consists of two components, luciferase buffer and lyophilized luciferase
489 substrate.

490 For long-term storage, unopened containers of the luciferase buffer and lyophilized luciferase
491 substrate can be stored at -70°C for up to six months.

492 To reconstitute luciferase reagent:

- 493 1. Remove luciferase buffer and luciferase substrate from -70°C freezer and allow
494 them to equilibrate to room temperature.
- 495 2. Add 10 mL of luciferase buffer solution to luciferase substrate container and swirl
496 or vortex to mix, the Luciferase substrate should readily go into solution.
- 497 3. Luciferase substrate should readily go into solution.
- 498 4. After solutions are mixed aliquot to a 15mL centrifuge tube.
- 499 5. Store complete solution at -20°C.

500 Reconstituted luciferase reagent is stable for 1 month at -20°C.

501 **8.6 Reconstituted CellTiter-Glo® Reagent**

502 CellTiter-Glo® reagent consists of two components, CellTiter-Glo® buffer and lyophilized
503 CellTiter-Glo® substrate.

504 For long term storage, unopened containers of the CellTiter-Glo® buffer and lyophilized
505 CellTiter-Glo® substrate can be stored at -70°C for up to one year.

506 To reconstitute CellTiter-Glo® reagent:

- 507 1. Remove CellTiter-Glo® buffer and CellTiter-Glo® substrate from -70°C freezer
508 and allow them to equilibrate to room temperature.
- 509 2. Add CellTiter-Glo® buffer solution to CellTiter-Glo® substrate container and
510 swirl or vortex gently to mix; the CellTiter-Glo® substrate should readily go into
511 solution.
- 512 3. CellTiter-Glo® substrate should readily go into solution.
- 513 4. After solutions are mixed aliquot to a 15mL centrifuge tube.
- 514 5. Store complete solution at -20°C.

515 Reconstituted CellTiter-Glo® reagent is stable for up to 3 months at -20°C.

516 **9.0 OVERVIEW OF PROPAGATION AND EXPERIMENTAL PLATING OF**
517 **BG1Luc4E2 CELLS**

518 The BG1Luc4E2 (BG-1) cells are stored in liquid nitrogen in 2 mL cryovials. BG-1 cells are
519 grown as a monolayer in tissue culture flasks in a dedicated tissue culture incubator at 37°C ±
520 1°C, 90% ± 5% humidity, and 5.0% ± 1% CO₂/air. The cells should be examined on a daily basis
521 during working days under an inverted phase contrast microscope, and any changes in
522 morphology and adhesive properties must be noted in the study notebook.

523 Two T150 flasks containing cells at 80% to 90% confluence will usually yield a sufficient
524 number of cells to fill four 96-well plates for use in experiments.

525 **9.1 Procedures for Thawing Cells and Establishing Tissue Cultures**

526 Warm all tissue culture media and solutions to room temperature by placing them under the
527 tissue culture hood several hours before use.

528 All tissue culture media, media supplements, and tissue culture plasticware must be quality
529 tested before use in experiments (**Section 17.0**).

530 **9.1.1 Thawing Cells**

- 531 1. Remove a cryovial of frozen BG-1 cells from the liquid nitrogen flask.
- 532 2. Facilitate rapid thawing by loosening the top slightly (do not remove top) to
533 release trapped gasses and retightening it. Roll vial between palms.
- 534 3. Use a micropipette to transfer cells to a 50 mL conical centrifuge tube.
- 535 4. Rinse cryovial twice with 1X PBS and add PBS rinse material to the conical tube.
- 536 5. Add 20 mL of RPMI to the conical tube.
- 537 6. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
538 for an additional 5 minutes.
- 539 7. Aspirate media from pellet and re-suspend it in 5 mL RPMI, drawing the pellet
540 repeatedly through a 2.0 mL serological pipette to break up any clumps of cells.
- 541 8. Transfer cells to a T25 flask, place them in an incubator (see conditions in
542 **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).

543 9.1.2 Establishing Tissue Cultures

544 Once cells have reached 80% to 90% confluence, transfer the cells to a T75 flask by performing,
545 for example, the following steps:

- 546 1. Remove the T25 flask from the incubator.
- 547 2. Aspirate the RPMI, then add 5 mL 1X PBS.
- 548 3. Aspirate 1X PBS, then add 1 to 2 mL 1X trypsin to the T25 flask, gently swirling
549 the flask to coat all cells with the trypsin.
- 550 4. Place the flask in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 551 5. Detach cells by hitting the side of the flask sharply against the palm or heel of the
552 hand.
- 553 6. Confirm cell detachment by examination under an inverted microscope. If cells
554 have not detached, return the flask to the incubator for an additional 2 minutes,
555 then hit the flask again.
- 556 7. After cells have detached, add 5 mL PBS, and transfer the suspended cells to a 50
557 mL centrifuge tube. Wash the flask one additional time with 5 mL PBS.
- 558 8. Immediately add 20 mL RPMI to the conical tube to inhibit further cellular
559 digestion by residual trypsin.
- 560 9. Pellet the cells by centrifugation, as described in **Section 9.1.1**, and re-suspend the
561 cells in 10 mL RPMI medium.
- 562 10. Draw the pellet repeatedly through a 25 mL serological pipette to break up
563 clumps of cells
- 564 11. Transfer cells to a T75 flask, then place the flask in an incubator (see conditions
565 in **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).

566 When cells have reached 80% to 90% confluency, transfer them into a T150 flask by performing,
567 for example, the following steps:

- 568 12. Remove the T75 flask from the incubator, aspirate the old media and add 5 mL
569 1X PBS.

- 570 13. Aspirate 1X PBS, add 2 mL of 1X trypsin to the flask, and place it in an incubator
571 (see conditions in **Section 9.0**) for 5 to 10 min.
- 572 14. Repeat **steps 5** through **11** in **Section 9.1.2**, re-suspending the pellet in 20 mL of
573 RPMI.
- 574 15. Transfer cells to a T150 flask and place it in the incubator (see conditions in
575 **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).
- 576 16. Remove the T150 flask from the incubator.
- 577 17. Aspirate the RPMI and add 5 mL 1X PBS.
- 578 18. Aspirate 1X PBS and add 3 mL 1X trypsin to the T150 flask, making sure that the
579 cells are coated with the trypsin.
- 580 19. Incubate cells in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 581 20. Detach cells by hitting the side of the flask sharply against the palm or heel of the
582 hand.
- 583 21. Confirm cell detachment by examination under an inverted microscope. If cells
584 have not detached, return the flask to the incubator for an additional 2 minutes,
585 then hit the flask again.
- 586 22. After cells have detached, add 5mL 1X PBS and transfer the suspended cells from
587 the T150 flask to a 50 mL conical tube. Add an additional 5 mL PBS to the flask,
588 then transfer to the 50 mL conical tube.
- 589 23. Immediately add 20 mL RPMI to the conical tube to inhibit further cellular
590 digestion by residual trypsin.
- 591 24. Centrifuge at 1000 x g for eight minutes. If a pellet of cells has not formed,
592 centrifuge for an additional 5 minutes.
- 593 25. Aspirate the media from the pellet and re-suspend it in 40 mL RPMI, drawing the
594 pellet repeatedly through a 25 mL serological pipette to break up any clumps of
595 cells.

596 26. Transfer 20 mL of cell suspension to each of two T150 flasks, place them in an
597 incubator (see conditions in **Section 9.0**) and grow to 80% to 90% confluence
598 (approximately 48 to 72 hrs).

599 **9.2 Ongoing Tissue Culture Maintenance, Conditioning in Estrogen-free Medium,**
600 **and Plating Cells for Experimentation**

601 The following procedure is used to condition the BG1Luc4E2 cells to an estrogen-free
602 environment prior to plating the cells in 96-well plates for analysis of estrogen dependent
603 induction of luciferase activity.

604
605 To start the tissue culture maintenance and estrogen-free conditioning, split the two T150 culture
606 flasks into four T150 flasks. Two of these flasks will be used for continuing tissue culture and
607 will use the RPMI media mentioned above. The other two flasks will be cultured in estrogen-free
608 DMEM for experimental use. Extra care must be taken to avoid contaminating the estrogen-free
609 cells with RPMI.

- 610 1. Remove both T150 flasks from the incubator.
- 611 2. Aspirate the medium and rinse the cells with 5 mL 1X PBS.
- 612 3. Aspirate 1X PBS, then add 3 mL 1X trypsin to the flasks, gently swirling the flask
613 to coat all cells with the trypsin.
- 614 4. Incubate cells in the incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 615 5. Detach cells by hitting the side of the flask sharply against the palm or heel of the
616 hand.
- 617 6. Confirm cell detachment by examination under an inverted microscope. If cells
618 have not detached, return the flask to the incubator for an additional 2 minutes,
619 then hit the flask again.
- 620 7. After cells have detached, add 5 mL 1X PBS to the first T150 flask and transfer
621 the suspended cells to the second T150 flask.
- 622 8. Transfer the contents of both flasks to a 50 mL conical tube. Repeat step 7 with an
623 additional 5 mL 1X PBS and transfer to the 50 mL conical tube.

- 624 9. Immediately add 20 mL estrogen-free DMEM to the 50 mL conical tube to inhibit
625 further cellular digestion by residual trypsin.
- 626 10. Centrifuge at 1000 x g for eight minutes. If a pellet of cells has not formed,
627 centrifuge for an additional 5 minutes.
- 628 11. Aspirate media from pellet and re-suspend it in 4 mL estrogen-free DMEM,
629 drawing the pellet repeatedly through a 25 mL serological pipette to break up
630 clumps of cells.

631 At this point, cells are ready to be divided into the ongoing tissue culture and estrogen-free
632 conditioning groups.

633 9.2.1 Ongoing Tissue Culture Maintenance

- 634 1. Add 20 mL RPMI to two T150 flasks.
- 635 2. Add 220 µL G418 to the RPMI in the T150 flasks
- 636 3. Add 1 mL of cell suspension from **Section 9.2 step 11** to each flask.
- 637 4. Place T150 flasks in tissue culture incubator (see conditions in **Section 9.0**) and
638 grow to 80% to 90% confluence (approximately 48 to 72 hrs).
- 639 5. Tissue culture medium may need to be changed 24 hours after addition of G418 to
640 remove cells that have died because they do not express reporter plasmid.
- 641 6. G418 does not need to be added to the flasks a second time.
- 642 7. Repeat **Section 9.2 steps 1-11** for ongoing tissue culture maintenance.

643 9.2.2 Conditioning in Estrogen-free Medium

- 644 1. Add 20 mL estrogen-free DMEM to two T150 flasks.
- 645 2. Add 150 µL G418 to the estrogen-free DMEM in the T150 flasks.
- 646 3. Add 1 mL of cell suspension from **Section 9.2 step 11** to each flask.
- 647 4. Tissue culture medium may need to be changed 24 hours after addition of G418 to
648 remove cells that have died because they do not express reporter plasmid.
- 649 5. G418 does not need to be added to the flasks a second time.

650 6. Place the T150 flasks in the incubator (see conditions in **Section 9.0**) and grow to
651 80% to 90% confluence (approximately 48 to 72 hrs).

652 9.2.3 Plating Cells Grown in Estrogen-free DMEM for Experimentation

- 653 1. Remove the T150 flasks that have been conditioned in estrogen-free DMEM for
654 48 to 72 hours from the incubator.
- 655 2. Aspirate the medium, then rinse the cells with 5 mL 1X PBS.
- 656 3. Aspirate 1X PBS, then add 3 mL 1X trypsin to the flasks, gently swirling the flask
657 to coat all cells with the trypsin.
- 658 4. Place the flasks in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 659 5. Detach cells by hitting the side of the flask sharply against the palm or the heel of
660 the hand.
- 661 6. Confirm cell detachment by examination under an inverted microscope. If cells
662 have not detached, return the flask to the incubator for 2 additional minutes, then
663 hit the flask again.
- 664 7. After cells have detached, add 5 mL 1X PBS and transfer the suspended cells
665 from the T150 flask to a 50 mL conical tube. Add an additional 5 mL PBS to the
666 flask, then transfer to the 50 mL conical tube.
- 667 8. Immediately add 20 mL estrogen-free DMEM to each conical tube to inhibit
668 further cellular digestion by residual trypsin.
- 669 9. Centrifuge at 1000 x g for eight minutes. If a pellet of cells has not formed,
670 centrifuge for an additional 5 minutes.
- 671 10. Aspirate off the media from the pellet and re-suspend it in 20 mL DMEM,
672 drawing the pellet repeatedly through a 25 mL serological pipette to break up any
673 clumps of cells.
- 674 11. Pipette 15 μ L of the cell suspension into the “v” shaped slot on the
675 hemocytometer. Ensure that the solution covers the entire surface area of the
676 hemocytometer grid, and allow cells to settle before counting.