


- 677 12. Using 100x magnification, view the counting grid.
- 678 13. The counting grid on the hemocytometer consists of nine sections, four of which
679 are counted (upper left, upper right, lower left, and lower right, see **Figure 9-1**).
680 Each section counted consists of four by four grids. Starting at the top left and
681 moving clockwise, count all cells in each of the four by four grids. Some cells
682 will be touching the outside borders of the square, but only count those that touch
683 the top and right borders of the square. This value is then used in the calculation
684 below to get to the desired concentration of 200,000 cells/mL.

Figure 9-1 Hemocytometer Counting Grid.



685

The volume of each square is 10^{-4} mL, therefore:

686

Cells/mL = (average number per grid) x 10^4 mL. x 1/(starting dilution).

687

688

Starting dilution: 20mL (for T150 flasks)

689

690 Harvested cells for a T150 flask are suspended in 20 mL of estrogen-free DMEM and sampled
691 for determination of concentration of cells/mL.

692

693 Example Calculation:

694

- Grids 1, 2, 3, and 4 are counted and provide the following data:

695

- 50, 51, 49, and 50: average number of cells per grid is equal to 50.

696 Cells/mL = 50 cells per grid ÷ 10⁻⁴ mL volume of grid = 50 X 10⁴ cells/mL (or 500,000
697 cells/mL)

698 Total # of Cells Harvested = 500,000 cells/mL x 20 mL

699 Desired Concentration (or Concentration_{Final})= 200,000 cells/mL

700 Formula: (Concentration_{Final} x Volume_{Final} = Concentration_{Initial} x Volume_{Initial})

701

702 Concentration_{Final} = 200,000 cells/mL

703 Concentration_{Initial} = 500,000 cells/mL

704 Volume_{Initial} = 20 mL

705 Volume_{Final} – to be solved for.

706

707 Therefore: 200,000 cells/mL x Volume_{Final} = 500,000 cells/mL x 20 mL

708 Solving for Volume_{Final} we find = 50 mL

709

710 Therefore, add 30 mL of DMEM Growth media to the cell suspension for a total volume of 50
711 mL, which will yield the desired concentration of 200,000 cells/mL for plating.

712 14. This dilution scheme will give a concentration of 200,000 cells/mL. 200 µL of
713 this cell suspension is used for each well of a 96-well plate (i.e., 40,000 cells per
714 well).

715 15. Remove a 96-well plate from its sterile packaging. Use a repeater pipetter to
716 pipette 200 µL of cell suspension into each well except the outside ring of wells.

717 16. Use a repeater pipetter to pipette 200 µL of estrogen-free DMEM to the outer
718 wells of the 96-well plate.

719 17. Incubate plate(s) in an incubator (see conditions in **Section 9.0**) for a minimum of
720 24 hours, but no longer than 48 hours before dosing.

721 Two T150 flasks containing cells at 80% to 90% confluence will typically yield sufficient cells
722 to fill four 96-well plates (not including the perimeter wells).

723 **10.0 PREPARATION OF TEST SUBSTANCES**

724 The solvent used for dissolution of test substances is 100% DMSO. All test substances should be
725 allowed to equilibrate to room temperature before being dissolved and diluted. Test substance
726 solutions (except for reference standards and controls) should not be prepared in bulk for use in
727 subsequent tests. Test substances are to be used within 24 hours of preparation. Solutions should
728 not have noticeable precipitate or cloudiness.

729 All information on weighing, solubility testing, and calculation of final concentrations for test
730 substances, reference standards and controls is to be recorded in the study notebook.

731 **10.1 Determination of Test Substance Solubility**

- 732 1. Prepare a 200 mg/mL solution of the test substance in 100% DMSO in a 15 mL
733 conical tube.
- 734 2. Vortex to mix.
- 735 3. If the substance does not go into solution after vigorous vortexing, sonicate for 10
736 min.
- 737 4. If the 200 mg/mL solution has visible precipitate or is cloudy, add an additional 9
738 mL of DMSO to the conical tube (20 mg/mL), then vortex and/or sonicate for 10
739 min. as before.
- 740 5. If the 20 mg/mL solution has visible precipitate or is cloudy, prepare a 2 mg/mL
741 solution in a 15 mL conical tube. Vortex and/or sonicate for 10 min. as before.
- 742 6. If the 2 mg/mL solution has visible precipitate or is cloudy, add an additional 9
743 mL of DMSO to the conical tube (200 µg/mL), then vortex and/or sonicate for 10
744 min. as before.
- 745 7. Continue testing, using 1/10 less substance in each subsequent attempt until the
746 test substance is solubilized in DMSO.

747 Once a solution of test substance has been obtained that does not have any visible precipitate or
748 cloudiness in 100% DMSO, the solubility of the test substance must be determined in the 1%
749 DMSO/99% estrogen-free DMEM mixture used for LUMI-CELL® ER testing.

750 8. Add 2 µL of the highest concentration of the test substance/DMSO solution to a
751 13 mm test tube.

752 9. Add 400 µL estrogen-free DMEM to the test tube and vortex gently,

753 10. If cloudiness or precipitate develop, vortex for up to 10 minutes.

754 11. If vortexing does not dissolve test substance, sonicate test substance for up to 10
755 minutes.

756 12. If test substance has visible precipitate or is cloudy return to **Section 10.1 step 7**
757 to try the next lower concentration for the test substance.

758 The Testing Facility shall forward the results from the solubility tests assay to the SMT through
759 the designated contacts in electronic format and hard copy upon completion of testing.

760 **11.0 PREPARATION OF REFERENCE STANDARD, CONTROL AND TEST**
761 **SUBSTANCE STOCK SOLUTIONS FOR RANGE FINDER AND**
762 **COMPREHENSIVE TESTING**

763 All information on preparation of test substances, reference standards and controls is to be
764 recorded in the study notebook.

765 **11.1 Preparation of Ral/E2 Stock Solutions**

766 E2 and raloxifene stocks are prepared separately and then combined into Ral/E2 stocks, which
767 are then used to prepare dosing solutions in **Section 12**.

768 **11.1.1 E2 Stock Solution**

769 The final concentration of the E2 stock solution is 5.0×10^{-3} µg/mL. Prepare the E2 stock as
770 shown in **Table 11-1**.

771

771 **Table 11-1 Preparation of E2 Stock Solution**

Step #	Action	DMSO	E2 Concentration
1	Make a 10 mg/mL stock solution in 100% DMSO in a 4mL vial.	-	10 mg/mL
2	Transfer 10 µL E2 solution from Step #1 to a new 4 mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	100 µg/mL
3	Transfer 10 µL E2 solution from Step #2 to a new 4mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	1 µg/mL
4	Transfer 100 µL E2 solution from Step #3 to a new glass container large enough to hold 15 mL.	Add 9.90 mL of 100% DMSO. Vortex to mix.	1.0×10^{-2} µg/mL
5	Transfer 5 mL E2 solution from Step #4 to a new glass container large enough to hold 15 mL	Add 5 mL of 100% DMSO. Vortex to mix.	5.0×10^{-3} µg/mL

772

773 11.1.2 Raloxifene Working Stock Solution

774 Prepare a 2.5 µg/mL raloxifene working stock solution as shown in **Table 11-2**.

775 **Table 11-2 Preparation of Raloxifene Working Stock Solution**

Step #	Action	DMSO	Raloxifene Concentration
1	Make a 10 mg/mL solution of raloxifene in a 4 mL glass vial.	-	1.0×10^4 µg/mL
2	Transfer 10 µL raloxifene solution from Step #1 to a new 4 mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	100 µg/mL
3	Transfer 150 µL raloxifene solution from Step #2 to a new 4 mL vial.	Add 2.850 mL of 100% DMSO. Vortex to mix.	5 µg/mL
4	Transfer 1.5 mL raloxifene solution from Step #3 to a new 13 mm test tube.	Add 1.5 mL of 100% DMSO. Vortex to mix.	2.5 µg/mL

776

777 **11.2 Ral/E2 Range Finder Stock**

778 11.2.1 Raloxifene Dilutions

779 Number three 4 mL vials with the numbers 1 to 3 and use the raloxifene solution prepared in

780 **Section 11.1.2** to make raloxifene dilutions as shown **Table 11-3**.

781

781 **Table 11-3 Preparation of Raloxifene Reference Standard for Range Finder Testing**

Step #	Action	DMSO	Raloxifene Concentration
1	Transfer 500 µL of the raloxifene working stock solution to a new 4 mL vial.	-	2.5 µg/mL
2	Transfer 50 µL of the raloxifene working stock solution to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	2.50×10^{-1} µg/mL
3	Transfer 2 µL of the raloxifene working stock solution to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	1.00×10^{-2} µg/mL

782

783 11.2.2 Ral/E2 Dilutions for Range Finder Stock:

784 Add 500 µL of the 5×10^{-3} µg/mL E2 solution prepared in **Section 11.1.1** to each of the
785 raloxifene dilution vials prepared in **Section 11.2.1**. Vortex each tube to mix. The final
786 concentrations for raloxifene and E2 are listed in **Table 11-4**.

787 **Table 11-4 Concentrations of Raloxifene and E2 in the**
788 **Ral/E2 Range Finder Stock Solution**

Tube #	Raloxifene (µg/ml)	E2 (µg/ml)
1	1.25	2.5×10^{-3}
2	1.25×10^{-1}	2.5×10^{-3}
3	5.00×10^{-3}	2.5×10^{-3}

789

790 **11.3 Ral/E2 Comprehensive Testing Stock**

791 11.3.1 Raloxifene Dilutions for Comprehensive Testing Stock

792 Use the raloxifene solution prepared in **Section 11.1.2** to make a 9 point serial dilution of
793 raloxifene as shown **Table 11-5**.

794 **Table 11-5 Preparation of Raloxifene 9-Point Serial Dilution**

Step #	Action	DMSO	Discard	Raloxifene Concentration
1	Transfer 500 µL of the raloxifene working stock solution to a new 4 mL vial.	-	-	2.5 µg/mL
2	Transfer 500 µL of the raloxifene working stock solution to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	1.25 µg/mL

Step #	Action	DMSO	Discard	Raloxifene Concentration
3	Transfer 500 µL raloxifene solution from Step #2 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	6.25×10^{-1} µg/mL
4	Transfer 500 µL raloxifene solution from Step #3 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	3.13×10^{-1} µg/mL
5	Transfer 500 µL raloxifene solution from Step #4 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	1.56×10^{-1} µg/mL
6	Transfer 500 µL raloxifene solution from Step #5 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	7.81×10^{-2} µg/mL
7	Transfer 500 µL raloxifene solution from Step #6 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	3.91×10^{-2} µg/mL
8	Transfer 500 µL raloxifene solution from Step #7 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	-	1.95×10^{-2} µg/mL
9	Transfer 500 µL raloxifene solution from Step #8 to a new 4 mL vial.	Add 500 µL of 100% DMSO. Vortex to mix.	Discard 500 µL from Tube #9	9.77×10^{-3} µg/mL

795

796 11.3.2 Ral/E2 Dilutions for Comprehensive Testing Stock:

797 Add 500 µL of the 5×10^{-3} µg/mL E2 solution prepared in Section 11.1.1 to each of the 9
798 raloxifene dilution vials (including the working stock solution in Tube #1). Vortex each tube to
799 mix. The final concentrations for raloxifene and E2 are listed in Table 11-6.

800 **Table 11-6 Concentrations of Raloxifene and E2 in the Ral/E2 Stock Solution**

Tube #	Raloxifene (µg/mL)	E2 (µg/mL)
1	1.25	2.5×10^{-3}
2	6.25×10^{-1}	2.5×10^{-3}
3	3.13×10^1	2.5×10^{-3}
4	1.56×10^{-1}	2.5×10^{-3}
5	7.81×10^2	2.5×10^{-3}
6	3.91×10^{-2}	2.5×10^{-3}
7	1.95×10^{-2}	2.5×10^{-3}
8	9.77×10^{-3}	2.5×10^{-3}
9	4.88×10^{-3}	2.5×10^{-3}

801

802 **11.4 Flavone/E2 Stock Solution**

803 To prepare the flavone/E2 stock solution, proceed as follows:

- 804 1. Prepare 1 mL of 10 mg/mL flavone
- 805 2. Add 1 mL of the 5×10^{-3} $\mu\text{g/mL}$ E2 (prepared as in **Section 11.1.1**) to the 10
- 806 mg/mL flavone. This will make a working solution of 5 mg/mL flavone with
- 807 2.5×10^{-3} $\mu\text{g/mL}$ E2.

808 **12.0 PREPARATION OF REFERENCE STANDARD, CONTROL AND TEST**

809 **SUBSTANCE DOSING SOLUTIONS FOR RANGE FINDER AND**

810 **COMPREHENSIVE TESTING**

811 **12.1 Preparation of Reference Standard and Control Dosing Solutions for Range**

812 **Finder Testing**

813 Range finder testing is conducted on 96-well plates using three concentrations of Ral/E2 in

814 duplicate as the reference standard. Three replicate wells for the DMSO, and E2 controls are

815 included on each plate.

816 All “dosing solutions” of test substance concentrations are to be expressed as $\mu\text{g/mL}$ in the study

817 notebook and in all laboratory reports.

818 Dosing solutions are to be used within 24 hours of preparation.

819 **12.1.1 Preparation of Ral/E2 Reference Standard Range Finder Dosing Solutions**

- 820 1. Label two sets of 13 mm glass tubes with the numbers 1 to 3.
- 821 2. Add 4 μL of Ral/E2 stock from tube #1 from **Section 11.2.2** to the two 13 mm
- 822 glass test tubes labeled #1.
- 823 3. Add 4 μL of Ral/E2 stock from tube #2 from **Section 11.2.2** to the two 13 mm
- 824 glass test tubes labeled #2. Repeat for tubes #3.
- 825 4. Add 400 μL of estrogen-free DMEM to each tube and vortex to mix.

826 **12.1.2 Preparation of DMSO Control Range Finder Dosing Solution**

- 827 1. Add 4 μL of 100% DMSO to three 13 mm glass test tubes.

828 2. Add 400 µL of estrogen-free DMEM to each tube and vortex to mix.

829 12.1.3 Preparation of E2 Control Range Finder Dosing Solution

830 1. Add 2 µL of the E2 stock from **Section 11.1.1** to three 13 mm glass test tube.

831 2. Add 400 µL of estrogen-free DMEM to each tube and vortex to mix.

832 **12.2 Preparation of Test Substance Dosing Solutions for Range Finder Testing**

833 Range finder experiments are used to determine the concentrations of test substance to be used
834 during comprehensive testing. Antagonist range finding for coded substances consists of six
835 point, logarithmic serial dilutions.

836 To prepare test substance dosing solutions:

837 1. Label two sets of six glass 13 mm test tubes with the numbers 1 through 6 and
838 place them in a test tube rack. Perform a serial dilution of test substance as shown
839 in **Table 12-1**.

840 **Table 12-1 Preparation of Test Substance Serial Dilution for Range Finder Testing**

Tube #	100% DMSO	Test Substance ¹	Final Volume
1	-	10 µL of test substance solution from Section 10.1	100 µL
2	90 µL	10 µL of test substance solution from Section 10.1	100 µL
3	90 µL	10 µL from Tube #2	100 µL
4	90 µL	10 µL from Tube #3	100 µL
5	90 µL	10 µL from Tube #4	100 µL
6	90 µL	10 µL from Tube #5	100 µL

841 ¹Vortex tubes #2 through 5 before removing test substance/DMSO solution to place in the next tube in the series.

842 ²Transfer test substance/DMSO solutions to a new set of 13 mm test tubes.

843 ²Vortex all tubes to mix media, test substance, and E2.

844

845 2. Transfer test substance/DMSO solutions to new tubes and add E2 as shown in
846 **Table 12-2**.

847

847 **Table 12-2 Addition of E2 to Test Substance Serial Dilution for Range Finder Testing**

Tube Number	Test Substance	E2	Remove	Estrogen-free DMEM ³	Final Volume
1	Transfer 4 µL of test substance from Tube #1 in Section 12.2 step 1 to a new tube	Add 4 µL of the 5 x 10 ⁻³ µg/mL E2 solution prepared in Section 11.1.1. Vortex to mix.	Remove and discard 4 µL from Tube #1	400 µL	404 µL
2	Transfer 4 µL of test substance from Tube #2 to a new tube	Add 4 µL of the 5 x 10 ⁻³ µg/mL E2 solution prepared in Section 11.1.1. Vortex to mix.	Remove and discard 4 µL from Tube #2	400 µL	404 µL
3	Transfer 4 µL of test substance from Tube #3 to a new tube	Add 4 µL of the 5 x 10 ⁻³ µg/mL E2 solution prepared in Section 11.1.1. Vortex to mix.	Remove and discard 4 µL from Tube #3	400 µL	404 µL
4	Transfer 4 µL of test substance from Tube #4 to a new tube	Add 4 µL of the 5 x 10 ⁻³ µg/mL E2 solution prepared in Section 11.1.1. Vortex to mix.	Remove 4 µL from Tube #4	400 µL	404 µL
5	Transfer 4 µL of test substance from Tube #5 to a new tube	Add 4 µL of the 5 x 10 ⁻³ µg/mL E2 solution prepared in Section 11.1.1. Vortex to mix.	Remove 4 µL from Tube #5	400 µL	404 µL
6	Transfer 4 µL of test substance from Tube #6 to a new tube	Add 4 µL of the 5 x 10 ⁻³ µg/mL E2 solution prepared in Section 11.1.1. Vortex to mix.	Remove 4 µL from Tube #6	400 µL	404 µL

848

849 Determination of whether a substance is positive in range finder testing and selection of starting
 850 concentrations for comprehensive testing will be discussed in Section 14.0.

851

851 **12.3 Preparation of Reference Standard and Control Dosing Solutions for**
852 **Comprehensive Testing**

853 Comprehensive testing is conducted on 96-well plates using nine concentrations of Ral/E2 in
854 duplicate as the reference standard. Three replicate wells for the DMSO, E2 and flavone/E2
855 controls are included on each plate.

856 All “dosing solutions” of test substance concentrations are to be expressed as $\mu\text{g/mL}$ in the study
857 notebook and in all laboratory reports.

858 Store dosing solutions at room temperature. Use within 24 hours of preparation.

859 **12.3.1 Preparation of Ral/E2 Reference Standard Dosing Solutions for Comprehensive**
860 **Testing**

861 In preparation for making Ral/E2 double serial dilutions, label two sets of nine glass 13 mm test
862 tubes with the numbers 1 through 9 and place them in a test tube rack. Tube number 1 will
863 contain the highest concentration of raloxifene (**Table 12-3**).

864 **Table 12-3 Preparation of Ral/E2 Reference Standard Dosing Solution**
865 **for Comprehensive Testing**

Tube Number	Ral/E2 Stock	Estrogen-free DMEM	Final Volume
1	4 μL of Tube #1 from Section 11.3.2	400 μL	404 μL
2	4 μL of Tube #2 from Section 11.3.2	400 μL	404 μL
3	4 μL of Tube #3 from Section 11.3.2	400 μL	404 μL
4	4 μL of Tube #4 from Section 11.3.2	400 μL	404 μL
5	4 μL of Tube #5 from Section 11.3.2	400 μL	404 μL
6	4 μL of Tube #6 from Section 11.3.2	400 μL	404 μL
7	4 μL of Tube #7 from Section 11.3.2	400 μL	404 μL
8	4 μL of Tube #8 from Section 11.3.2	400 μL	404 μL
9	4 μL of Tube #9 from Section 11.3.2	400 μL	404 μL

866

867 **12.3.2 Preparation of DMSO Control Comprehensive Testing Dosing Solution**

868 1. Add 4 μL of 100% DMSO to three 13 mm glass test tubes.

869 2. Add 400 µL of estrogen-free DMEM to each tube and vortex to mix.

870 12.3.3 Preparation of E2 Control Comprehensive Testing Dosing Solution

871 1. Add 2 µL of the E2 stock from **Section 11.1.1** to three 13 mm glass test tube.

872 2. Add 400 µL of estrogen-free DMEM to each tube and vortex to mix.

873 12.3.4 Preparation of Flavone/E2 Control Comprehensive Dosing Solution

874 1. Add 4 µL of flavone/E2 from **Section 11.4** to three 13 mm glass test tubes.

875 2. Add 400 µL of estrogen-free DMEM to each tube and vortex to mix.

876 **12.4 Preparation of Test Substance Dosing Solutions for Comprehensive Testing**

877 Comprehensive testing experiments are used to determine whether a substance possesses ER
878 antagonist activity in the LUMI-CELL® ER test method. Antagonist comprehensive testing for
879 coded substances consists of 11 point, double serial dilutions, with each concentration tested in
880 triplicate wells of the 96-well plate.

881 To prepare test substance dosing solutions for comprehensive testing:

882 1. Determine the concentration at which maximal antagonism occurs in the range
883 finder experiment. Start an 11-point serial dilution curve at a 20-fold higher
884 concentration than the concentration causing a maximal antagonist response (i.e.,
885 if the maximum antagonist response occurred at 0.01 mg/mL, start the serial
886 dilution curve at 0.2 mg/mL).

887 2. If there is no observable antagonist response in the range finder experiment, start
888 the serial dilution at the maximum soluble concentration.

889 3. Label three sets of 11 glass 13 mm test tubes with the numbers 1 through 11 and
890 place them in a test tube rack. Prepare three serial dilutions as indicated in **Table**
891 **12-4**.

892

892 **Table 12-4 Preparation of Test Substance Dosing Solutions for Comprehensive Testing**

Tube Number	100% DMSO	Test Substance ¹	Discard	E2 Testing Stock	Discard	Estrogen-free DMEM ²	Final Volume
1	-	4 µL of test substance solution from Section 10.2.4 step 1	-	4 µL	4 µL	400 µL	404 µL
2	4 µL	4 µL of test substance solution from Section 10.2.4 step 1	-	4 µL	4 µL	400 µL	404 µL
3	4 µL	4 µL from Tube #2	-	4 µL	4 µL	400 µL	404 µL
4	4 µL	4 µL from Tube #3	-	4 µL	4 µL	400 µL	404 µL
5	4 µL	4 µL from Tube #4	-	4 µL	4 µL	400 µL	404 µL
6	4 µL	4 µL from Tube #5	-	4 µL	4 µL	400 µL	404 µL
7	4 µL	4 µL from Tube #6	-	4 µL	4 µL	400 µL	404 µL
8	4 µL	4 µL from Tube #7	-	4 µL	4 µL	400 µL	404 µL
9	4 µL	4 µL from Tube #8	-	4 µL	4 µL	400 µL	404 µL
10	4 µL	4 µL from Tube #9	-	4 µL	4 µL	400 µL	404 µL
11	4 µL	4 µL from Tube #10	4 µL	4 µL	4 µL	400 µL	404 µL

893 ¹Vortex tubes #2 through 10 before removing test substance/DMSO solution to place in the next tube in
894 the series.

895 ²Vortex all tubes to mix media, test substance, and E2.

896

897 **13.0 GENERAL PROCEDURES FOR THE TESTING OF CODED SUBSTANCES**

898 Range finder experiments are used to determine the concentrations of test substance to be used
899 during comprehensive testing. Comprehensive testing experiments are used to determine whether
900 a substance possesses ER antagonist activity in the LUMI-CELL® ER test method.

901 General procedures for range finder and comprehensive testing are nearly identical. For specific
902 details (such as plate layout) of range finder testing see Section 14.0. For specific details of
903 comprehensive testing, see Section 15.0.

904

904 **13.1 Application of Reference Standard, Control and Test Substances**

- 905 1. Remove the 96-well plates (from **Section 9.2.3 step 18**) from the incubator;
906 inspect them using an inverted microscope. Only use plates in which the cells in
907 all wells receive a score of 1 according to **Table 11-1**.
- 908 2. Remove medium by inverting the plate onto blotter paper. Gently tap plate against
909 the bench surface to remove residual liquid trapped in the wells.
- 910 3. Add 200 µL of medium, reference standard, control or test substance to each well
911 (see **Sections 14.0** and **15.0** for specific plate layouts).
- 912 4. Return plates to incubator (see **Section 9.0** for details) for 19 to 24 hours to allow
913 maximal induction of luciferase activity in the cells.

914 **13.1.1 Preparation of Excel® Data Analysis Template**

- 915 1. In Excel®, open a new “AntICCVAMTemplate” and save it with the appropriate
916 project name as indicated in the NICEATM Style Guide.
- 917 2. Add appropriate information regarding the assay to the “Compound Tracking”
918 tab.
- 919 3. Enter substance testing information to the “List” page (i.e. Project /Sample ID,
920 Concentration, and Comments (or compound name). This should populate the
921 “Template”, “Compound Mixing” and “Visual Inspection” tabs with the
922 appropriate information for the experiment.
- 923 4. Save the newly named project file
- 924 5. Print out either the “List” or “Template” page for help with dosing the 96-well
925 plate. Sign and date the print out and store in study notebook.

926 **13.2 Visual Evaluation of Cell Viability**

- 927 1. 19 to 24 hours after dosing the plate, remove the plate from the incubator and
928 remove the media from the wells by inverting the plate onto blotter paper. Gently
929 tap plate against the bench surface to remove residual liquid trapped in the wells.

- 930 2. Use a repeat pipetter to add 50 µL 1X PBS to all wells. Immediately remove PBS
931 by inversion.
- 932 3. Using an inverted microscope, inspect all of the wells used in the 96-well plate
933 and record the visual observations using the scores in **Table 13-1**.

934 **Table 13-1 Visual Observation Scoring**

Viability Score	Brief Description¹
1	Normal Cell Morphology and Cell Density
2	Altered Cell Morphology and/or Small Gaps between Cells
3	Altered Cell Morphology and/or Large Gaps between Cells
4	Few (or no) Visible Cells
1P	Score of 1 with Precipitate
2P	Score of 2 with Precipitate
3P	Score of 3 with Precipitate
4P	Score of 4 with Precipitate
5P	Unable to View Cells Due to Precipitate

935 ¹Reference micrographs will be provided by NICEATM.

936

937 **13.3 Lysis of Cells for LUMI-CELL® ER**

- 938 1. Apply the reflective white backing tape to the bottom of the 96-well plate (this
939 will increase the effectiveness of the luminometer).
- 940 2. Add 30µL 1X lysis reagent to the assay wells and place the 96-well plate on an
941 orbital shaker for one minute.
- 942 3. Remove plate from shaker and measure luminescence (as described in **Section**
943 **11.5**).

944 **13.4 CellTiter-Glo® Assessment of Cell Viability**

945 When considered necessary, a quantitative evaluation of cell viability will be performed with the
946 Promega CellTiter-Glo® assay system. CellTiter-Glo® uses luminescence as an indicator of the
947 number of cells per plate and therefore must be conducted in parallel with the LUMI-CELL® ER
948 test method (i.e., both test methods cannot be conducted on the same plate).

- 949 1. Dose and incubate cells under the same conditions as for LUMI-CELL® ER.

- 950 2. Remove plates from incubator and discard the medium by inverting the plate onto
951 blotter paper. Gently tap plate against the bench surface to remove residual liquid
952 trapped in the wells.
- 953 3. Use a repeat pipetter to add 50 µL 1X PBS to all assay wells. Immediately
954 remove PBS by inversion.
- 955 4. Examine all wells used under an inverted microscope. Make notes of any well
956 with codes described in **Table 13-1**.
- 957 5. Place white backing tape on the bottom of the 96-well plate.
- 958 6. Add 100 µL estrogen-free DMEM to each well containing cells.
- 959 7. Add 100 µL CellTiter-Glo® reagent to each well containing cells.
- 960 8. Place plate on an orbital shaker for 1 minute to induce cell lysis.
- 961 9. Incubate (see **Section 9.0** for details) for 10 min.
- 962 10. Measure luminescence promptly. Do *not* add luciferase reagent to CellTiter-Glo®
963 plates.

964 **13.5 Measurement of Luminescence**

965 Luminescence is measured in the range of 300 to 650 nm, using an injecting luminometer and
966 with software that controls the injection volume and measurement interval. Light emission from
967 each well is expressed as relative light units (RLU) per well. The luminometer output is saved as
968 raw data in an Excel® spread sheet. A hard copy of the luminometer raw data should be signed,
969 dated and stored in the study notebook.

970 **13.6 Data Analysis**

971 LUMI-CELL® ER uses an Excel® spreadsheet to collect and adjust the RLU values obtained
972 from the luminometer and a GraphPad Prism® template to analyze and graph data. Plate
973 reduction is calculated using unadjusted RLU values.

974 The Excel® spreadsheet subtracts background luminescence (average DMSO solvent control
975 RLU value) from test substance, reference standard and control RLU values. Test substance,
976 reference standard, and control RLU values are then adjusted relative to the highest Ral/E2

977 reference standard RLU value, which is set to 10,000. After adjustment, values are transferred to
978 GraphPad Prism® for data analysis and graphing.

979 13.6.1 Correction and Adjustment of Luminometer Data

980 The following steps describe the procedures required to populate the Excel® spreadsheet that has
981 been configured to collect and adjust the RLU values obtained from the luminometer.

- 982 1. Open the raw data file and the corresponding experimental Excel® spreadsheet
983 from **Section 13.1.1**.
- 984 2. Copy the raw data using the Excel® copy function, then paste the copied data into
985 cell C22 of the “RAW DATA” tab in the experimental Excel® spreadsheet using
986 the Paste Special – Values command. This position corresponds to position B2 in
987 the table labeled Table 1 in this tab.
- 988 3. Examine the DMSO data in Table 1 of the Excel® spreadsheet to determine
989 whether there are any potential outliers. See **Section 13.6.2** for further explanation
990 of outlier determinations.
- 991 4. If an outlier is identified, perform the following steps to remove the outlier from
992 calculations:
 - 993 ▪ correct the equation used to calculate DMSO background in Table 1,
994 e.g., if outlier is located in cell G24, adjust the calculation in cell H42 to
995 read =AVERAGE(F24,H24,I24)
 - 996 ▪ then correct the equation used to calculate the average DMSO value in
997 Table 2, e.g., following the above example, adjust cell M44 to read
998 =AVERAGE(F36,H36,I36)
 - 999 ▪ then correct the equation used to calculate the standard deviation of the
1000 DMSO value in Table 2, e.g., following the above example, adjust cell
1001 M45 to read =STDEV(F36,H36,I36)
- 1002 5. Excel® will automatically subtract the background (the average DMSO control
1003 value) from all of the RLU values in Table 1 and populate Table 2 with these
1004 adjusted values.

- 1005 6. To calculate plate reduction, identify the cell containing the Ral/E2 replicate that
1006 has the lowest RLU value and the cell containing the RLU value for the same
1007 concentration in the corresponding Ral/E2 replicate (e.g., if the lowest RLU value
1008 for Ral/E2 is located in cell E23, the corresponding cell would be E22).
- 1009 7. Identify the cell containing the Ral/E2 replicate that has the highest RLU value
1010 and the cell containing the RLU value for the same concentration in the
1011 corresponding Ral/E2 replicate
- 1012 8. Click into cell D16 and enter the cell numbers from **Section 13.6.1 step 6** into the
1013 numerator and the cell numbers from **step 7** into the denominator.
- 1014 9. Click on the “ER Antagonist Report” tab.
- 1015 10. The data for the Ral/E2 reference standard, DMSO, E2, and Flavone/E2 replicates
1016 populate the left portion (columns A-F) of the spreadsheet. The data is
1017 automatically placed into an Excel® graph.
- 1018 11. To set the highest RLU value for the reference standard to 10,000 RLU, go to cell
1019 D2 of “ER Antagonist Report” tab and check the formula contained within that
1020 cell. The divisor should be the cell number of the cell containing the highest
1021 averaged Ral/E2 RLU value (column E).
- 1022 12. Copy the data into GraphPad Prism® for graphing and analysis as indicated in the
1023 NICEATM Prism® users guide.

1024 13.6.2 Determination of Outliers

1025 The Study Director will use good statistical judgment for determining “unusable” wells that will
1026 be excluded from the data analysis and will provide an explanation in the study notebook for any
1027 excluded data. This judgment for data acceptance will include Q-test analysis.

1028 The formula for the Q test is:

1029
$$\frac{\text{Outlier} - \text{Nearest Neighbor}}{\text{Range (Highest} - \text{Lowest)}}$$

1030 where the outlier is the value proposed for exclusion, the nearest neighbor is the value closest to
1031 the outlier, and the range is the range of the three values.

1032 If the value of this ratio is greater than 0.94 (the Q value for the 90% confidence interval for a
1033 sample size of three), the outlier may be excluded from data analysis. For Ral/E2 reference
1034 standard replicates (sample size of two), any adjusted RLU value for a replicate at a given
1035 concentration of Ral/E2 is considered and outlier if its value is more than 20% above or below
1036 the adjusted RLU value for that concentration in the historical database.

1037 13.6.3 Acceptance Criteria

1038 Acceptance or rejection of a test is based on evaluation of reference standard and control results
1039 from each experiment conducted on a 96-well plate. Results are compared to quality controls
1040 (QC) for these parameters derived from the historical database, which are summarized below.

- 1041 • Reduction: Plate reduction, as measured by dividing the averaged highest Ral/E2
1042 reference standard RLU value by the averaged lowest Ral/E2 reference standard
1043 RLU value, must be greater than three fold.
- 1044 • Reference standard results: Calculated Ral/E2 reference standard IC₅₀ values must
1045 be within 2.5 times the standard deviation of the historical database IC₅₀ mean
1046 value.
- 1047 • Solvent control results: DMSO control RLU values must be within 2.5 times the
1048 standard deviation of the historical solvent control mean RLU value.
- 1049 • E2 control results: E2 control RLU values must be within 2.5 times the standard
1050 deviation of the historical E2 control mean RLU value.
- 1051 • Positive control results: Flavone/E2 control RLU values must be within 2.5 times
1052 the standard deviation of the historical database flavone/E2 control mean RLU
1053 value.


1054 An experiment that fails any single acceptance criterion will be discarded and repeated.

1055 13.6.4 Calculation of Relative IC₅₀ Values

1056 Relative IC₅₀ values are calculated as directed in the NICEATM Prism® users guide. The
1057 concentration curve must have a sigmoidal shape and reach saturation at the highest and lowest
1058 concentrations tested (**Figure 13-1**) for Prism® to calculate a valid relative IC₅₀ value.

1059

1059 **Figure 13-1 Example Concentration Curve for Calculation of Relative IC₅₀ Value**



1060

1061 The mathematical model used by Prism® to calculate an IC₅₀ value uses the Hill function as
1062 described in Section 6.0.

1063 13.6.5 Calculation of Absolute IC₅₀ Values

1064 Calculate an absolute IC₅₀ (Eli Lilly 2005) for all substances that are positive for antagonism at
1065 one or more concentrations that cause a 50% reduction in response compared to the E2 control
1066 (Figure 13-2).