## 6.3.6. Chemical Exposure to Cells

- (1) Mix well the medium in the deep-well plate.

**Note: Confirm that the chemicals in DCC-FBS-EMEM have not precipitated. If precipitation is observed, it must be recorded.**

- (2) Add 50  $\mu$ L of the medium solution to both plates for antagonist and cytotoxicity assays prepared in "6.1." in triplicate.
- (3) Record if there is precipitation or not in wells by visually comparing the spike-in control wells in the cytotoxicity assay plate using a microscope.

(The plate for cytotoxicity assay has a clear bottom to allow this observation.)

- (4) Incubate for 20-24 hr in CO<sub>2</sub> incubator at 37°C.

**Note: The same medium dilution should be used to prepare the plate for cytotoxicity testing.**

|   | Chem.-A            |   |   | Chem.-B            |   |   | Chem.-C            |   |   | Chem.-D            |    |    |
|---|--------------------|---|---|--------------------|---|---|--------------------|---|---|--------------------|----|----|
|   | 1                  | 2 | 3 | 4                  | 5 | 6 | 7                  | 8 | 9 | 10                 | 11 | 12 |
| A | 10 <sup>-3</sup> M | → | → | 10 <sup>-3</sup> M | → | → | 10 <sup>-2</sup> M | → | → | 10 <sup>-3</sup> M | →  | →  |
| B | 10 <sup>-4</sup> M | → | → | 10 <sup>-4</sup> M | → | → | 10 <sup>-4</sup> M | → | → | 10 <sup>-4</sup> M | →  | →  |
| C | 10 <sup>-5</sup> M | → | → | 10 <sup>-5</sup> M | → | → | 10 <sup>-5</sup> M | → | → | 10 <sup>-5</sup> M | →  | →  |
| D | 10 <sup>-6</sup> M | → | → | 10 <sup>-6</sup> M | → | → | 10 <sup>-6</sup> M | → | → | 10 <sup>-6</sup> M | →  | →  |
| E | 10 <sup>-7</sup> M | → | → | 10 <sup>-7</sup> M | → | → | 10 <sup>-7</sup> M | → | → | 10 <sup>-7</sup> M | →  | →  |
| F | 10 <sup>-8</sup> M | → | → | 10 <sup>-8</sup> M | → | → | 10 <sup>-8</sup> M | → | → | 10 <sup>-8</sup> M | →  | →  |
| G | 0.1%<br>DMSO       | → | → | 0.1%<br>DMSO       | → | → | 1 $\mu$ M<br>OHT   | → | → | 100 $\mu$ M<br>Dig | →  | →  |
| H | 0.1%<br>DMSO       | → | → | 0.1%<br>DMSO       | → | → | 1 nM of<br>E2      | → | → | 1 nM of<br>E2      | →  | →  |

**Fig. 7 Assay plate layout for estrogenic assay**

Medium to be added: 10% DCC-FBS-EMEM containing 75 pM E2 (=  $7.5 \times 10^{-11}$ M), resulting in a final concentration of 25pM E2.

Medium to be added: 10% DCC-FBS-EMEM containing 0.3 % DMSO, resulting in a final DMSO concentration of 0.2%.

#### 6.4. LUCIFERASE ASSAY

After 20-24 hr incubation, luciferase activity can be measured with the luciferase assay reagent and a luminometer in accordance with the manufacturer's instructions.

If Steady-Glo Luciferase Assay System is used, the following procedure can be used to minimize the use of substrate solution.

- (1) Newly prepare or thaw the substrate solution
- (2) Mix 3 mL of the substrate solution and 3 mL of PBS (+) in a (for 1-plate use) in a reagent tray.
- (3) Remove all the exposure medium.
- (4) Add 50  $\mu$ L of substrate solution with PBS (+) using a multi-channel pipettor.
- (5) Leave the plate for 10-min at room temperature in dark area to protect it from light.
- (6) Read plates on a Chemiluminescence plate reader.

#### 6.5. CYTOTOXICITY (CELL VIABILITY) ASSAY USING CELL COUNTING KIT-8 (CCK-8)

After 20-24 hr incubation, cell viability can be measured with CCK-8 assay kit.

- (1) Mix 16.5 mL of 10% DCC-FBS-EMEM and 880  $\mu$ L of CCK-8 reagent in a reagent tray (for one plate) [CCK-8 solution].

*Note: CCK-8 solution should be prepared for each plate since this solution will increase in colour as time advances.*

- (2) Remove exposure medium.
- (3) Add 100  $\mu$ L of PBS (-) to each well and remove it to wash cells.
- (4) Add 150  $\mu$ L of CCK-8 solution to each well.
- (5) Incubate for 90-min in CO<sub>2</sub> incubator at 37°C.
- (6) Measure the absorbance at 450 nm using a microplate reader.

## 7. DATA ANALYSIS

### 7.1. DATA ANALYSIS FOR ESTROGENIC AND ANTI-ESTROGENIC ASSAYS

Using the spreadsheet provided, PC10 and PC50 for estrogenic assay and linIC30 and linIC50 for anti-estrogenic assay should be calculated. Additionally EC50 and IC50 for estrogenic and anti-estrogenic assays, respectively, should be calculated with the appropriate nonlinear curve-fitting software (GraphPad Prism (GraphPad Software Inc., San Diego, CA) is recommended.) using the following Hill's logistic equation;

$$Y = \text{Bottom} + (\text{Top} - \text{Bottom}) / (1 + 10^{((\text{LogE}(\text{IC}50 - X) * \text{HillSlope}))})$$

\* Where X is the logarithm of concentration, Y is the response and Y starts at the Bottom and goes to the Top with a sigmoid shape.


Fig. 8 Description of PC10, PC50 lin IC30 and IC50

### 7.2. CYTOTOXICITY DATA

In the cytotoxicity test, if the cell viability is less than 80%, the test chemical is cytotoxic at "that" test concentration and such data for evaluating anti-estrogenic assay should be omitted for the calculations of lin.IC30, lin.IC50 or IC50.

Cell viability should be calculated as below;

$$\text{Cell viability (\%)} = \frac{[(\text{Abs}_{.450} \text{ of test well}) - (\text{Mean of Abs}_{.450} \text{ of cytotox. control})]}{[(\text{Mean of Abs}_{.450} \text{ of VC}) - (\text{Mean of Abs}_{.450} \text{ of cytotox. control})]} \times 100$$

### 7.3. PERFORMANCE CRITERIA

In order to guarantee the assay performance, the performance criteria provided in this section should be fulfilled.

#### 7.3.1. Performance Criteria For Estrogenic Assay for Test System Setup (for Task-1)

Before entering task-2 testing, all laboratories should demonstrate the following performance criteria (Table 2) by using the qualified data that meets the quality criteria in

Table 3 to confirm whether the test system is correctly set up at each laboratory.

**Table 2 Performance Criteria for agonistic assay (Acceptable range of reference chemicals)**

|                    | Log[EC50]     | Log[PC50]     | Log[PC10]    | Hill slope |
|--------------------|---------------|---------------|--------------|------------|
| 17β-Estradiol (E2) | -11.3 ~ -10.1 | -11.4 ~ -10.1 | <-11         | 0.7 ~ 1.5  |
| 17α-Estradiol      | -9.6 ~ -8.4   | -9.6 ~ -8.1   | -10.7 ~ -9.3 | 0.9 ~ 2.0  |
| Corticosterone     | -             | -             | -            |            |

**Table 3 Quality Criteria for each plate for agonistic assay (control chemical)**

|  |   |
|--|---|
| Fold-induction <sup>*1</sup> of 1 nM of E2 | > = 4   |
| 10% fold-induction of 1 nM E2              | > 1 ± 2SD of fold-induction of vehicle control<br>(Mean luminescence intensity of 1 nM of E2) |
| *1: Fold-induction =                       | (Mean luminescence intensity of vehicle control)  |

### 7.3.2. Performance Criteria For Anti-estrogenic Assay (for Task-2 and Task-3)

Criteria for each plate with control chemicals (Table 4) and criteria for reference chemicals (Table 5) should be fulfilled.

**IMPORTANT NOTE: The final criteria may be modified during the study.**

**Table 4 Quality Criteria for each plate for anti-estrogenic assay**

|  |          |
|--|----------|
| Fold-induction of Spike-in Control (25 pM of E2) <sup>*1</sup> | >= 4     |
| RTA <sup>*2</sup> of 1 nM E2                                   | >= 100%  |
| RTA of 1 μM OHT  | =< 39.4% |
| RTA of 100 μM Dig.   | =< 0%    |

\*1: Fold-induction of Spike-in Control =  $\frac{(\text{Mean luminescence intensity of Spike-in Control})}{(\text{Mean luminescence intensity of vehicle control})}$

\*2: Relative Transcriptional Activation (RTA) =  $\frac{[(\text{luminescence intensity of a well}) - (\text{Mean luminescence intensity of vehicle control})]}{[(\text{Mean luminescence intensity of Spike-in Control}) - (\text{Mean of luminescence intensity of vehicle control})]}$

**Table 5 Performance Criteria for anti-estrogenic assay (Acceptable range for reference chemicals)**

|           | log [lin.IC30] | log [lin.IC50] | log [var.IC50] |
|-----------|----------------|----------------|----------------|
| OHT       | -9.62 ~ -8.73  | -9.46 ~ -8.16  | -9.32 ~ -8.20  |
| TAM       | -7.55 ~ -6.84  | -7.08 ~ -6.26  | -7.02 ~ -6.32  |
| RU486     | -6.18 ~ -5.41  | -5.61 ~ -5.08  | -5.53 ~ -4.86  |
| Flutamide | -              | -              | -              |

#: Reference chemicals should be tested in at least one run in each experiment (conducted under the same conditions including the materials, source of cells (passage), and technicians).

## Appendix-1

### Preparation of Serum treated with Dextran Coated Charcoal (DCC)

The treatment of serum with Dextran-coated charcoal (DCC) is a generally used methodology for the removal of estrogenic compounds from serum. It is added to the cell medium in order to exclude the biased response associated with residual estrogens in serum.

The following materials and equipments will be required;

#### Materials

- Activated charcoal (Sigma, Catalog# C9157)
- Dextran (MW 64,000~76,000, Sigma, Catalog# D4751)
- Magnesium chloride hexahydrate ( $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$ , Wako, Catalog# 135-00165,  $\geq 98\%$  or its equivalent)

Prepare 1 M  $\text{MgCl}_2$  aq. by dissolving 20.3 g of  $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$  in 100 mL of Milli-Q and filtering it with sterile filter.

- Sucrose (Wako, Catalog# 196-00015 or its equivalent)
- 1 M HEPES buffer solution (pH 7.4) (Gibco, Catalog# 15630)
- Ultrapure water produced from a filter system

#### Equipment

- Autoclaved glass container (size should be adjusted as appropriate)
- General Laboratory Centrifuge (that can set temperature at 4°C.)

The following procedure is adjusted for the use of 50 mL centrifuge tubes.

[Day-1] Prepare 1 litre of dextran coated charcoal suspension by adding the following reagents in the autoclaved glass container and stir it at 4°C, overnight.

- |                           |        |
|---------------------------|--------|
| • 1 M $\text{MgCl}_2$ aq. | 1.5 mL |
| • Sucrose                 | 85.5 g |
| • Activated charcoal      | 2.5 g  |
| • Dextrane                | 0.25 g |
| • 1 M HEPES               | 5 mL   |

[Day-2] Dispense the suspension in 50 mL centrifuge tubes and centrifuge at 10,000 rpm at 4°C for 10 minutes. Remove the supernatant and store half of the charcoal sediment at 4°C for the use on Day-3. Suspend the other half of the charcoal with fetal bovine serum (FBS) that is thawed at 42°C and left for 30 minutes at 56°C for heat inactivation, and transfer into the autoclaved glass container such as an Erlenmeyer flask. Stir this suspension gently at 4°C, overnight.

[Day-3] Dispense the suspension with FBS into centrifuge tubes for centrifugation at 10,000 rpm at 4°C for 10 minutes. Collect FBS and transfer into the new charcoal sediment prepared and stored on Day-2. Suspend the charcoal sediment and stir this suspension gently at 4°C, overnight.

[Day-4] Dispense the suspension for centrifugation at 10,000 rpm at 4°C for 10 minutes and sterilise the supernatant by filtration through 0.22 µm sterile filter. This DCC treated FBS should be stored at -20°C.

## HeLa ATG Validation Study

### Test Plate Quality Criteria and Assay Performance Criteria

#### Amendment:

Reason: The prior criteria were provisionally defined based solely upon the preliminary assay data from the lead laboratory (CERI), and it was expected that these might need minor modification following data analyses from the participating laboratories.

This amendment modifies the Quality Criteria and the Performance Criteria for the HeLa ATG assay based on the Task2 results from three Japanese laboratories. These modifications were agreed upon with the biostatistician at the data analysis meeting of HeLa ATG validation study held on 15 Dec 2008 at CERI.

#### Quality Criteria for each plate for anti-estrogenic assay:

The newly defined quality criteria and background calculated data are shown in Table 1. Modified sections are indicated in red characters.

# Detailed individual data were shown in Appendix. Table 1 ~4.

Table 1. Quality Criteria for each plate for anti-estrogenic assay

|  | Prior   | New            | mean±2SD      | mean±2.5SD    | mean±3SD     |
|--|---------|----------------|---------------|---------------|--------------|
| Fold-induction of Spike-in Control (25 pM of E2) | > 6     | ≥ 4            | 6.0 ~ 14.7    | 4.9 ~ 15.8    | 3.9 ~ 16.8   |
|  |         | (>6 recommend) | 7.5 ~ 12.8    | 6.8 ~ 13.4    | 6.2 ~ 14.1   |
| RTA of 1 nM E2                                   | > 100%  | ≥ 100%         | 114.9 ~ 204.5 | 103.7 ~ 215.7 | 92.5 ~ 226.9 |
| RTA of 1 μM OHT                                  | < 16.9% | ≤ 39.4%        | 0.5 ~ 39.4    | -4.4 ~ 44.3   | -9.3 ~ 49.1  |
|  |         |                | -0.2 ~ 38.7   | -5.0 ~ 43.6   | -9.9 ~ 48.4  |
| RTA of 100 μM Dig.                               | < 0%    | ≤ 0%           | -14.2 ~ -5.9  | -15.2 ~ -4.8  | -16.3 ~ -3.8 |
|  |         |                | -13.8 ~ -6.9  | -14.6 ~ -6.1  | -15.5 ~ -5.2 |

RTA : Relative transcriptional activation

#### Rationale for definition of Quality Criteria:

##### Fold induction:

While all data from three Japanese laboratories met the prior provisional criteria, two additional laboratories have so far had difficulty in achieving >6 fold induction (FI).

This suggests that the prior provisional criteria may be too stringent, and that for practical reasons the range could be extended a little, whilst also maintaining adequate performance, thereby accommodating all participating laboratories. It is expected that sufficient performance will be achieved if the FI is ≥4, based upon the experience of the lead laboratory with the agonist assay (the FI Criteria of the agonist assay is ≥4).

However, it is still recommended that the FI >6 remains as the preferred FI recommended value to obtain optimal results

RTA of 1  $\mu$ M OHT:

One laboratory (Kaneka) could not meet this criteria within five assay runs (see Appendix Table 3). However, their data was reproducible and concordant results (i.e. IC30, IC50) were observed for the reference chemicals. We therefore concluded that this was inter-laboratory variation. New criteria was defined by mean+2SD (39.4%) of all data (n=13).

RTA of 1 nM E2 and RTA of 100  $\mu$ M Dig.:

There was no change in prior criteria.

**Performance Criteria for anti-estrogenic assay (Acceptable range for reference chemicals):**

The newly defined criteria range for reference chemicals are shown in Table 2.

Modifications are indicated in red characters.

# Detailed individual data are shown in Appendix. Table 5 ~7.

Table.2 Performance Criteria for anti-estrogenic assay

|       |                | Prior         | New           | n* |
|-------|----------------|---------------|---------------|----|
| OHT   | log [lin.IC30] | -9.86 ~ -8.76 | -9.62 ~ -8.73 | 13 |
|       | log [lin.IC50] | -9.79 ~ -8.28 | -9.46 ~ -8.16 | 13 |
|       | log [var.IC50] | -9.15 ~ -8.94 | -9.32 ~ -8.20 | 13 |
| TAM   | log [lin.IC30] | -7.88 ~ -6.99 | -7.55 ~ -6.84 | 13 |
|       | log [lin.IC50] | -7.48 ~ -6.50 | -7.08 ~ -6.26 | 13 |
|       | log [var.IC50] | -7.17 ~ -6.77 | -7.02 ~ -6.32 | 13 |
| RU486 | log [lin.IC30] | -6.20 ~ -5.32 | -6.18 ~ -5.41 | 10 |
|       | log [lin.IC50] | -5.70 ~ -5.09 | -5.61 ~ -5.08 | 10 |
|       | log [var.IC50] | -6.22 ~ -5.32 | -5.53 ~ -4.86 | 10 |

n\*: Number of data used for calculation.

lin.IC30: 30% inhibition concentration by linear fitting.

lin.IC50: 50% inhibition concentration by linear fitting.

var.IC50: 50% inhibition concentration by non-linear fitting to Hill equation using Prism ver.4.

**Rationale for definition of Performance Criteria:**

The criteria ranges are defined based on mean $\pm$ 2SD as calculated from the results obtained from 13 assay runs during Task 2 from three Japanese laboratories (CERI: 5 assays, Otsuka: 3 assays, Kaneka: 5 assays). There were no large differences between the prior provisional range and the newly proposed provisional ranges on each parameter for each compound. For the parameter of RU486, three assay data sets (CERI: #5, Kaneka: #2, #5) were excluded from the calculation due to cytotoxicity, and at the recommendation of the biostatistician (see follows).

**Exclusion of data based upon chemical-induced cytotoxicity:**

The limit for cytotoxicity was set at <20% in the protocol for this STTA. Thus, no data were used for the calculation of the criteria parameters from any chemical concentration for which cell viability was less than 80% (e.g., cell cytotoxicity > 20%). In addition, non-linear curve fitting was not performed in cases where there were no data points were observed for more than 50% inhibition of ER-mediated gene transcriptional activation, after exclusion of the data points for cytotoxicity.

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表2-2 Quality Criteria for each plate for agonistic assay (control chemical)

|  |  |
|--|--|
| Fold-induction <sup>†1</sup> of 1 nM of E2 | >= 4   |
| 10% fold-induction of 1 nM E2              | > 1 ± 2SD of fold-induction of vehicle control |

$$^{\dagger 1}: \text{Fold-induction} = \frac{(\text{Mean luminescence intensity of 1 nM of E2})}{(\text{Mean luminescence intensity of vehicle control})}$$

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表4-1 Quality Criteria for each plate for anti-estrogenic assay

|  |         |
|--|---------|
| Fold-induction of Spike-in Control (25 pM of E2) <sup>*1</sup> | > 6     |
| RTA <sup>*2</sup> of 1 nM E2                                   | > 100%  |
| RTA of 1 μM OHT  | < 16.9% |
| RTA of 100 μM Dig.   | < 0%    |

$$*1: \text{Fold-induction of Spike-in Control} = \frac{(\text{Mean luminescence intensity of Spike-in Control})}{(\text{Mean luminescence intensity of vehicle control})}$$

$$[(\text{luminescence intensity of a well}) - (\text{Mean luminescence intensity of vehicle control})]$$

$$*2: \text{Relative Transcriptional Activation (RTA)} = \frac{[(\text{Mean luminescence intensity of Spike-in Control}) - (\text{Mean of luminescence intensity of vehicle control})]}{(\text{Mean luminescence intensity of vehicle control})}$$

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