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Mycobacteriology PT Procedual Manual (New York State Department of Health)

# Mycobacteriology PT Procedural Manual

## Table of Contents

Sample Set Description.....	p. 2
NYS Laboratory Permit Categories.....	p. 3
Lab Notification.....	p. 4
Organism Preparation.....	p. 4
Labeling and Tracking Organisms.....	p. 4
Transferring Specimens to 16 X 125 vials.....	p. 5
Cytospin Procedure.....	p. 5
Slide Staining procedure.....	p. 6
Transferring to Wheaton Bottles .....	p. 6
Contaminant Preparation.....	p. 7
Slide Preparation.....	p. 7
Protocol for Specimen Preparation .....	p. 9
Protocol for Isolate Preparation.....	p. 11
Shipping PT Samples.....	p. 12

## Sample Set Description

**A. Slides Categories: All Slides are numbered as follows:**

- # 0351 (first 2 digits are the year, third is sample type, last is sample #)
  - # 0352
  - # 0353
  - # 0354
  - # 0355
- 

**B. Specimens: Categories: All except Smears Only**

Specimens are numbered as follows:

- # SP 0361 (first 2 digits are the year, third is sample type, last is sample #)
  - # SP 0362
  - # SP 0363
  - # SP 0364
  - # SP 0365
- 

**C. Isolates Categories: Restricted, Restricted-S and Reference**

Isolates are numbered as follows:

- # IS 0371 (first 2 digits are the year, third is sample type, last is sample #)
- # IS 0372

**Categories: General, General-S, and Reference**

- # IS 0373 (first 2 digits are the year, third is sample type, last is sample #)
  - # IS 0374
- 

**D. Isolates Categories: Restricted-S, General S and Reference**

Isolates are numbered as follows:

- # IS 0381 (first 2 digits are the year, third is sample type, last is sample #)
- # IS 0382
- # IS 0383

## New York State Laboratory Permit Categories:

### MYCOBACTERIOLOGY

MYCOBACTERIOLOGY - Smears Only	(5 SLIDES)
MYCOBACTERIOLOGY - Restricted	(5 SLIDES, 5 SPECIMENS, 2 ISOLATES)
MYCOBACTERIOLOGY - Restricted-S	(5 SLIDES, 5 SPECIMENS, 2 ISOLATES, 3 ISOLATES)
MYCOBACTERIOLOGY - General	(5 SLIDES, 5 SPECIMENS, 2 ISOLATES)
MYCOBACTERIOLOGY - General-S	(5 SLIDES, 5 SPECIMENS, 2 ISOLATES, 3 ISOLATES)
MYCOBACTERIOLOGY - Reference	(5 SLIDES, 5 SPECIMENS, 6 ISOLATES, 3 ISOLATES)

---

**SMEARS ONLY.** This category is for laboratories that only examine smears for acid-fast bacilli. Laboratories holding this category must submit all specimens for growth detection and identification to a laboratory holding a New York State permit in the appropriate Mycobacteriology category.

**CATEGORY RECEIVES: 5 SLIDES**

---

**RESTRICTED and RESTRICTED-S** This category is for laboratories that examine smears for acid-fast bacilli, attempt to isolate acid-fast bacilli and identify *Mycobacterium tuberculosis*, *Mycobacterium bovis*, *Mycobacterium bovis* BCG, and *Mycobacterium africanum* as members of the *Mycobacterium tuberculosis* complex, and for those that perform *in vitro* susceptibility studies (**Restricted-S**)

**CATEGORY RECEIVES: 5 SLIDES**

**5 SIMULATED CLINICAL SPECIMENS  
2 ISOLATES FOR IDENTIFICATION**

**RESTRICTED-S-**

**3 ISOLATES FOR DRUG SUSCEPTIBILITY TESTING**

---

**GENERAL OR GENERAL-S.** These categories are for laboratories that conduct all clinical tests for detection, isolation and identification of all mycobacteria (**General**), and for those that perform *in vitro* susceptibility studies (**General-S**).

**CATEGORY RECEIVES: 5 SLIDES**

**5 SIMULATED CLINICAL SPECIMENS  
2 ISOLATES FOR IDENTIFICATION**

**GENERAL-S**

**3 ISOLATES FOR DRUG SUSCEPTIBILITY TESTING**

---

**REFERENCE:** These labs serve to validate test analytes. Reference Labs test all analytes.

**CATEGORY RECEIVES: 5 SLIDES**

**5 SIMULATED CLINICAL SPECIMENS  
4 ISOLATES FOR IDENTIFICATION  
3 ISOLATES FOR DRUG SUSCEPTIBILITY TESTING**

## Lab Notification

Mail out notification letter to all of the participating labs stating the date that they should expect to receive their samples and what they should do with those samples. Include a copy of updated mycobacteriology laboratory standards and any additional readings that max would like to send.

## Organism Preparation

### MATERIALS

- |                                    |                  |
|------------------------------------|------------------|
| ➤ ATCC strains                     | Chocolate Plates |
| ➤ ATCC list with freezer locations | 7H10 Plates      |
| ➤ Media order forms                | 3CC needles      |
| ➤ 12B vials                        | TB Data Base     |
| ➤ PAPR                             |                  |

- 1) Meet with Max to select the organisms for the event. Do not discuss the chosen specimens with others in the lab.
- 2) If possible use the ATCC list to find samples, and allow the frozen samples to thaw. If necessary use the database to find clinical specimens that meet the criteria of the event.
- 3) Following BSLIII protocols: Transfer the selected organisms into 12B vials. Depending on availability dispense .2ml into each 12B, .1ml onto each a chocolate and 7 H 10 plate.
- 4) If there is insufficient amounts of the original specimen inoculate only one 12 B, or if the organism is slow growing you should also inoculate only one 12B.
- 5) The ideal growing temperature of each organism should be considered in order to ensure proper growth. (Certain organisms grow better at alternate temperatures.)

## Labeling and Tracking Organisms

### MATERIALS

- |                                 |                         |
|---------------------------------|-------------------------|
| ➤ SATO thermal transfer printer | Thermal transfer labels |
| ➤ Label Maker program           | BACTEC                  |

- 1) Using the label maker program, label the 12B vials along with the chocolate, and 7H10 plates with the specimen identification number, last name, or ATCC strain, and the date.
- 2) Designate each vial containing the same specimen as either vial I or vial II.
- 3) Each specimen is to be tracked on a separate protocol sheet so that the GI and chocolate plate activity can be recorded daily.
- 4) **Select the needed tests for each individual organism (growth detection, identification, and susceptibility) on their respective protocol sheet so that the clinical lab pre-tests the samples to ensure that the correct organisms are prepared.**

- 5) Specimens should be placed in "Box two" and protocol sheets should be placed in the corresponding notebook so the growth indices can be monitored and recorded daily.
- 6) While the organisms are incubating, order 16x125mm vials containing Dubois Davis +tween 80 from the Media department. Be sure to mark "Sterile" on the order form. Order 6 tubes for every specimen used during the event)

## Transferring Organisms to 16X125 Vials

### MATERIALS

- 16x125mm vials containing Dubois Davis +tween 80 ordered from Media
- 1cc needles
- SATO thermal transfer printer      P APR
- Label Maker program                      7H10 Plates
- Thermal transfer labels                  Chocolate Plates
- Color dot labels                              37 Degree incubator
- Glas-Col rotator

- 1) For each time the organism is needed during the PT event prepare 5 16X125 vials containing Dubois Davis + tween 80 to account for any contamination that may occur. For each specimen choose the 12B vial that has grown to 999, or if both have grown to 999 choose the one that appears clearest, or the one that has grown closest to a doubling rate.
- 2) Each top of the 16X125 vials should be color-coded. The tube itself should be labeled as to which 12B vial the contents were taken, the specimen number, and the date the tube was inoculated.
- 3) The tubes should also be labeled sequentially
- 4) TB slow growing mycobacterium and slow growing non-TB mycobacterium require inoculation of 1ml. For all others types of bacteria, both TB and non-TB, you should inoculate with .5 ml. The 16x125mm vials are then placed in the Glas-Col rotator in the 37-degree incubator and allowed to grow until McFarland Standard # 2.
- 5) Chocolate plates should be done at the time of transfer to ensure that there was no contamination.
- 6) After two days of incubation a cytospin should be done on a 12x125 vial to ensure that there was no contamination during the transfer. Once McFarland two is suspected a final cyto spin can be done to check for purity. (Please see Cyto Spin Procedure below)
- 7) Order from the media department 20, 250ml Wheaton bottles with 200 ml of Dubois Davis + tween 80. Be sure to mark "sterile" on the order form.

## Cytospin Procedure

### MATERIALS

- Clear Frost slides                          Centrifuge
- Cyto Funnels                                  3CC Needles
- Cyto Funnel Caps                              P APR
- Cyto Spin Clamps                              Bovine Protein
- Cytospin Disk                                  Pencil

- 1) The GI of the organism will determine the amount of the organism that is needed to properly transfer bacterial colonies to the slide. The higher the GI the lower the amount of the organism needed, no more than .5 and no less than .2
- 2) Using a pencil, label clear frosted slides with the information that appears on the 12B vial so that there is no chance for cross contamination.
- 3) Set slides along with cyto funnels into slide holders. Be sure both are completely in place and clamped down. It is important to make sure that the slide and the funnel are locked tightly into place. When installing the funnels and the funnel caps it is crucial that the cap and the opening to the funnel are not contaminated.
- 4) Under BSLIII conditions transfer .2 ml of a Bovine protein into the funnel, this will help adhere the bacteria to the slide.
- 5) Maintaining BSL III conditions, extract between .2 and .5ml of the desired organism from the 12B vial and place in the funnel.
- 6) Place funnels into centrifuge disk. Make sure that there is a balanced load of specimens within the disk; it may be necessary to use an empty funnel and slide set.
- 7) Air seal and bring to centrifuge.

### Slide Staining Procedure

#### MATERIALS

- |                              |              |
|------------------------------|--------------|
| ➤ Hot Rack TB                | Methyl Blue  |
| ➤ Ziehl Neelsen              | Acid Alcohol |
| ➤ Distilled H <sub>2</sub> O | Rinse Rack   |

- 1) Place slides on hot rack.
- 2) Cover slides with Ziehl Neelsen stain for five minutes
- 3) Place Slides on rinse racks and flood with H<sub>2</sub>O.
- 4) Decolorize completely with acid alcohol usually 1min 30 seconds.
- 5) Flood Slide with H<sub>2</sub>O
- 6) Cover Slides with TB Methyl Blue for 45 Seconds.
- 7) Rinse with H<sub>2</sub>O. Let air dry and read.

### Transfer to Wheaton Bottles

#### MATERIALS

- 250 ml Wheaton Bottles containing 230 ml Dubois Davis
- Electric pipette dispenser
- 2 ml pipette tubes
- Chocolate plates
- PAPR
- 7H10 plates
- 37 degree incubator

- 1) Once McFarland 2 standard is assumed in the 16X125 vials they must be transferred into 250 ml Wheaton Bottles containing 230 ML of Dubois Davis+ tween 80.

- 2) Under BSL III standards use an electric pipette and 2ml pipette tubes, transfer the entire contents (10ml) of the 16X125 tube into the 250-ml Wheaton bottle containing 230 ml Dubois Davis and a sterile stir bar. A chocolate and 7H10 plate should be made at this time.
- 2) Each Wheaton bottle is to be labeled with the specimen number and also the tube from which it was inoculated.
- 3) The Wheaton bottles are to be placed on magnetic stir plates to ensure equal consistency within the bottle.
- 4) A cyto spin should be done after two days of incubation to ensure the presence of AFB and to rule out contamination.
- 5) When McFarland two is assumed smears can be made to check for slide concentration. If a numerous slide concentration is present slides can be prepared for the event. (See slide Protocol)
- 6) Order (15) 125-ml Wheaton bottles with 90 ml mucin from media department.
- 7) Order (8) ONE-liter bottles with 700 ml Mucin from media department.
- 8) Order (10) 16x125mm tubes with 10-ml tripticase soy broth for contaminate preparation

## Contaminant Preparation

### MATERIALS

- 10 – 16x125 mm tubes with tripticase Soy Broth
  - Three different contaminants from freezer
  - 1CC needles
  - Chocolate plates
1. Remove 3 of the 10 16x125mm tubes of tripticase soy broth (TSB) from the cooler. Along with the contaminants.
  2. Under BSLIII conditions add .2ml of one contaminant into each TSB tube.
  3. Remove the contaminates from the incubator at the 24 hour mark and place in the cooler so that McFarland 2 is not surpassed.
  4. Add 1ml of McFarland #1 contaminants to 90 ml H<sub>2</sub>O for proper dilution of contaminants.
  5. 1ml of properly diluted contaminants into specimen.

## Slide Preparation


### MATERIALS

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>➤ Color-frost slides (previously labeled)</li> <li>➤ Samples grown to McFarland 2 in 250 ml Wheaton's</li> <li>➤ Electric pipette</li> <li>➤ Pipette tubes 10-ml</li> <li>➤ 125-ml Wheaton bottles with 90 ml Mucin in each</li> <li>➤ Stir Bars</li> </ul> | <ul style="list-style-type: none"> <li>Forceps</li> <li>Bunsen burner</li> <li>Pipette tips</li> <li>Repeating pipette</li> <li>Slide trays</li> <li>Stir plates</li> </ul> |
|--|---|

- 1) Once a numerous concentration is reached within the 250 ml Wheaton bottle, and there is no Contamination in the sample, slides can be made.
- 2) Transfer enough of the sample from the McFarland 2, 250-ml Wheaton bottles into a 125 ml Wheaton with 90 ml mucin.

- 3) The procedure should start by transferring 10 ml into the mucin and should go up by 5 ml at a time until a numerous count is reached within the 125-ml mucin bottles.
- 4) The 125 ml Wheaton containing the mucin and a sterile stir bar should be placed on a stir plate and allowed to spin for a few minutes before preparing sample slides, to ensure a homogenous sample.
- 5) Once the numerous bottle is finished it should be labeled as numerous, and the specimen number, name, and the date the bottle was created.
- 6) Use the numerous bottle to create the remainder of your slides.
- 7) To create a moderate sample transfer the same amount of the numerous specimen into a 125-ml mucin as you did from the 250 ml Wheaton to the numerous.
- 8) To create a few sample transfer the same amount of the moderate sample into a 125 ml Mucin bottle as you did from the numerous to the moderate.
- 9) To create a rare sample transfer the same amount of the few sample into a 125ml mucin bottle as you did from the moderate to few.
- 10) Negative slides get contaminants ONLY.
- 11) Negative slides should be made first.
- 12) Set pre-labeled slides onto slide trays and arrange under the hood. It is important that the slides sit flat so that the Mycobacterium or contaminants do not run to one side of the slide and an even consistency remains in the center of the slide.
- 13) Using the rapid dispenser place 20 micro liters onto the slide and spread with the tip to a nickel sized spot.
- 14) Allow all of the slides to dry under the hood, and then heat fix.

After dilutions are made, organisms are dispensed on slides, air dried, heat fixed, and packed

## Protocol for Specimen Preparation


### MATERIALS

- Culture of AFB grown on rotary mixer to at least McFarland Standard No. 2.
- Mucin
- Contaminants
- Vacuum Pipette and Tips
- Vortex & magnetic stirrer
- Pre-labeled Nunc tubes
- Para film cut into 1 inch squares

### Specimen Preparation

1. Use selected organism that has grown to McFarland #2 in 200 ml of Dubois Davis. This is the undiluted specimen
2. **NUMEROUS** = Pipette needed amount (X) to of the McFarland #2 Specimen into the 1L Wheaton bottle containing 700 ml mucin so that upon smear a numerous count is seen.
3. **MODERATE** = Pipette X amount from the Numerous Wheaton Bottle into a different 1L wheaton bottle with 700ml mucin.
4. **FEW** = Pipette X amount of the moderate suspension into 1L Wheaton Bottle containing 700 ml mucin.
5. **RARE** = Pipette X amount of the few suspension into 1L Wheaton Bottle containing 700 ml mucin.
6. Negative specimens get only the contaminants.
7. Before Dispensing allow Dispensing Flasks to mix on magnetic stirrer for ~30 minutes.
8. Using 5 ml pipette tips and vacuum pipette aid dispense 3.0 ml into each labeled and color coded nunc tube.
9. Dispense Negative Specimens first
10. When completed, Para film the nunc tubes and store Specimens in the walk in cooler.

# Specimen Protocol


## Protocol for Isolate Preparation

### MATERIALS

- 250-ml Wheaton bottles with samples grown to McFarland 2
- Electric Pipette
- 5ml pipette tubes
- Chocolate Plates
- 7 H 10 plates
- Para film cut into 1inch squares
- Pre-labeled Nunc tubes
- Stir Plate
- Stir Bars

1. Once McFarland number 2 is achieved in the 250 ml Wheaton the samples that are to be used for growth detection, and Identification can be dispensed.
2. Use an electric pipette and 5ml pipette tubes to transfer 3 ml from the McFarland 2 bottle to individual pre-labeled nun tubes. A chocolate plate should be made at this
3. The top of each nunc tube should be wrapped with Para film to ensure closure and stored in the cold room until shipment.

## Shipping PT Samples

### Smears Only Laboratory

1. Into a cushioned mailer, place one (1) slide sample set.
2. Place paperwork into 7X10 envelope.
3. Deposit cushioned envelope and paperwork envelope into FedEx shipping envelope.
4. Affix corresponding FedEx mailing label to shipping envelope.

### Other Laboratories

Procedure:

1. Organize samples: retrieve the exact numbers needed for each category.

<b>Restricted</b>	<b>Restricted-S</b>	<b>Reference</b>	<b>General</b>	<b>General-S</b>
Slides	Slides	Slides	Slides	Slides
Isolate 1 & 2	ID Isolates 1 & 2	ID Isolates 1 - 4	ID Isolates 3 & 4	ID Isolate 3 & 4
5 Specimens	5 Specimens	5 Specimens	5 Specimens	5 Specimens
	Susceptibility Isolates 1,2 & 3	Susceptibility Isolates 1,2 & 3		Susceptibility Isolates 1,2 & 3

2. Retrieve Fed-EX mailing labels from LABLIC program and place on disk. Bring disk to mail room to print out dangerous goods labels and mailing labels.
3. Place canisters and appropriate paper work inside PT mailer box. Affix all necessary labels to box and seal.
4. Ship samples via FedEx overnight mail.

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