Realtime RT-PCR assays for the detection of pandemic (H1N1) 2009 virus, seasonal influenza viruses or avian (only group protocol No.2) influenza A viruses

Group protocol No.1²⁰

This protocol describes realtime RT-PCR procedures for the detection of:

1. Influenza type A viruses (M gene)
2. Pandemic (H1N1) 2009 virus (H1 gene)
3. Seasonal influenza A (H1N1) and A (H3N2) viruses (H1h, N1h, H3h, N2h genes)

This protocol is to be used in case of suspicion of human infection by the pandemic (H1N1) 2009 virus.

The following testing strategy is recommended:

- RNA extraction.
- Amplification in parallel of M, H1 pandemic and GAPDH (to assess quality of the specimen and extraction procedure) genes.
- In a separate set, amplification of the H1h, N1h, H3h and N2h genes.
- All the above genes are tested in a single plex and not in a multiplex format.

Materials required

- QIAamp Viral RNA (QIAGEN mini Kit 50) (Qiagen®, Cat. No. 52904)
- SuperScript™ III Platinum® One-Step qRT-PCR System (Invitrogen, Cat. No. 11732-020)
- Superscript™ III Platinum® One-Step qRT-PCR system (Invitrogen, Cat. No. 11732-088)
- Non acetyled BSA 10% (Invitrogen, Cat. No. P2046)
- LightCycler 1.5 or 2.0 (Capillaries)
- LightCycler 480
- 7500 Real-Time PCR System, Applied® Biosystems
- Smart Cycler® Cepheid

Nucleic acid extraction

RNA is extracted from specimens using the QIAamp Viral RNA kit (QIAGEN® Mini Kit 50 ref 52904). RNA extracted from 200 µl of the original sample is eluted in 60 µl of elution buffer.

All primers and probes described below were validated under the following conditions using the above equipment and reagents mentioned in the materials required section above.

Adjustments may be required for the use of other kits or other real-time PCR instruments.

Primers and probes for the detection of influenza A viruses (M gene), GAPDH and the influenza A (H1N1) pandemic 2009 virus (H1swl gene) were also validated under the following conditions:

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- 7500 Applied® Biosystems
- SmartCycler® Cepheid

Primers and probes

If the sample is positive for M and negative for H1sw, we suggest use of the H1N1h and H3N2h sets of primers.

Type/subtype	Gene	Name	Sequences	Bases	PCR Product size	Reference
Seasonal H1N1	HA	H1h-678Fw	CACCCAGAAATAGCCAAAA	20		1
Seasonal H1N1	HA	H1h-840Rv	TCCTGATCCAAAGCCTCTAC	20	163 bp	1
Seasonal H1N1	HA	H1h-715probe	5'-Fam-CAGGAAGGAAGAATCAACTA 3'-BHQ-1	20		1
Seasonal H3N2	HA	H3h-177Fw	GAGCTGGTTCAGAGTTCCTC	20		1
Seasonal H3N2	HA	H3h-388Rv	GTGACCTAAGGGAGGCATAATC	22	211 bp	1
Seasonal H3N2	HA	H3h-306Probe	5'-Fam-TTTTGTGAACGCAGCAAAG 3'-BHQ-1	20		1
Seasonal H3N2	NA	N2h-1150 Fw	GTCCAMACCTAAYTCAA	18		1
Seasonal H3N2	NA	N2h-1344 Rv	GCCACAAAACACAACAATAC	20	194 bp	1
Seasonal H3N2	NA	N2h-1290 probe	5'-Fam-CTTCCCCTTATCAACTCCACA-3'-BHQ-1	21		1
Seasonal H1N1	NA	N1h-1134 Fw	TGGATGGACAGATACCGACA	20		1
Seasonal H1N1	NA	N1h-1275 Rv	CTCAACCCAGAAGCAAGGTC	20	142 bp	1
Seasonal H1N1	NA	N1h-1206 probe	5'-Fam-CAGCGGAAGTTTCGTTCAACAT 3'-BHQ-1	22		1
Pandemic H1N1	HA	GRswH1-349Fw	GAGCTAAGAGAGCAATTGA	19		1
Pandemic H1N1	HA	GRswH1-601Rv	GTAGATGGATGGTGAATG	18	253bp	1
Pandemic H1N1	HA	GRswH1-538Probe(-)	Fam 5'-TTGCTGAGCTTTGGGTATGA -3'-BHQ-1	20		1
Influenza type A		GRAM/7Fw	CTTCTAACCGAGGTCGAAACGTA	23		2
Influenza type A		GRAM/161Rv	GGTGACAGGATTGGCTTGCTTTA	25	202 bp	2
Influenza type A		GRAM probe/52/+	5'-Fam-TCAGGCCCTCAAAGCCGAG 3'-BHQ-1	21		2
Human internal control	GAPDH	GAPDH-6Fw	GAAGGTGAAGTCCGGAGT	18		3
Human internal control	GAPDH	GAPDH-231Rv	GAAGATGGTGATGGGATTC	20	226 bp	3
Human internal control	GAPDH	GAPDH-202Probe(-)	5'-Fam-CAAGCTTCCCGTTCTCAGCC-3'-BHQ-1	20		3

1. National Influenza Center (Northern-France), Institut Pasteur, Paris.

2. Wong et al., 2005, J. Clin. Pathol. 58;276-280.

3. National Influenza Center (Southern-France), CHU, Lyon.

Procedure

LightCycler 1.5 or 2.0 (Capillaries) (ROCHE)

Reagent	Volume (μ l)	Final concentration
H ₂ O PPI	1.06	
Reaction mix 2x	10.0	3.0 mM Mg
MgSO ₄ (50mM)	0.24	0.6 mM Mg
Forward Primer (10 μ M)	1.0	0.5 μ M
Reverse Primer (10 μ M)	1.0	0.5 μ M
Probe (10 μ M)	0.4	0.2 μ M
BSA non acetylated (10mg/ml)	0.5	0.25 mg/ml
Superscript III RT/Platinum Taq Mix	0.8	
Total	15.0	

LightCycler 480 (96-well format) (ROCHE)

Reagent	Volume (μ l)	Final concentration
H ₂ O PPI	1.56	
Reaction mix 2x	10.0	3.0 mM Mg
MgSO ₄ (50mM)	0.24	0.6 mM Mg
Forward Primer (10 μ M)	1.0	0.5 μ M
Reverse Primer (10 μ M)	1.0	0.5 μ M
Probe (10 μ M)	0.4	0.2 μ M
Superscript III RT/Platinum Taq Mix	0.8	
Total	15.0	

15 μ l of reaction mix + 5 μ l of RNA samples.

7500 Real-Time PCR System (Applied Biosystems) or SmartCycler (Cepheid) Real-Time PCR machines

Reagent	Volume (μ l)	Final concentration
H ₂ O PPI	1.5	
Reaction mix 2x	12.5	3.0 mM Mg
Forward Primer (10 μ M)	2.0	0.8 μ M
Reverse Primer (10 μ M)	2.0	0.8 μ M
Probe (10 μ M)	1.0	0.2 μ M
ROX reference dye (diluted 1/10)	0.5	
Superscript III RT/Platinum Taq Mix	0.5	
Total	20.0	

20 μ l of reaction mix + 5 μ l of RNA samples

Controls

Each real-time RT-PCR assay includes additional unknown samples:

- Two negative samples bracketing unknown samples during RNA extraction (negative extraction controls);
- Positive controls (in duplicate); when using *in vitro* synthesized transcripts as controls include five quantification positive controls (in duplicate) including 10^4 , 10^3 and 10^2 copies of *in vitro* synthesized RNA transcripts; and
- One negative amplification control.

LightCycler System

Amplification cycles:	Temperature (°C)	Time (minute:second)	No. of cycles
Reverse transcription	45	15:00	1
Denaturation	95	3:00	1
Amplification	95	0:10	} 50
	55	0:10	
	72	0:20	
Cooling	40	0:30	1

7500 Applied or Smartcycler System

Amplification cycles:	Temperature (°C)	Time (minute:second)	No. of cycles
Reverse transcription	50	2:00	1
Denaturation	95	15:00	
Amplification	95	0:15	50
	60	0:40	

Sensitivity

For the M real-time RT-PCR

Sensitivity is about 100 copies of RNA genome equivalent per reaction (95% confidence level). This amount of target sequences is always detected, however, the probability to detect lower amounts of virus decreases accordingly. In our settings, samples containing 10 copies could be detected.

For the H1 pandemic real-time RT-PCR

Sensitivity is comparable to that of the M real-time RT-PCR and comparable to the sensitivity of the CDC kit (Cp <36 for all positive specimens tested so far).

For the H3h and N2h real-time RT-PCR

Sensitivity of the H3h real-time RT-PCR is equivalent to that of the M real-time PCR (Cp H3h ≈ Cp M) but the sensitivity of the N2h real-time RT-PCR is lower (Cp N2h ≈ Cp M + 5 Cp)

For the H1h and N1h real-time RT-PCR

Sensitivity of the N1h real-time RT-PCR is equivalent to that of the M real-time PCR (Cp N1h ≈ Cp M) but the sensitivity of the H1h real-time RT-PCR is lower (Cp H1h ≈ Cp M + 4 Cp).

Specificity

For the H1 pandemic real-time RT-PCR

So far, limited testing showed no detection for seasonal influenza viruses (influenza A(H1N1), A(H3N2), B) nor for specimens known to be positive for other respiratory viruses (influenza C, RSV A, B, hBoV, hPIV1,3, hMPV, HRV, enterovirus, adenovirus, CMV, HSV, VZV).

For swine influenza viruses, detection was positive for A/sw/England/117316/86 (classical swine lineage) and negative for A/sw/England/502321/94 (H3N2).

For pandemic (H1N1) 2009 viruses, detection was positive for A/California/4/2009, as well as for more than 10 specimens positive for the novel A(H1N1) pandemic virus.

NOTE: The H1 pandemic real-time RT-PCR does not detect the positive control from the CDC kit.

Positive control for M and GAPDH real-time RT-PCR

Positive control for M real-time RT-PCR is an *in vitro* transcribed RNA derived from strain A/Paris 650/06(H1N1). The transcript contains the Open Reading Frame of the M gene (from the ATG to nt 982) as negative standard. Each microtube contains 10^{11} copies of target sequences diluted in yeast tRNA and lyophilised.

Positive control for GAPDH real-time RT-PCR is an *in vitro* transcribed RNA. The transcript contains the Open Reading Frame of the M gene (from nt 6 (ATG = 1) to nt 231) as a negative strand. Each microtube contains 10^{11} copies of target sequences diluted in yeast tRNA and lyophilised.

Reconstitution of transcribed RNA

- Add 100 μ l of distilled water to obtain a solution at a concentration of 10^9 copies/ μ l. Store at -80°C .
- Dilute in H_2O to prepare a master bank at 2×10^6 copies/ μ l. Store at -80°C .
- From this, prepare a working bank of reagent at 2×10^4 copies/ μ l in order to avoid freeze/thaw cycles. Working tubes may be stored at -20°C for less than one week.
- Positive controls are available upon request at grippe@pasteur.fr.

Interpretation of results

GAPDH reactions should give a $C_p < 35$; if higher and otherwise negative results are obtained this may result from:

- Poor quality of the specimen with insufficient number of cells; obtain a new specimen for the same patient.
- Presence of inhibitors; repeat the procedure with dilutions of the extracted RNA (e.g. 1:10, 1:100) and/or repeat RNA extraction.

Positive reactions for M and H1 pandemic and negative reactions for H1h, N1h, H3h, N2h: confirmed case for pandemic (H1N1) 2009 virus.

Positive reaction for M and negative for H1 pandemic and for H1h, N1h, H3h, N2h (usually seen for low virus load in specimen), repeat reactions and/or repeat RNA extraction.

Positive reaction for M and either N1h or H3h but negative for H1 pandemic and for H1h and N2h (usually seen for low virus load in specimen) may indicate an infection with seasonal virus; in this case, repeat reactions and/or repeat RNA extraction to determine subtypes.

Positive reaction for M and H1 pandemic virus and positive reaction for either N1h, H3h may reflect a cross-contamination or a possible co-infection with both the pandemic (H1N1) 2009 virus and a seasonal virus; repeat RNA extraction and repeat reactions with all necessary precautions to avoid cross-contamination.

Group protocol No. 2²¹

This protocol describes realtime RT-PCR procedures for the detection of:

1. Influenza type A viruses (M gene)
2. Pandemic influenza A (H1N1) 2009 virus (N1 gene)
3. Pandemic influenza A (H1N1) 2009 virus (HA gene)
4. Seasonal influenza A (H1N1) and A (H3N2) viruses (H1h, H3h, genes)
5. Avian influenza A (H5) viruses, including highly pathogenic H5N1 viruses
6. Avian influenza A (H7) viruses
7. Influenza type B viruses

Materials required

1. ABI7500 machine
2. TaqMan EZ RT-PCR Core Reagents kit (Applied Biosystems, Cat. No. N8080236)
3. Amplification and detection is performed on ABI7500 machine with the TaqMan EZ RT-PCR Core Reagents kit (Applied Biosystems), Using 20 µl of RNA extract in an end volume of 50 µl

Test validation

All tests have been optimized and validated to work in the platform described above. They have not been optimized for other platforms and other reagents. For instance, some of these tests do not work well with other reagents (e.g. Ampli-Taq Gold).

The influenza B virus test, the influenza A virus M-gene test, the human influenza A/H1 virus test, and the human influenza A/H3 virus test have been used successfully for several years in a diagnostic virology laboratory with human specimens. The influenza A virus M-gene test, the influenza A/H5 virus tests and the influenza A/H7 test have been used successfully for several years in a diagnostic virology laboratory with human and animal specimens.

The seasonal influenza A/H1 test does not detect the pandemic H1N1 2009 virus H1 sequences, only seasonal influenza A (H1N1) virus H1 sequences. The pandemic H1N1 2009 NA and HA tests do not detect seasonal influenza A (H1N1) NA or HA sequences, only N1 and H1 of the pandemic H1N1 2009 virus sequences. This has been validated using sets of numerous recent human H1N1 viruses and swine and the pandemic H1N1 viruses. In addition, the tests have proven successful in detecting sporadic cases of the pandemic H1N1 2009 viruses in The Netherlands.

²¹ Dept Virology, Erasmus MC Rotterdam, Netherlands. (www.virology.nl)

Reaction mixture (EZ RT-PCR)

Reagent	Volume (µl)
DEPC water	4.5 µl
5x Taqman EZ buffer	10.0 µl
Manganese acetate (25 mM)	6.0 µl
dNTP's	6.0 µl
Primer/probe mix (1 µl/reaction)	1.0 µl
rTth DNA Polymerase (2,5U/µl)	2.0 µl
Amperase UNG (1U/µl)	0.5 µl
Total mix volume	30.0
Input RNA	20.0 µl
Total volume	50.0

Thermal Cycling conditions

Temperature (°C)	Time (minute:second)	No. of cycles
50	02:00	1
60	30:00	1
95	05:00	1
94	00:20	40[a]
62	01:00	

[a] For pandemic (H1N1) 2009 NA and HA gene specific test, 45 cycles are used.

Primers and probes

Type/subtype	Primer/probe ²²	Gene fragment	Oligonucleotide sequence	Quantity/reaction (pmol)
Flu B	Primer 1	M	GAGACACAATTGCCTACCTGCTT	10
Flu B	Primer 2	M	TTCTTCCCACCGAACCAAC	40
Flu B	Probe	M	AGAAGATGGAGAAGGCAAAGCAGAAGCTAGC	5
Flu type A	Primer 1	M	AAGACCAATCCTGTACCTCTGA	40
Flu type A	Primer 2	M	CAAAGCGTCTACGCTGCAGTCC	30
Flu type A	Probe	M	TTGTGTTCACGCTCACCGTGCC	5
Flu A (H1N1) seasonal	Primer 1	HA	GAATAGCCCCACTACAATTGGGTAA	10
Flu A (H1N1) seasonal	Primer 2	HA	GTAATTCGCATTCTGGGTTTCCT	10
Flu A (H1N1) seasonal	Probe	HA	AAGATCCATCCGGCAACGCTGCA	10
Flu A (H3) seasonal	Primer 1	HA	GATGTGTACAGAGATGAAGCATTAAACA	40
Flu A (H3) seasonal	Primer 2	HA	TAGGATCCAATCTTTGTATCCTGACTT	15
Flu A (H3) seasonal	Probe 1	HA	AGCTCAACACCTTTGATCTGGAACCGG	10
Flu A (H3) seasonal	Probe 2	HA	AGCTCAACGCCTTTGATCTGGAACCGG	10
Flu A (H5N1) HPAI, Set 1	Primer 1	HA	GAGAGGAAATAAGTGGAGTAAATGGGA	30
Flu A (H5N1) HPAI, Set 1	Primer 2	HA	AAGATAGACCAGCTACCATGATTGC	2.5
Flu A (H5N1) HPAI, Set 1	Probe	HA	TTTATTCAACAGTGCGGAGTCCCTAGCACT	10
Flu A (H5N1) HPAI, Set 2	Primer 1	HA	GGAACCTACCAAATACTGTCAATTTATTCA	40
Flu A (H5N1) HPAI, Set 2	Primer 2	HA	CCATAAAGATAGACCAGCTACCATGA	2.5
Flu A (H5N1) HPAI, Set 2	Probe	HA	TTGCCAGTGCTAGGGAACCTGCCAC	10
Flu A (H7) AI viruses	Primer 1	HA	GGCAACAGGAATGAAGAATGTTCC	40
Flu A (H7) AI viruses	Primer 2	HA	AATCAGACCTTCCCATCCATTTTC	15
Flu A (H7) AI viruses	Probe	HA	AGAGGCCTATTTGGTGCTATAGCGGGTTTCA	10
Pandemic 2009 virus*	Primer 1	HA	GGAAAGAAATGCTGGATCTGGTA	30
Pandemic 2009 virus*	Primer 2	HA	ATGGGAGGCTGGTGTATTATAGC	30
Pandemic 2009 virus*	Probe	HA	TGCAATACAACCTGTCAGACACCCAAGGG	10
Pandemic 2009 virus	Primer 1	NA	ACATGTGTGTGCAGGGATAACTG	30
Pandemic 2009 virus	Primer 2	NA	TCCGAAAATCCCACTGCATAT	30
Pandemic 2009 virus	Probe	NA	ATCGACCGTGGGTGCTTTCAACCA	10


* A new protocol which did not exist in the 18 August 2009 version of this document.

Test result interpretation


M gene pos, seasonal H1 pos, pandemic 2009 H1 and N1 neg: Human seasonal H1 influenza

M gene pos, seasonal H1 neg, pandemic 2009 H1 or N1 pos: pandemic (H1N1) 2009 virus

Annex 3:
PCR testing algorithm and results interpretation

In the event that consent for a full autopsy is not forthcoming, the attending physicians should consider a para-mortem needle biopsy approach. In Hong Kong, during the SARS outbreak in 2003, this was used in 24 of the 39 post-mortem cases. These biopsies can be obtained percutaneously and apart from minor puncture wounds there is no disfigurement to the body and infectious risk to autopsy personnel is low.

Trucut biopsy needle that can be used for needle sampling of organs.

LUNG

The area of the lung chosen for Trucut ideally should correspond to radiological signs of pneumonia. For negative control, a biopsy from a radiological normal area is also advised.

A minimum of four cores of tissue should be sampled, five if liquid nitrogen is available. It should be noted that samples involving fixation should be performed last so as to avoid potential cross contamination of specimen bottles.

Biopsy 1 (fresh): Half is put into viral transport medium, the other half is placed in OCT embedding media and frozen in liquid nitrogen, if it is available.

Biopsy 2 (fresh): *Placed in RNAlater Tissue Protect Tube (Qiagen®, Cat. No. 76163) ,if it is available.

Biopsy 3 – 70% Ethanol: Ethanol is a superior fixative to formalin for RNA extraction but is inferior for morphology. This fixative can be used if the tissue needs to be transported over

Annex 4:
Collection of autopsy tissues in resource-limited settings

This document should be read as a supplement to existing autopsy guidelines recommended by the US Centers for Disease Control and Prevention (CDC) for submission procedures (see: <http://www.cdc.gov/h1n1flu/tissuesubmission.htm>) and safety requirements (see: http://www.cdc.gov/h1n1flu/post_mortem.htm.)

The information below specifically supplements the above CDC information when rapid freezing and storage methods such as dry ice are not readily available, and also, if consent has only been given for a limited or para-mortem examination.

If consent has been given for a full autopsy/postmortem, then the relevant tissues should be ideally divided to be placed in appropriate containers for each of the following: (1) frozen tissue, (2) fixed in 70% alcohol (this gives better preservation of RNA than formalin), (3) placed in RNAlater Tissue Protect Tube (Qiagen®, Cat. No. 76163) for RNA preservation (see below), (4) virus transport medium for virus isolation and (5) 10% neutral buffered formalin.

For long term storage and transport, RNAlater (Ambion) Tissue Protect Tubes provide pre-measured volumes of RNAlater RNA Stabilization Reagent in re-closable tubes for convenient handling and sample storage. The reagent is **not** a fixative, but rather preserves RNA for up to one day at 37°C, seven days at 18–25°C, or four weeks at 2–8°C, allowing processing, transportation, storage, and shipping of samples without liquid nitrogen or dry ice. Alternatively, the samples can also be placed at –20°C or –80°C for archival storage. Morphological analysis, however, will be very difficult to interpret from tissues placed in this reagent.

Tissues to be sampled:

- Lung: all lobes, including affected areas
- Proximal and distal trachea
- Liver
- Brain
- Kidneys
- Nasopharynx (this is best done from a supratentorial approach with removal en bloc).

distance and liquid nitrogen is not available. (If the same Trucut needle is used for all biopsies fixation in ethanol should be done first as even a small amount of formalin in ethanol will impair RNA extraction).

Biopsy 4 and 5 – 10% neutral buffered formalin: This will allow good morphology to determine the extent of lung damage, secondary bacterial infection, etc. It will also enable immunohistochemistry to be performed for viral antibodies and co-localization studies to determine the cells infected with the infectious agent. Formalin fixed tissue can be used for RNA extraction but is less optimal as a fixative in this respect than alcohol. One of the biopsy samples should be processed for paraffin embedding; the other will remain in formalin for future extraction of RNA/DNA or electron microscope examination. It should be recognized that the RNA that is from formalin fixed tissue extracted will be degraded and one should not aim at amplifying fragments larger than 60-120 bp. If freshly prepared, 2.5% glutaraldehyde in cacodylate buffer (0.1 M sodium cacodylate-HCl buffer pH 7.4) is available; a 1-2 mm sample from the core can be placed directly into this fixative for electron microscope examination.

OTHER ORGANS

Sampling of other organs, especially the spleen, bone marrow and liver, would be important (at least for standard histology), and for detecting the presence of virus (by culture, RT-PCR and viral antigen by immunohistology), and to look for haemophagocytosis or changes in these organs. To avoid possible PCR carryover contamination a new needle should be used for each organ. This is also important for determining the extent of disease dissemination.

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1. Outbreak of swine-origin influenza A (H1N1) virus infection - Mexico, March-April 2009. *MMWR Morb Mortal Wkly Rep* 58, 467-470 (2009).
2. Update: drug susceptibility of swine-origin influenza A (H1N1) viruses, April 2009. *MMWR Morb Mortal Wkly Rep* 58, 433-435 (2009).
3. Baz, M., et al. Emergence of oseltamivir-resistant pandemic H1N1 virus during prophylaxis. *N Engl J Med* 361, 2296-2297 (2009).
4. Lackenby, A., et al. Emergence of resistance to oseltamivir among influenza A(H1N1) viruses in Europe. *Euro Surveill* 13(2008).
5. NITE/NIID. NITE/NIID Protocol for Sequencing Influenza A (H1N1) SWL Viral Genome Segments version 1.2 (30 May 2009). (2009).
6. 仙台市衛生研究所. 仙台市インフルエンザ等流行情報. (2009).
7. Nelson, M., et al. The early diversification of influenza A/H1N1pdm. *PLoS Curr Influenza*, RRN1126 (2009).
8. Hurt, A.C., Holien, J.K., Parker, M., Kelso, A. & Barr, I.G. Zanamivir-Resistant Influenza Viruses with a Novel Neuraminidase Mutation. *J Virol* (2009).
9. 国立感染症研究所感染症情報センター. <速報> 新型インフルエンザ (A/H1N1pdm) オセルタミビル耐性株 (H275Y) の国内発生状況 [第1報]. (2009).
10. Chen, H., et al. Oseltamivir-resistant influenza A pandemic (H1N1) 2009 virus, Hong Kong, China. *Emerg Infect Dis* 15, 1970-1972 (2009).
11. Hsu, S.B. & Hsieh, Y.H. On the role of asymptomatic infection in transmission dynamics of infectious diseases. *Bull Math Biol* 70, 134-155 (2008).
12. Cheng, V.C., et al. Prevention of nosocomial transmission of swine-origin pandemic influenza virus A/H1N1 by infection control bundle. *J Hosp Infect* (2010).
13. Williams, C.J., et al. Seasonal influenza risk in hospital healthcare workers is more strongly associated with household than occupational exposures: results from a prospective cohort study in Berlin, Germany, 2006/07. *BMC Infect Dis* 10, 8 (2010).

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Furuse Y, et al. Reversion of Influenza A(H3N2) from Amantadine-resistant to Amantadine-sensitive by Further reassortment in Japan during the 2006-2007 Influenza Season. *Journal of Clinical Microbiology* 2009, vol.47, No3, p841-844

Furuse Y, et al. Large-Scale Sequence Analysis of M Gene of Influenza A Viruses from Different Species: Mechanisms for Emergence and Spread of Amantadine Resistance. *Antimicrobial Agents and Chemotherapy* 2009 ;53(10):4457-63

Furuse Y, et al. Evolution of the M gene of the influenza A virus in different host species: large-scale sequence analysis. *Virology Journal* 2009 29;6:67

Furuse Y, et al. Reassortment between Swine Influenza A Viruses Increased their Adaptation to Humans in Pandemic H1N1/09. *Infection, Genetics & Evolution* (In Press)

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Mixed infection of influenza A viruses is common, 14th International Congress on Infectious Disease, March 2010, Miami, Florida

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