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Training course on genotyping and molecular analyses of antimicrobial resistant bacteria
Time Table

2015/02/02 Mon.		02/03 Tue.		02/04 Wed.		02/05 Thu.		02/06 Fri.	
Gram-negative bacteria	<i>C. difficile</i>	Gram-negative bacteria	<i>C. difficile</i>	Gram-negative bacteria	<i>C. difficile</i>	Gram-negative bacteria	<i>C. difficile</i>	Gram-negative bacteria	<i>C. difficile</i>
9		9	presentaion & discussion of participants	9	presentaion & discussion of participants	9	presentaion & discussion of participants	9	Lab
10		10	Lecture Antimicrobial sysceptibility testing & defection method for CRE	10	Lab results of disk diffusion test	10	Lab Electrophoresis	10	Lecture plasmid sequence GPAT
11		11	Lab Antimicrobial susceptibility testing	11	Lab DNA extraction	11	Lab Electrophoresis	11	Culture of <i>Clostridium difficile</i> from stool specimens
12		12	Lunch	12	Lunch	12	Lunch	12	Lunch
13	Overview Epidemiology and Surveillance of Antimicrobial resistant bacteria in Japan (Dr. S Suzuki)	13	Lunch	13	Lecture phenotype&genotype, typing method	13	Lab	13	Discussion
14	Lecture Biosafety (Dr. Tanabayashi) room 204	14	Lab disk diffusion test CarbaNP test	14	Lab PCR results of disk diffusion test	14	Lab comparison of pgenotype & genotype. PFGE, sequence (Demo)	14	Free
15		15		15		15		15	
16	Lab Labo Tour & Automatic system for bacteria identification (MALDI- TOF MS) (Dr. M Suzuki & Matsui)	16	PCR	16	Making gel for PCR ribotyping	16		16	
17		17		17		17		17	

Antimicrobial Susceptibility Testing & Detection Methods for Gram-Negative Rods

Objectives of AMR test at NIID

- 1 To reveal mechanisms of resistance
- 2 To confirm nosocomial outbreak

Before starting the tests...

1. Confirm the isolate information
 1. Number of isolates
 2. Species
 3. Resistance

2. Confirm Epidemiological information for each isolate

1. Date of isolation (sampling), at least "Year"
 2. Place (hospital, ward)
 3. Type of specimen (Blood, Urine etc.)
 4. Patient ID
3. Inoculation of the isolate
- One isolate per one agar plate
 - Use antimicrobial drug containing disk
 - Check the single colony, antimicrobial resistance
 - Confirm species identification and antimicrobial resistance
4. Preserve the isolates and record the data into your database.

Let's Start!

Antimicrobial Susceptibility Testing and Detection at reference lab


At clinical lab
Predict/speculate resistance mechanism from antimicrobial susceptibility pattern

At reference lab
Confirm the resistance mechanisms by molecular method (usually detection of resistance gene by PCR)
Check the whether the phenotype (antimicrobial susceptibility pattern, disk diffusion test, Carba NP test etc.) of isolates are consistent with the results of genotype.

- > Molecular typing
- ✓ Pulsed Field Gel Electrophoresis (PFGE)
 - ✓ Multilocus Sequence Typing (MLST)
 - ✓ Plasmid analysis


Antimicrobial susceptibility Test (AST)

- In Japan...
 - AST is usually conducted at clinical level
 - Middle to large hospitals, large commercial labs
⇒ microdilution method with automated system
 - Small hospitals and small commercial labs
⇒ disk diffusion test

Disk diffusion test


Relatively cheap
Continuous data available
Accurate MIC is not available
Difficulty in automation

Etest

Expensive
MIC is available

Microdilution Broth method

Use 96 well plate
Add 100µl of bacterial suspension
↓
Incubation
Read the plate for bacterial growth

Automation system is available.
Dry plate and frozen plate are commercially available.

96 well antimicrobial locating design

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Break point panel

- When the S/I/R category is defined as
Susceptible $\leq 1\mu\text{g/ml}$
Intermediate $= 2\mu\text{g/ml}$
Resistant $\geq 4\mu\text{g/ml}$

Meropenem	54	32	16	8	4	2	1	0.5	0.25	0.125

EUCAST screening cut-off for detecting
carbapenemase-producing Enterobacteriaceae
MIC $>0.12 \mu\text{g/ml}$ (meropenem)
Disk diffusion $<25\text{mm}$

Detection of
acquired β -lactam resistance among GNR

- Disk diffusion test with β -lactamase inhibitor
- Carba NP test
- Detection of β -lactamase genes by PCR

Check the consistency

This is the main objectives of this training course


Additional analysis

- Plasmid analysis, Whole genome analysis


Classification of β -lactamase

Ambler class	genotype	resistance phenotypes					Inhibition by				
		2 nd generation carbapenemase (KPC, CTX)	4 th generation carbapenemase	carbapenemase (KPC, CTX)	carbapenemase (KPC, CTX)	carbapenemase (KPC, CTX)	clavulanic acid	EDTA	DPA	PBA	clavulanic acid
A	TEM SHV (ESBL)	R (S CTX)	R	R	S	S	+	-	-	-	
	CTXM (ESBL)	R (CAZ)	R	R	S	S	+	-	-	-	
	KPC	R	R	R	R	R	-	-	+	-	
B	Bla _{TEM} SHV	R	R	S	R	R	-	+	-	-	
C	CMY etc...	R	I/R	R	I/R	S	-	-	+	+	
	OXA-48	S	S	S	S	R	-	-	-	-	

SMA: sodium mercaptoacetate acid, DPA: dipicolinic acid
ACV: amoxicillin/clavulanate, CAZ: ceftazidime, CTX: cefotaxime, CMZ: cefmetazole, CMN: cefminox, IPM: imipenem, MEPM: n

Carba NP test

Why it is important to know the genotype of resistance?

- Knowing the molecular epidemiology enables
 - Easy understanding of resistant mechanisms at clinical lab
 - Easy selection of resistant gene to do the PCR detection
 - Easy understanding of nosocomial outbreak of resistant bacteria

Nationwide survey in 2010

Between September and December in 2010
Enterobacteriaceae isolates resistant to carbapenems, aminoglycosides and fluoroquinolones were collected from hospitals across Japan.

	IMP-1 type	Genotype of Carbapenemase			
		KPC	NDM-1	OXA-48 like	SMB-1*
<i>Escherichia coli</i>	23				
<i>Klebsiella pneumoniae</i>	19	2	2	1	
<i>Enterobacter cloacae</i>	22				
**Other	8				1
Total	72	2	2	1	1

*Newly identified metallo-β-lactamase (Vlachino et al. AAC 2011)
** *Providencia* spp., *Citrobacter freundii*, *Serratia marcescens*, *Morganella morganii*, *Proteus mirabilis*, *Klebsiella oxytoca*

•IMP-1 type metallo-β-lactamase producing Enterobacteriaceae is already endemic in Japan.
•Most of KPC, NDM and OXA-48 producers are imported cases.

IMP-type carbapenemase domestic thread of Japan

- Metallo-β-lactamase, belong to Class B of the structural classification of β-lactamase
- First detected in clinical *Serratia marcescens* isolate from Japan in the late 1980s. The first acquired metallo-β-lactamase
- More than 30 different allotypes have been described. (IMP-1, IMP-2 . . .)
- The different IMP variants often have a defined geographical distribution.
- In Japan, IMP-1 metallo-β-lactamase have become endemic among carbapenem-resistant *Pseudomonas aeruginosa* since late 1990s.
- Most of the hospital microbiologist are aware of detecting metallo-β-lactamase by disk containing sodium mercaptoacetic acid (SMA)


CPE is not always CRE

MIC distribution of IMP-1-type metallo-β-lactamase producing Enterobacteriaceae

In Japan, MICs are more frequently measured for imipenem

MIC distribution of IMP-1/IMP-6 metallo-β-lactamase producers.

Classification of β -lactamase genes

Ambler class	gene	resistance phenotypes						Inhibition by			
		penicillin + β -lactamase inhibitor (e.g. ACV)	3 rd generation cephalosporin (e.g. CAZ, CTX)	4 th generation cephalosporin	monobactam (e.g. aztreonam)	cephamycin (e.g. CMZ, CMN)	carbapenem (e.g. IPM, MEPM)	clavulanate	EDTA /SMA /DPA	PBA	cloxacillin
A	TEM SHV (ESBL)	S	R (S: CTX)	R	R	S	S	+	-	-	-
	CTX-M (ESBL)	S	R (S: CAZ)	R	R	S	S	+	-	-	-
	KPC	R	R	R	R	R	R	-	-	+	-
B (metallo- β -lactamase)	NDM IMP VIM TMB	R	R	R	S	R	R	-	+	-	-
C (AmpC β -lactamase)	CMY etc..	R	R	S	R	I/R	S	-	-	+	+
D	OXA-48	R	S	S	S	S	R	-	-	-	-

SMA: sodium mercaptoacetic acid, DPA: dipicolinic acid

ACV: amoxicillin/clavulanate, CAZ: ceftazidime, CTX: cefotaxime, CMZ: cefmetazole, CMN: cefminox, IPM: imipenem, MEPM: meropenem

microbroth dilution

colony suspension in 3ml Mueller-Hinton broth (ca. McFarland=1 / OD₆₀₀= 0.257)

↓

dilution (add 25 μ l to 12ml Mueller-Hinton broth)

↓


to plate (100 μ l/well)

↓


incubate at 37°C overnight

	1	2	3	4	5	6	7	8	9	10	11	12
A	CPFX 32	IPM 64	32	16	8	4	2	1	0.5	0.25	0.125	S/A 32/64
B	16	MEPM 64	32	16	8	4	2	1	0.5	0.25	0.125	16/32
C	8	LVFX 32	16	8	4	2	1	0.5	0.25	CAZ 32	CFPM 32	8/16
D	4	AMK 128	64	32	16	8	4	2	PL 4	16	16	4/8
E	2	CL 4	MINO 32	16	8	4	2	1	2	8	8	2/4
F	1	2	PIP/TAZ 256/4	128/4	64/4	32/4	16/4	8/4	1	4	4	1/2
G	0.5	1	S/C 32/64	16/32	8/16	4/8	2/4	1/2	0.5	2	2	0.5/1
H	0.25	0.5	GM 64	32	16	8	4	2	1	1	1	GC


Plasmid
The drive force of resistance gene dissemination


Inter-species transmission of resistant plasmid

“Typical” nosocomial transmission of resistant bacteria



Nosocomial outbreak of CRE ?

Is this a nosocomial transmission??

S1 PFGE Southern blotting

Carba NP test

Aim

To detect a carbapenemase by change a color of phenol red

Materials

- bacteria on a plate (overnight culture)
- Lysis buffer (B-PERII, Bacterial Protein Extraction Reagent; Thermo)
- Imipenem-phenol red solution
- 96-well plate or 1.5 ml tube

Preparation of imipenem-phenol red solution

1. Dissolve 50 mg phenol red in 1 ml pure water after adding 1 drop 1N NaOH.
2. Add 9 ml pure water (0.5% phenol red solution).
3. Mix 2 ml 0.5% phenol red solution and 16.6 ml pure water, and adjust the pH 7.8 with 1N NaOH (drop each 3 μ L)
4. Dissolve 288 mg ZnSO₄ in 10 ml pure water.
5. Mix 2ml phenol red solution prepared in step 3 and 2 μ L ZnSO₄ solution prepared in step 4.
6. Dissolve 3 mg imipenem in 1 ml solution prepared in step 5.

Methods

1. Dispense 100 μ L lysis buffer in 1.5 ml tube or 96-well plate.
2. Scrape off with a 10- μ L loop and suspend in lysis buffer (two times).
3. Vortex for 5 min.
4. Incubate for 30 min, RT
5. Dispense 100 μ L imipenem-phenol red solution in 1.5 ml tube or 96-well plate.
6. A 30 μ L-volume of this lysate is mixed with imipenem-phenol red solution.
7. Incubate for 60-120 min, at 37 degree.
8. The color check with visual observation.

Reference

1. Patrice Nordmann *et al.* 2012. *Emerg Infect Dis.* 18:1503-7
2. Laurent Dortet *et al.* 2012. *Antimicrob. Agents Chemother.* 56:6437-6440
3. Monica O' sterblad *et al.* 2014. *Antimicrob. Agents Chemother.* 58:7553-7556

	1	2	3	4	5	6	7	8	9	10	11	12
A	CPFX 32	IPM 64	32	16	8	4	2	1	0.5	0.25	0.125	S/A 32/64
B	16	MEPM 64	32	16	8	4	2	1	0.5	0.25	0.125	16/32
C	8	LVFX 32	16	8	4	2	1	0.5	0.25	CAZ 32	CFPM 32	8/16
D	4	AMK 128	64	32	16	8	4	2	PL 4	16	16	4/8
E	2	CL 4	MINO 32	16	8	4	2	1	2	8	8	2/4
F	1	2	PIPC/TAZ 256/4	128/4	64/4	32/4	16/4	8/4	1	4	4	1/2
G	0.5	1	S/C 32/64	16/32	8/16	4/8	2/4	1/2	0.5	2	2	0.5/1
H	0.25	0.5	GM 64	32	16	8	4	2	1	1	1	GC

isolate No. _____

		MIC(μ g/ml)
IPM	imipenem	
MEPM	meropenem	
LVFX	levofloxacin	
CPFX	ciprofloxacin	
AMK	amikacin	
CL	colistin	
MINO	minocycline	
PIPC /TAZ	piperacillin/tazobactam	
S/C	sulbactam/cefoperazone	
GM	gentamicin	
PL	polymyxin B	
CAZ	ceftazidime	
CFPM	cefepime	
S/A	sulbactam/ampicillin	

isolates (AMR course at NIID)

	isolates No.	species	Memo (resistance genes etc...)
MDRP	10-654	<i>Pseudomonas aeruginosa</i>	
	10-655	<i>Pseudomonas aeruginosa</i>	
	10-747	<i>Pseudomonas aeruginosa</i>	
MDRA	12-277	<i>Acinetobacter baumannii</i>	
	09-775	<i>Acinetobacter baumannii</i>	
CRE	13-273	<i>Klebsiella pneumoniae</i>	
	10-810	<i>Klebsiella pneumoniae</i>	
	10-811	<i>Escherichia coli</i>	
	10-822	<i>Klebsiella pneumoniae</i>	
	11-327	<i>Klebsiella pneumoniae</i>	
	12-143	<i>Escherichia coli</i>	
	13-367	<i>Klebsiella pneumoniae</i>	
	06-358	<i>Escherichia coli</i>	
	06-359	<i>Klebsiella pneumoniae</i>	
	06-355	<i>Klebsiella pneumoniae</i>	
	08-1508	<i>Escherichia coli</i>	

PCR

1. DNA extraction

plate culture (incubated overnight at 37°C)

↓

colony suspension in 500µl sterile water (ca. McFarland 0.5)

↓

boil 100°C, 10min

↓

centrifugation 13000rpm, 4°C, 5min

↓

collect the supernatant (100µl) to other tube

2. reaction mix

PCR sheet

Date _____

		×	(number of samples+3)	Total volume
PCR beads*	1			
sterile water	22µl			
Forward primer (20µM)	1µl			
Reverse primer (20µM)	1µl			
template DNA	1µl			
Total	25µl			

*illustra Hot Start Mix RTG Bead, which includes: dNTPs, PCR buffer, Mg⁺⁺, Taq(enzyme)

primer type (gene type) _____

#1	13-273 <i>K. pneumoniae</i>
2	10-810 <i>K. pneumoniae</i>
3	10-822 <i>K. pneumoniae</i>
4	11-327 <i>K. pneumoniae</i>
5	10-654 <i>P. aeruginosa</i>
6	positive control
7	negative control

primer list (carbapenemase genes)

gene type	primer sequence (5' →3')	amplicon size
NDM type	F: TTGCCAATATTATGCACCC R: ATTGGCATAAGTCGCAATCC	420bp
KPC type	F: ATGTCACTGTATCGCCGTCT R: TTTTCAGAGCCTTACTGCC	893bp
IMP-1 type	F: ACCGCAGCAGAGTCTTTGCC R: ACAACCAGTTTTGCCTTACC	587bp
IMP-2 type	F: GTTTTATGTGTATGCTTCC R: AGCCTGTTCCCATGTAC	678bp
VIM-2 type	F: ATGTTCAAACCTTTGAGTAAG R: CTAACAACGACTGAGCG	801bp
OXA-48-like	F: TTGGTGGCATCGATTATCGG R: GAGCACTTCTTTGTGATGGC	744bp

PCR program

94°C	2min	} 30 cycle
94°C	1min	
55°C	1min	
72°C	1min30sec	
72°C	5min	
4°C	∞	

Reference

- Shibata N *et al.* 2003. J Clin Microbiol. 41(12):5407-13
- Poirel L *et al.* 2004. Antimicrob Agents Chemother. 48(1): 15-22