

Fig. 2. RET Week 4: *Pig-a* MF dose response data following exposure to ENU from the different laboratories was analysed using the BMD covariate approach, using BMRs of (a) 10% and (b) 100%.

One MMS dose response data set from Zeller et al. [17] was used to increase the precision of the BMD estimates [7]. The 4 parameter exponential (top horizontal line) and Hill (bottom horizontal line) models provided a suitable fit to the data, with 'laboratory' used as covariate. The width of the horizontal lines represents the BMDL-BMDU, which are ranked from lowest to highest concentration by BMD. During this combined analysis, the maximum response (parameter c) and log-steepness (d) parameters were assumed equal for all response curves, while parameters a (background response), b (potency) and var (within group variation) were covariate dependant. The use of constant 'shape' parameters (parameters c, d) still provided a strong description of the individual response curves. Overlapping lines show equipotency, with potency decreasing from top left to bottom right. CED-0.1 is equivalent to, or another name for, BMD₁₀. X-axis are Log₁₀.dose (mg/kg/day); Y-axis are laboratory.

203 results, although in three datasets (labs) the confidence intervals
204 related to the BMD for RET were relatively large (probably due
205 to all observations in the controls being zero, *i.e.* below limit of
206 quantification, LOQ). In both figures the confidence intervals do
207 not overlap among all labs (datasets), indicating differences among
208 some labs. Moreover, based on visual inspection of the figures, these
209 differences are roughly estimated to be within only one order of
210 magnitude.

211 3.1. Discussion

212 The BMD potency ranking plots provide information on the
213 reproducibility of the *Pig-a* dose response data. The overlapping
214 confidence intervals established from mutant RBC and RET fre-
215 quencies at week 4 show evidence that the mutant phenotype
216 population can be reasonably well reproduced in different datasets,
217 *i.e.* within one order of magnitude. Although the differences may be
218 due to the labs, they could also just be replication error. This could
219 not be established here, as no replicate studies within labs were
220 available.

221 This is the first instance where potency estimates using BMD
222 covariate analysis has been used to examine interlaboratory repro-
223 ducibility. Non-overlapping BMD metrics and relatively large
224 differences between the potential values of the true BMDs would
225 indicate that the level of reproducibility was low. In the present
226 ring study, we found reasonable reproducibility, but it would be
227 worthwhile to improve it. This could be achieved by increasing the
228 number of cells scored to a minimum of $1-5 \times 10^6$ RETs or RBCs
229 per animal, as this would reduce the number of zero's within the
230 data [20]. In this regard, laboratory A was the only laboratory not
231 to contain zero data points in the control *Pig-a* replicates of RETs,
232 and laboratory A also has the smallest BMDL-BMDU width. For the
233 RBC *Pig-a* MF, the laboratories that did not contain any zero's in
234 control were B, G, L2 and O, and the BMDL-BMDU for each of these
235 laboratories overlaps very well with each other (Fig. 1).

236 Another example of the BMD covariate analysis approach in
237 genetic toxicology is provided in the recent paper by Wills et al.,
238 who found no significant differences between BMD potency esti-
239 mates in different experimental replicates [7]. As another example
240 the approach was used to show the effect of sampling day in *in vivo*
241 *Pig-a* data. The approach is robust and provides a suitable way of
242 comparing potencies across covariates (*e.g.* laboratories, sampling
243 day, compound, sex, species *etc.*).

244 Acknowledgements

245 We thank Robert Heflich and the ILSI/HESI Genetic Toxicology
246 Technical Committee (GTTC) for initiating this collaboration.


247 Appendix A. Supplementary data

248 Supplementary data associated with this article can be found, in
249 the online version, at <http://dx.doi.org/10.1016/j.mrgentox.2016.04.004>.

251 References

- 252 [1] S.D. Dertinger, S. Phonethpawth, P. Weller, J. Nicolette, J. Murray, P. Sonders,
253 H.W. Vohr, J. Shi, L. Krsmanovic, C. Gleason, L. Custer, A. Henwood, K. Sweder,
254 L.F. Stankowski Jr., D.J. Roberts, A. Giddings, J. Kenny, A.M. Lynch, C. Defrain, F.
255 Nesslany, B.J. van der Leede, T. Van Doninck, A. Schuermans, K. Tanaka, Y.
256 Hiwata, O. Tajima, E. Wilde, A. Elhajouji, W.C. Gunther, C.J. Thiffeault, T.J.
257 Shutsky, R.D. Fiedler, T. Kimoto, J.A. Bhalli, R.H. Heflich, J.T. MacGregor,
258 International *Pig-a* gene mutation assay trial: evaluation of transferability
259 across 14 laboratories, *Environ. Mol. Mutagen.* 52 (2011) 690-698.
- 260 [2] B.B. Gollapudi, G.E. Johnson, L.G. Hernandez, L.H. Pottenger, K.L. Dearfield,
261 A.M. Jeffrey, E. Julien, J.H. Kim, D.P. Lovell, J.T. Macgregor, M.M. Moore, J. van
262 Benthem, P.A. White, E. Zeiger, V. Thybaud, Quantitative approaches for
263 assessing dose-response relationships in genetic toxicology studies, *Environ.*
264 *Mol. Mutagen.* 54 (2013) 8-18.
- 265 [3] G.E. Johnson, L.G. Soeteman-Hernandez, B.B. Gollapudi, O.G. Bodger, K.L.
266 Dearfield, R.H. Heflich, J.G. Hixon, D.P. Lovell, J.T. MacGregor, L.H. Pottenger,
267 C.M. Thompson, L. Abraham, V. Thybaud, J.Y. Tanir, E. Zeiger, J. van Benthem,
268 P.A. White, Derivation of point of departure (PoD) estimates in genetic


- 269 toxicology studies and their potential applications in risk assessment,
270 Environ. Mol. Mutagen. 55 (2014) 609-623.
- 271 [4] J.T. MacGregor, R. Frötschl, P.A. White, K.S. Crump, D.A. Eastmond, S.
272 Fukushima, M. Guérard, M. Hayashi, L.G. Soeteman-Hernandez, G.E. Johnson,
273 T. Kasamatsu, D. Levy, T. Morita, L. Müller, R. Schoeny, M.J. Schuler, V.
274 Thybaud, IWGT report on quantitative approaches to genotoxicity risk
275 assessment II. Use of point-of-departure (PoD) metrics in defining acceptable
276 exposure limits and assessing human risk, Mutat. Res. 783 (2015) 66-78.
- 277 [5] J.T. MacGregor, R. Frötschl, P.A. White, K.S. Crump, D.A. Eastmond, S.
278 Fukushima, M. Guérard, M. Hayashi, L.G. Soeteman-Hernandez, T. Kasamatsu,
279 D. Levy, T. Morita, L. Müller, R. Schoeny, M.J. Schuler, V. Thybaud, G.E. Johnson,
280 IWGT report on quantitative approaches to genotoxicity risk assessment I.
281 Methods and metrics for defining exposure-response relationships and points
282 of departure (PoDs), Mutat. Res. 783 (2015) 55-65.
- 283 [6] J.W. Wills, G.E. Johnson, S.H. Doak, L.G. Soeteman-Hernandez, W. Slob, P.A.
284 White, Empirical analysis of BMD metrics in genetic toxicology part I: *in vitro*
285 analyses to provide robust potency rankings and support MOA
286 determinations, Mutagenesis 31 (2016).
- 287 [7] J.W. Wills, A.S. Long, G.E. Johnson, J.C. Bemis, S.D. Dertinger, W. Slob, P.A.
288 White, Empirical analysis of BMD metrics in genetic toxicology part II: *in vivo*
289 potency comparisons to enhance the utility of experimental animals for
290 genetic toxicity assessment, Mutagenesis 31 (2016).
- 291 [8] L.G. Hernández, J. Van Benthem, W. Slob, Estimating the carcinogenic potency
292 of chemicals from the *in vivo* micronucleus test: RIVM Report
293 340700007/2012 RIVM Report (2012).
- 294 [9] L.G. Hernández, W. Slob, H. van Steeg, J. van Benthem, Can carcinogenic
295 potency be predicted from *in vivo* genotoxicity data? a meta-analysis of
296 historical data, Environ. Mol. Mutagen. 52 (2011) 518-528.
- 297 [10] L.G. Hernández, G.E. Johnson, L.H. Pottenger, J. van Benthem, Analysis of
298 low-dose mutagenic responses and the applicability of genotoxicity tests for
299 carcinogen potency prediction, Environ. Mol. Mutagen. 52 (2011) S26.
- 300 [11] L.G. Hernández, W. Slob, H. van Steeg, J. van Benthem, Comparison of
301 carcinogenic potency estimates to *in vivo* genotoxic potencies from the
302 micronucleus, transgenic rodent mutation and comet assay using the
303 benchmark dose approach, Environ. Mol. Mutagen. 51 (2010) 707.
- 304 [12] L. Soeteman-Hernandez, M. Fellows, G.E. Johnson, W. Slob, Carcinogenic
305 potency estimation from the *in vitro* micronucleus test in TK6 cells: a pilot
306 study, Toxicol. Sci. (2015), <http://dx.doi.org/10.1093/toxsci/kfv189>.
- 307 [13] G.E. Johnson, W. Slob, S.H. Doak, M.D. Fellows, B.B. Gollapudi, R.H. Heflich, B.J.
308 Rees, L.G. Soeteman-Hernandez, J.R. Verma, J. Wills, G.J.S. Jenkins, P.A. White,
309 New approaches to advance the use of genetic toxicology analyses for human
310 health risk assessment, Toxicol. Res. 4 (2015) 667-676.
- 311 [14] L.G. Soeteman-Hernandez, G.E. Johnson, W. Slob, Estimating the carcinogenic
312 potency of chemicals from the *in vivo* micronucleus test, Mutagenesis (2015),
313 <http://dx.doi.org/10.1093/mutage/gev043>.
- 314 [15] W. Slob, R.W. Setzer, Shape and steepness of toxicological dose-response
315 relationships of continuous endpoints, Crit. Rev. Toxicol. 44 (2014) 270-297.
- 316 [16] D. Miura, et al., Potential value of the *Pig-a* gene mutation assay in
317 reticulocytes (PIGRET assay) using rats as a short-term genotoxicity test:
318 Summary of the collaborative study across 16 laboratories using 24 chemicals
319 by the Mammalian Mutagenicity Study (MMS) Group of the Japanese
320 Environmental Mutagenicity Society (JEMS), Mutat. Res. (2016) () in this issue.
- 321 [17] A. Zeller, L. Tang, S.D. Dertinger, J. Funk, G. Duran-Pacheco, M. Guerdard, A
322 proposal for a novel rationale for critical effect size in dose-response analysis
323 based on a multi-endpoint *in vivo* study with methyl methanesulfonate,
324 Mutagenesis 31 (2016).
- 325 [18] EFSA, Use of benchmark dose approach in risk assessment: Guidance of the
326 Scientific Committee, 2009, 1150 (2009) 1-72.
- 327 [19] J.C. Bemis, J.W. Wills, S.M. Bryce, D.K. Torous, S.D. Dertinger, W. Slob,
328 Comparison of *in vitro* and *in vivo* clastogenic potency based on benchmark
329 dose analysis of flow cytometric micronucleus data, Mutagenesis 31 (2016).
- 330 [20] B.B. Gollapudi, A.M. Lynch, R.H. Heflich, S.D. Dertinger, V.N. Dobrovolsky, R.
331 Froetschl, K. Horibata, M.O. Kenyon, T. Kimoto, D.P. Lovell, L.F. Stankowski Jr.,
332 P.A. White, K.L. Witt, J.Y. Tanir, The *in vivo Pig-a* assay: a report of the
333 international workshop on genotoxicity testing (IWGT) workgroup, Mutat.
334 Res. 783 (2015) 23-35.



Contents lists available at ScienceDirect

Mutation Research/Genetic Toxicology and Environmental Mutagenesis

journal homepage: www.elsevier.com/locate/gen tox
 Community address: www.elsevier.com/locate/mutres



The *in vivo* Pig-a assay: A report of the International Workshop On Genotoxicity Testing (IWGT) Workgroup[☆]



B. Bhaskar Gollapudi^{a,1}, Anthony M. Lynch^{b,1}, Robert H. Heflich^{c,*,2},
 Stephen D. Dertinger^d, Vasily N. Dobrovolsky^c, Roland Froetschl^e, Katsuyoshi Horibata^f,
 Michelle O. Kenyon^g, Takafumi Kimoto^h, David P. Lovellⁱ, Leon F. Stankowski Jr.^j,
 Paul A. White^k, Kristine L. Witt^l, Jennifer Y. Tanir^m

^a Exponent, Midland, MI, USA^b GlaxoSmithKline, Hertfordshire, UK^c US Food and Drug Administration, National Center for Toxicological Research, Jefferson, AR, USA^d Litron Laboratories, Rochester, NY, USA^e Federal Institute for Drugs and Medical Devices, Bonn, Germany^f National Institute of Health Sciences, Tokyo, Japan^g Pfizer Global Research and Development, Groton, CT, USA^h Teijin Pharma, Tokyo, Japanⁱ St. George's, University of London, UK^j BioReliance, Rockville, MD, USA^k Environmental Health Science and Research Bureau, Health Canada, Ottawa, ON, Canada^l National Institutes of Health, National Institute of Environmental Health Sciences, Division of the National Toxicology Program, Research Triangle Park, NC, USA^m Health and Environmental Sciences Institute, Washington, DC, USA

ARTICLE INFO

Article history:

Received 13 September 2014

Accepted 15 September 2014

Available online 26 September 2014

Keywords:

Red blood cells

Reticulocytes

Glycosylphosphatidylinositol

Flow cytometry

CD59

Mutation

ABSTRACT

The *in vivo* Pig-a assay uses flow cytometry to measure phenotypic variants for antibody binding to cell surface glycosylphosphatidylinositol (GPI)-anchored proteins. There is good evidence suggesting that the absence of antibody binding is the result of a mutation in the endogenous X-linked *Pig-a* gene, which forms the rationale for the assay. Although the assay has been performed with several types of hematopoietic cells and in a variety of mammalian species, including humans, currently it is optimized only for measuring CD59-deficient (presumed *Pig-a* mutant) erythrocytes in the peripheral blood of rats. An expert workgroup formed by the International Workshop on Genotoxicity Testing considered the state of assay development and the potential of the assay for regulatory use. Consensus was reached on what is known about the *Pig-a* assay and how it should be conducted, and recommendations were made on additional data and refinements that would help to further enhance the assay for use in hazard identification and risk assessment.

Published by Elsevier B.V.

[☆] The views expressed in this report are not necessarily shared by the institutions at which the workgroup members are employed. The information in these materials is not a formal dissemination of information by the U.S. Food and Drug Administration and does not represent agency position or policy. Data presented in this communication were partially funded by a grant from the National Institute of Health/National Institute of Environmental Health Sciences (NIEHS; grant no. R44ES021973). The contents of this report are solely the responsibility of the authors, and do not necessarily represent the official views of the NIEHS.

* Corresponding author at: U.S. Food and Drug Administration, NCTR, 3900 NCTR Road, Jefferson, AR 72079, USA. Tel.: +1 870 543 7493.

E-mail address: robert.heflich@fda.hhs.gov (R.H. Heflich).

¹ Co-Chairs.

² Rapporteur.

<http://dx.doi.org/10.1016/j.mrgentox.2014.09.007>

1383-5718/Published by Elsevier B.V.

1. Introduction

1.1. Background to this report

The *in vivo* *Pig-a* assay was first described for rats and mice in 2008 [1–4], and received almost immediate interest as a potential assay for evaluating the *in vivo* mutagenicity of new and existing substances. Subsequently, formal and informal workshops and presentations on the assay were held at meetings of the US Environmental Mutagenesis and Genomics Society and the Health and Environmental Sciences Institute (HESI) of the International Life Sciences Institute [5]. A Workgroup, made up of experts from academic, regulatory, and industrial laboratories, was formed in 2012 under the auspices of the International Workshop on Genotoxicity Testing (IWGT) to review the development of the *Pig-a* assay in the context of safety assessment strategies. The remit of the Workgroup was to consider the underlying science of the *Pig-a* assay and technical considerations for the assay protocol with a view to assay acceptance in a regulatory context and to recommend where and how further progress on developing the

assay can be made. Subsequently, two meetings of the Workgroup were held in 2013, one in conjunction with the HESI Genetic Toxicology Technical Committee annual meeting in Washington, DC, on April 15 and the second at the IWGT meeting in Foz do Iguaçu, Brazil, on October 31–November 1. This report describes consensus statements on the *Pig-a* assay developed by the IWGT Workgroup.

1.2. Principle of the assay

The *Pig-a* assay is based on the identification of mutant cells that have an altered repertoire of cell surface markers (Fig. 1). The assay was developed from an understanding of the molecular nature of a rare human acquired genetic disorder, paroxysmal nocturnal hemoglobinuria (PNH) (reviewed [6–8]). The *Pig-a* gene (phosphatidylinositol glycan, class A gene) codes for the catalytic subunit of an *N*-acetyl glucosamine transferase that is involved in an early step of glycosylphosphatidylinositol (GPI) biosynthesis. GPI anchors an assortment of protein markers (e.g., CD59, CD24, CD55) to the exterior surface of the cytoplasmic membranes of higher eukaryotes. In mammals, *Pig-a* (nomenclature: *PIG-A*




Fig. 1. Principle of the rat erythrocyte *Pig-a* assay. In wild-type nucleated erythroid precursor cells (upper left), glycosylphosphatidylinositol (GPI) anchors are synthesized in a series of steps at the endoplasmic reticulum (ER). The first step is catalyzed by the product of the X-linked *Pig-a* gene; all other steps of the GPI synthesis pathway are catalyzed by enzymes encoded by autosomal genes. The GPI core structure consists of 3 mannoses, 1 glucosamine, and 1 phosphatidylinositol, with 2 or 3 hydrophobic tails imbedded in the ER lipid bilayer membrane (a more detailed description of GPI structure and synthesis can be found in [7]). In the final stage of biosynthesis, a surface marker (e.g., CD59 protein, which is synthesized independently) is attached to the GPI anchor and the assembly is transported to the surface of the cell. *Pig-a* mutant cells (upper right) are deficient in GPI anchor synthesis, and thus GPI-anchored markers are not transported to the cytoplasmic membrane. The *Pig-a* assay is conducted with peripheral blood erythrocytes having no nuclei (lower, left and right). At the surface of the cytoplasmic membrane, GPI anchored CD59 reacts with fluorescent anti-CD59 antibodies which make the whole wild-type cell fluoresce (lower left). *Pig-a* mutant cells do not express CD59 on their surface and thus are not fluorescently labeled (lower right). Flow cytometry is used to quantify the number of wild-type and mutant cells and estimate a *Pig-a* mutant frequency.

in humans and other primates; *Pig-a* in mice and rats) is an X-linked gene present as a single functional copy in cells from both males and females. Other genes involved in GPI biosynthesis are autosomal and have two functional copies. A single inactivating mutation in the *Pig-a* gene is sufficient to make a cell deficient in GPI anchors and, as a consequence, deficient in surface-bound GPI-anchored markers. Since it is exceedingly unlikely that anchor deficiency would occur due to inactivating mutations in both copies of the autosomal genes involved in GPI synthesis, measuring GPI deficiency is considered 'virtually equivalent' to measuring *Pig-a* mutation.

The assay distinguishes the *Pig-a* mutant phenotype from the wild-type phenotype by labeling cells with fluorescent antibodies against a GPI-anchored marker or markers [8] (Fig. 1). Wild-type cells bind marker-specific antibodies and fluoresce, while mutant cells do not bind the antibodies and do not fluoresce. Cells are analyzed using high throughput flow cytometry to quantify wild-type (*i.e.* labeled) and rare mutant (*i.e.* unlabeled) cells. There are multiple commercially available antibodies, stains, and reagents suitable for designing *Pig-a* assays in different cell types and in different mammalian species.

The use of flow cytometry for detecting cells deficient in surface markers imposes certain constraints on the tissues amenable for the *Pig-a* assay. For instance, the samples must be prepared as single-cell suspensions, and the cells should not be subjected to treatments that may alter the cell surface membrane (*e.g.*, proteolytic digestion). Currently, these two requirements have limited the *in vivo* application of the *Pig-a* assay to hematopoietic tissue (peripheral blood erythrocytes and white blood cells, and, to a lesser degree, bone marrow).

A combination of theory and experimental evidence indicates that *Pig-a* mutations occurring in bone marrow cells result in the emergence of the mutant cells that are measured in peripheral blood [9,10]. The appearance of *Pig-a* mutant cells in peripheral blood exhibits a time delay that is dependent upon the cell and GPI-anchor turnover rates, as well as the trafficking time from the bone marrow to the periphery for the specific cell type. Red blood cells (RBCs) are the most practical cell type for performing the assay since microliter volumes of blood contain sufficient quantities of RBCs for the enumeration of the rare mutants required for a successful *Pig-a* assay. The abundance of RBCs in peripheral blood allows serial blood collections from even small laboratory animals and thus permits longitudinal studies to be conducted on the same set of animals. In addition, the small fraction of newly formed erythrocytes, *i.e.*, reticulocytes (RETs), can be distinguished by appropriate staining. RETs have a rapid turnover rate (few days), and express the mutant phenotype originating in the bone marrow faster than the population of total RBCs in peripheral blood. Finally, at least for the rat assay, *Pig-a* mutant erythrocytes act as though they have a neutral phenotype [9], meaning that mutant cells are not subjected to negative selection and can accumulate with repeated dosing, potentially increasing the sensitivity of the assay. Observations concerning *Pig-a* mutant manifestation in the peripheral blood and mutant accumulation with repeat doses are consistent with what occurs in other *in vivo* gene mutation assays [11,12], and lend support to the mutational origin of the phenotype measured in the *Pig-a* assay.

1.3. Potential for translation of the endpoint from experimental models to humans

The pathway for GPI biosynthesis is conserved in most mammalian species, including common laboratory animals, such as mice, rats and monkeys, and in humans. In fact, approaches for scoring *Pig-a*/*PIG-A* mutant cells have been described for all these

species [1–3,13–16]. Although credit for the original flow cytometric methodology for measuring mutant cells goes to David Araten and Lucio Luzzatto, who described *PIG-A* assays for human granulocytes and RBCs [17], human and monkey assays are now much less well developed than the rat assay. Also, while comparable methodology to the rat assay is available for the mouse, at present there is considerably less experience conducting the *Pig-a* assay in mice.

The ability to evaluate *Pig-a* mutants in both humans and laboratory animals means that hypotheses about the responsiveness of humans to potential mutagen exposures can be tested in animal models. Also, the human *PIG-A* assay may have value in clinical settings, where *Pig-a* assays conducted in laboratory animals could provide a translational biomarker for endogenous mutation *in vivo*. For example, the *Pig-a* and *PIG-A* assays could be used for monitoring the long-term effects of genotoxic chemotherapy [15] and aid in estimating the likelihood of developing secondary malignancies. Additionally, the *PIG-A* assay could be used in epidemiology studies for monitoring the health status (mutation accumulation) in populations exposed to potentially adverse environments, including occupational or accidental exposures to hazardous chemicals.

2. Topics discussed by the Workgroup

2.1. Placement of the assay within genotoxic testing strategies

The *Pig-a* assay should be considered an appropriate *in vivo* follow-up to positive results in bacterial and *in vitro* mammalian cell gene mutation assays. However, unless bone marrow exposure to the parent compound or to its metabolite(s) can be demonstrated directly or indirectly *via* plasma or by a reduction in the percentage of RETs, caution must be exercised in interpreting negative results. Moreover, based on analysis of the testing performed to date (see Section 2.2.3), the assay is not limited by any requirement of the test agent for metabolic activation. Although *Pig-a* is an endogenous gene, it is located on the X chromosome, and there is only one functional X chromosome. Thus, the Workgroup recognizes that the assay may not be suitable as a follow-up to either *in vitro* clastogenicity or aneugenicity findings since such events could potentially lead to cell death rather than mutation [18]. However, there currently is limited experience on the ability of the assay to detect clastogens or aneugens and the Workgroup considered the available information insufficient to justify a data-driven recommendation on this issue. Results from the *Pig-a* assay also could be used to build weight of evidence for the *in vivo* mutagenicity of a test agent, as might occur in cancer mode-of-action evaluations. Finally, the Workgroup noted the potential of using the *Pig-a* assay, as a measure of gene mutation, to complement the micronucleus (MN) assay, which measures clastogenicity/aneugenicity. Both can be readily included in routine *in vivo* safety evaluations, especially when the assays are integrated into subchronic (28 or 90 day) treatment protocols.

2.2. State of validation

Multi-laboratory trials initiated by Litron Laboratories [19], the Japanese Research Group [20], and the HESI initiative [5] have contributed to establishing protocols for the assay, testing the inter-laboratory reproducibility of the assay, and in expanding the number of agents tested.

2.2.1. Species, strain, and cell type

The assay has been investigated extensively only in rats. While the strain of rat may influence the absorption, distribution, metabolism, and elimination (ADME) of a test material and thus

affect the assay outcome, this is not anticipated to influence the biology of the endpoint; however, strain variation, if any, deserves further investigation.

Mammalian species used in other toxicological/investigational studies are amenable for use, theoretically, and there are a number of publications on *Pig-a* (or *PIG-A*) assays using other species (e.g., human [15,17], mouse [2,13]). However, a recommendation for their routine use in safety assessment studies cannot be made at this time because standard protocols have not yet been evaluated extensively or published.

The assay has been extensively investigated using RBCs and RETs. Measuring the *Pig-a* mutant phenotype in RETs has the advantage of observing an induced response more quickly than in RBCs and has the theoretical advantage of mitigating any effects of immune lysis (which does not seem to be a major factor in rats [9]). For routine screening, however, it is recommended that RBCs also should be investigated for GPI-anchored protein deficiency because of the ability to score a much larger sample size of cells. Additional studies have used other cell types in the *Pig-a* assay (e.g., lymphocytes [3,4], bone marrow erythroid cells [10]); but while these cell types are potentially useful, a recommendation for their routine use in safety assessment studies cannot be made at this time.

2.2.2. Intra/inter-laboratory reproducibility

The Workgroup could not identify any studies specifically intended to test the intra-laboratory reproducibility of the assay beyond assaying technical replicates (e.g., [21]; also see results in Fig. 2 and Table 3). There is, however, good evidence for a high degree of inter-laboratory transferability and reproducibility based on the results with several potent and weak mutagens (and low doses of potent mutagens) that were tested in a systematic manner as part of the Litron trial [22–27]. This trial employed common protocols (refined for each stage of the trial), and common reagents, although the rat strains and flow cytometers used were at the discretion of the participating laboratory. Stage I consisted of information gathering, and Stage II tested the transferability of an anti-CD59 antibody-based method using data from Litron Laboratories as the comparator. Stage II was based on assays employing acutely administered doses of *N*-ethyl-*N*-nitrosourea (ENU), and involved 14 laboratories [22]. The results demonstrated a high level of agreement, and provided confidence in extending the studies to investigate the performance characteristics of the assay further. Stage III utilized a 28-day repeat-dose design, and experiments were performed with ENU [23], dimethylbenz[*a*]anthracene (DMBA; [24]), *N*-methyl-*N*-nitrosourea [25], benzo[*a*]pyrene (BaP; [26]), and 4-nitroquinoline-1-oxide (4-NQO; [27]) in at least 2 laboratories per agent. Results from these studies also showed a high degree of concordance among laboratories. In addition, these data provided the impetus for establishing a method to dramatically increase the number of cells evaluated per sample through the use of immunomagnetic separation technology [21]. Stage IV studies have been conducted with these updated methods and demonstrate improved statistical power [21,28,29], but only a fraction of the completed work has been published or available to the Workgroup in a form that could be used to assess intra- and inter-laboratory reproducibility.

Results from the Japanese trial, using a different staining protocol, have begun to appear, and demonstrate a similar degree of inter-laboratory reproducibility [20]. Therefore, the Workgroup concluded that the available data show that the assay is both robust and demonstrably reproducible within and across laboratories for test agents with a range of mutagenic potency. The primary variable influencing the assay outcome across laboratories appears to be sample labeling and processing for flow cytometry.

Regarding assessments of sample labeling and processing, the Workgroup acknowledged that performing independent studies with the same test agent is not always the most effective way of evaluating the impact of these factors on mutant frequency measurements. Rather, results from blood samples processed at one site on multiple occasions, and blood samples split among multiple sites, can often be more informative. For instance, Fig. 2 shows the results from 15 Sprague Dawley rats treated with either the vehicle (i.e., negative control) or 1 of 2 dose levels of thiothepa. Four weeks after cessation of treatment for 28-days, blood samples were collected and processed by Litron personnel (at Rochester, NY) for same-day determination of mutant RET and RBC frequencies. A second aliquot of each whole blood sample was maintained in a refrigerator for next-day labeling and analysis at Litron, while a third set of coded aliquots was shipped overnight to Groton, CT to collaborators at Pfizer for labeling and flow cytometric analysis. All “next day” blood samples were maintained as whole blood in EDTA-coated vials until processing occurred. For the overnight shipment, the samples were kept cold, but not frozen. Both next-day analyses produced *Pig-a* mutant frequencies that were similar to those of the assays conducted on fresh samples (Fig. 2). Another study, using a different labeling protocol, demonstrated that the shipment of refrigerated blood samples had little effect on the measurement of *Pig-a* mutant frequencies [20]. These results strongly support the Workgroup’s conclusion that, where there are sufficiently trained personnel, the sample labeling and flow cytometric analysis procedures are highly reproducible within and across laboratories.

2.2.3. Test agent coverage

A list of 41 agents tested in the rat *Pig-a* assay was compiled by the Workgroup (Table 1). Most of the agents identified by the Workgroup were genotoxic in one or more tests, including 26 that were Ames’ test positive. The testing that has been conducted identified most of the agents expected to be positive as positive in the assay. Besides direct-acting simple alkylating agents, positive responses were detected from a number of chemicals that require metabolic activation in order to manifest their genotoxicity, including 2-acetylaminofluorene, aflatoxin B1, aristolochic acids, BaP, cyclophosphamide, diethylnitrosamine (DEN), dibenzo[*a,l*]pyrene, DMBA, and urethane. It is apparent from these observations that sufficient reactive metabolites from these agents damage the DNA of erythroid precursor cells to produce a response in the assay, even when other tissues may be responsible for the metabolism of these compounds. None of the unanticipated negative responses in the assay were clearly associated with a requirement for metabolic activation. Of note, the magnitude of positive responses indicates that the dynamic range of the *Pig-a* assay is at least 2–3 orders of magnitude: while background frequencies were $\leq 5 \times 10^{-6}$ (see Section 2.5.1), induced frequencies for the most potent mutagen tested to date, ENU, were (depending on the dose) in excess of 1000×10^{-6} (e.g., [9]).

The Workgroup found a few examples of responses produced by genotoxic agents that illustrated both the strengths and the limitations of the assay. Urethane, a genotoxicant requiring metabolic activation not present in many *in vitro* test systems [30], was positive in the assay. DEN, which is considered to be a liver-specific genotoxicant that is often negative for assays of bone marrow genotoxicity (like the erythrocyte MN assay), initially tested negative in the assay using the ‘basic’ protocol [24]. Subsequent testing using an immunomagnetic separation protocol produced a positive response [31]. Etoposide, azidothymidine (AZT), and 5-fluorouracil tested negative in the assay but all have mutagenic mechanisms which would be anticipated to be difficult to detect with an X-linked reporter gene, and only the basic protocol and/or acute dosing have been used for testing to date. Finally, aflatoxin B1 was positive in the assay, but only following subchronic, and not acute dosing,

Table 1
Agents tested in the rat *Pig-a* assay.

Compound	CAS no.	Chemical class	Expected outcome ^a	Rat <i>Pig-a</i> result		References	Comments
				Acute exposure ^b	Subchronic exposure ^c		
2-Acetylaminofluorene	53-96-3	Aromatic amine	+	+	+	[33]; Novartis (unpublished)	Requires metabolic activation; 2 labs
Acetaminophen	103-90-2	Hydroxylaniline	–	–HT	–HT	Janssen (unpublished)	Potent liver toxicant
Acrylamide	79-06-1	α,β-unsaturated amide	–/Weak+	–	–/? HT	FDA-NCTR (unpublished)	TGR, <i>Hprt</i> lymphocyte positive (2 mo treatment)
Aflatoxin B1	1162-65-8	Mycotoxin	+	–HT	Weak+	Janssen (unpublished)	Requires metabolic activation by CYP3A4 (not highly expressed in rats compared to humans)
o-Anthranilic acid	118-92-3	Aromatic amine	–	–HT	–HT	[33]	Non-alerting structure
Aristolochic acids	313-67-7	Aromatic nitro	+	+	+	[32]; Novartis (unpublished)	2 labs; weak <i>in vivo</i> MN
Azathioprine	446-86-6	Aromatic nitro	–/Weak +	Weak+	Weak+	[33]	Immunotoxicant
Azidothymidine	30516-87-1	Nucleoside analog	–/Weak+	–Basic (7 days)		[60]	<i>In vivo</i> MN positive; <i>Tk</i> mutation positive (mouse)
Benzo[<i>a</i>]pyrene	50-32-8	PAH	+	+	+	[26,29,37,52]	Requires metabolic activation; 3 labs
2-Butoxyethanol	111-76-2	Ethylene glycol derivative	–	–HT	–HT	Pfizer (unpublished)	Hemolytic agent
Chlorambucil	305-03-3	Alkylator; nitrogen mustard	+	+	+	[33,53]; Novartis (unpublished)	2 labs; twice at one lab
4-Chloro-1,2-diaminobenzene	95-83-0	Aromatic amine	–(Male)	–HT (male)	–HT (male)	GSK (unpublished)	<i>In vivo</i> MN positive acute exposure only; Females not tested
Cisplatin	15663-27-1	Antineoplastic; crosslinker	+(Female)	+	+	[54,58]	2 labs
Cyclophosphamide	6055-19-2	Alkylator, nitrogen mustard	Weak+	+	+	[33,56,58]	Requires metabolic activation; strong <i>in vivo</i> MN positive; 3 labs
Dibenzo[<i>a,l</i>]pyrene	191-30-0	PAH	+	+		Roche (unpublished)	Requires metabolic activation
Diethylnitrosamine	55-18-5	Nitrosamine	–(Blood)		–Basic; weak + HT	[24,31]	Requires metabolic activation; weak <i>in vivo</i> MN
[1-(3-dimethylaminopropyl)-3-ethylcarbodiimide hydrochloride, (EDAC)]	25952-53-8	Carbodiimide	–		–HT	[61]	Ames and <i>in vitro</i> MN positive; Tested up to MTD, rapid degradation into non-mutagenic EDAU in acid conditions and in aqueous solutions (SD rats)
7,12-Dimethylbenz[<i>a</i>]anthracene	57-97-6	PAH	+	+	+	[1,20,24,37,52]; Janssen (unpublished)	Requires metabolic activation; 5 labs; twice at one lab (SD & Wistar rats)
1,2-Dimethylhydrazine hydrochloride	306-37-6	Hydrazine	–(Blood)	–Basic		Pfizer (unpublished)	
Ethylmethane sulfonate	62-50-0	Alkylator	+	+	+	[56,62]; FDA-NCTR (unpublished)	3 labs
<i>N</i> -Ethyl- <i>N</i> -nitrosourea	759-73-9	Alkylator	+	+	+	[1,3,9,20,22,23,37,51,56,57,62]; PBR Laboratories (unpublished)	Multiple labs (13+ for 3 day; 5 for 28-day)

Table 1 (Continued)

Compound	CAS no.	Chemical class	Expected outcome ^a	Rat <i>Pig-a</i> result		References	Comments
				Acute exposure ^b	Subchronic exposure ^c		
Etoposide	33419-42-0	Topoisomerase inhibitor	–	–Basic	–Basic	Pfizer (unpublished)	Bacterial specific mutagen
5-Fluorouracil	51-21-8	Pyrimidine analog	–/?	–HT		NCDSE (unpublished)	Toxicity to red blood cells observed; positive in concurrent <i>in vivo</i> MN test
Furan	110-00-9	Heterocycle	–		–Basic (56 days)	[63]	Requires metabolic activation by CYP2E1; concurrent spleen lymphocyte <i>Hprt</i> and <i>Pig-a</i> and <i>cII</i> liver mutation assays were negative
Hydroxyurea	127-07-1	Inhibits DNA replication	–	–HT	–HT	[33]	
Ionizing radiation	N/A		+	+(1 and 4 days)		[58]; Litron Laboratories (unpublished)	2 labs
Isopropylmethane sulfonate	926-06-7	Alkylator	+	+(1 and 3 days)	+	Pfizer (in preparation)	
Isopropyltoluene sulfonate	2307-69-9	Alkylator	+	+(1 day)		Pfizer (unpublished)	
Melamine	108-78-1	Triaminotriazine	–	–HT		NCDSE (in preparation)	
Melphalan	148-82-3	Alkylator; nitrogen mustard	+	+	+	[33,53]; Novartis (unpublished)	2 labs; twice at one lab
Methylmethane sulfonate	66-27-3	Alkylator	Weak+	+	+	[33]; Roche (unpublished)	2 labs
Methylphenidate	113-45-1	Phenylethylamine derivative (psycho-stimulant)	–		–Basic (21 days)	[55]	
4-Nitroquinoline-1-oxide	56-57-5	Aromatic <i>N</i> -oxide	+	+	+	[20,27,37,52,56,64]; AbbVie (unpublished)	4 labs; twice at one lab
<i>N</i> -Nitroso- <i>N</i> -methylurea	684-93-5	Alkylator	+	+	+	[25,37,52,65]	3 labs; twice at one lab (basic and HT)
1,3-Propane sultone	1120-71-4	Alkylator	+	+	+	[28,53]; Novartis (unpublished)	2 labs; twice at one lab
Pyrene	129-00-0	PAH	–	–HT	–HT	[29]	
Sodium chloride	7647-14-5	Salt	–	–HT	–HT	[33]	
Sulfisoxazole	127-69-5	Aryl sulfonamide	–	–HT	–HT	[33]	
Temozolomide	85662-93-1	Imidazotetrazine; (Alkylating agent)	+	+(5 days)		Roche (unpublished)	
Thiotepa	52-24-4	Aziridine	+	+	+	[33,53]; Novartis (unpublished); Litron/Pfizer (summarized within)	2 labs; twice at one lab and a shipment study with a third lab
Urethane	51-79-6	Carbamate	+		+	[66]	Requires metabolism that is not available <i>in vitro</i> ; TGR positive; weak effects

^a Based on expert judgment of Ames and *in vitro* mammalian gene mutation data, mechanism of genotoxicity, and tissue distribution *in vivo*.

^b Acute is 3 consecutive days dosing unless otherwise noted.

^c Subchronic is 28 consecutive days dosing unless otherwise noted.

HT, immunomagnetic separation; only noted when a compound is negative; basic, no immunomagnetic separation; only noted when a compound is negative; TGR, transgenic rodent; –, negative; +, positive; ?, equivocal; N/A, not applicable; ND, no data; PAH, polycyclic aromatic hydrocarbon; NCDSE, National Center for Drug Safety Evaluation and Research.




Fig. 2. *Pig-a* data from 15 individual male Sprague Dawley rats treated with vehicle (C) or 1 of 2 dose levels of thiotepa (T). All animals were treated daily for 28 consecutive days. 'C' animals were treated with the vehicle, T13–T18 animals were given 0.492 mg thiotepa/kg/day, while animals T19–T24 were given 10 mg thiotepa/kg/day for 17 days and 5 mg/kg/day for the remaining 11 days. Blood samples were collected 28 days after the last treatment and processed by Litron Laboratories' personnel for same-day determination of *Pig-a* mutant RET and RBC frequencies. A second aliquot of each whole blood sample was maintained in a refrigerator overnight prior to labeling and analysis at Litron (Rochester, NY), while a third set of aliquots was shipped to collaborators at Pfizer (Groton, CT) for labeling and flow cytometric analysis. All "next day" blood samples were maintained as whole blood in EDTA-coated vials until processing occurred. In the case of the overnight shipment, vials were kept cold, not frozen. Each bar represents the *Pig-a* mutant frequency for each individual rat's blood sample at one of the 3 conditions tested. Note that some bars for the control and low-dose-treated rats are hardly visible on the figure. The results, which show extreme rat-to-rat variability among the high-dose animals at this late blood collection time point, indicate that across a wide dynamic range, *Pig-a* mutant frequencies are highly reproducible within and between laboratories.

while acrylamide was at best inconclusive after subchronic testing; acrylamide was not tested using acute treatment.

Based on these observations, the Workgroup concluded that definitive negative calls should be based on data from the most sensitive protocols (*i.e.*, using immunomagnetic separation), and tests conducted to the maximum tolerated dose (MTD) or the limit dose (see Section 2.4.2). This recommendation extends to agents anticipated to be negative, including Ames' negative agents. Of the 8 agents tested that generally are considered non-genotoxicants, 5 have been tested to the above standards (*o*-anthranilic acid, 2-butoxyethanol, pyrene, sodium chloride, sulfisoxazole), while 3 have been tested incompletely (furan,

melamine, methylphenidate). All 8 have tested negative in the assay. Of these 8 chemicals, 2-butoxyethanol is of particular significance because it caused marked intravascular lysis of RBCs and induced a strong compensatory erythropoiesis, yet it did not affect *Pig-a* mutant frequencies.

Within the group of agents that have been tested, the Workgroup noted two reports of assays (using 3- and 28-day treatments with aristolochic acids and 28-day treatments with 4-NQO) that were negative for erythrocyte MN induction, but positive for *Pig-a* mutation [27,32]. There also were 3 agents identified by the Workgroup (AZT, hydroxyurea and 5-fluorouracil) that induced micronucleated RETs but did not increase *Pig-a* mutant

erythrocyte frequencies. The experiments with hydroxyurea employed both 3- and 28-day study designs with immunomagnetic separation (Table 1) [33]. Similarly, etoposide was negative in limited *Pig-a* testing, and previous work with this non-DNA-reactive compound indicates that it would have tested positive for MN induction (e.g. [34]). It was noted that these observations support the routine use of the *Pig-a* assay in conjunction with the *in vivo* MN assay, especially when subchronic treatments are conducted, as the information from these 2 assays could be complementary and/or confirmatory.

The Workgroup recommended expanding the number of tested agents in general and Ames' negative agents in particular (see Section 3.1.3). The Workgroup further recommended testing more Ames' negative agents that might interfere with an erythrocyte-based assay (e.g., those that might induce erythropoiesis and/or lyse RBCs, like 2-butoxyethanol [Table 1]).

2.3. Evaluating mutations responsible for *Pig-a* mutants

Guidance provided by the Organization for Economic Cooperation and Development (OECD; [35]) indicates that assay validation should include a demonstration that the assay measures what it is designed to measure. In this case, it is important that a *Pig-a* gene mutation assay measures mutation induction in the *Pig-a* gene rather than, for instance, gene silencing or the clonal expansion of pre-existing *Pig-a* mutations. Establishing the mutational basis of the response is a challenge, however, because the assay measures the mutant phenotype in enucleated erythrocytes. To date, there is no direct evidence that *Pig-a* mutation is responsible for producing *Pig-a* mutant erythrocytes. However, there is good evidence from nucleated bone marrow and spleen cells from rodents, and extensive experience from human PNH patients (see below), that supports the conclusion that GPI-anchored protein-deficient cells, in fact, have mutations in the *Pig-a* gene.

Studies indicate that proaerolysin (a bacterial toxin that potentiates its toxicity through binding to GPI) selects T-lymphocyte clones from ENU-treated rats that contain mutations in the *Pig-a* gene [4,36]. In addition, flow-sorted erythroid mouse bone marrow cells that are deficient in GPI-anchored markers also contain mutations in the *Pig-a* gene [10]. In addition, it has been demonstrated that rat CD59-negative erythrocytes are indeed GPI anchor-deficient (*i.e.*, resistant to aerolysin) [37]. Furthermore, Lemieux et al. [38] reported similar dose-response relationships for BaP-induced *Pig-a* mutant erythrocytes compared with bone marrow *lacZ* mutations using a transgenic mouse model. The Workgroup noted, however, that these data fall short of demonstrating that the individual mutants identified by flow cytometry (from either lymphocytes or erythrocytes) contain (or contained) mutations in the *Pig-a* gene. There is also an extensive literature indicating that PNH patients contain large clonal expansions of lymphocytes with *Pig-a* mutations and large numbers of phenotypically mutant erythrocytes [39–41].

The Workgroup concluded that the weight of evidence was consistent with a direct association between *Pig-a* mutation (rather than gene silencing or enzymatic inactivation) and the phenotype measured in the assay, and indicated that lack of absolute proof should not preclude use of the *Pig-a* assay for regulatory purposes. The Workgroup, however, encourages further work in this area to better characterize the mutations presumably measured by the assay (see Section 3.1.1). Until this proof is provided, the Workgroup advised that the endpoint measured by the assay should be referred to as 'GPI-anchored protein-deficient erythrocytes' or, more simply, '*Pig-a* mutant erythrocytes'.

2.4. Assay protocol

2.4.1. Test animals/sex

The assay has been investigated extensively and standardized only in the rat. Methods exist, however, for measuring *Pig-a* mutation in other species, in particular, the mouse. However, a recommendation for their routine use in safety assessment studies cannot be made at this time as standard protocols have not yet been published.

So far, nearly all the rat studies have been conducted in males, with only limited data in females. Consequently, no definitive statement can be made as to whether X-chromosome silencing in females (*viz.*, "Lyonization") has any influence on assay sensitivity. Nevertheless, evidence from the X-linked hypoxanthine-guanine phosphoribosyltransferase (*Hprt*) locus indicates the frequency of *Hprt* mutants in circulating T-lymphocytes is similar in male and female humans [42] and mice (e.g. [43,44]). Female human carriers of *HPRT* germ line mutations are healthy but their somatic cells are mosaics in terms of enzyme activity because of the random inactivation of the X-chromosome; in contrast *HPRT*-mutant males develop one of two clinical syndromes: Lesch-Nyhan syndrome or the Kelley–Seegmiller syndrome (reviewed by [45]). The incidence of Lesch–Nyhan syndrome (*HPRT* nullizygous) is very rare in females, with only a handful of women described [46]; the mechanisms responsible are complex and non-random X-inactivation has been observed in a single female patient [47]. There are no notable publications indicating whether Lyonization has any influence on the *HPRT* cell mutation assay. Nevertheless, the Workgroup suggested establishing the relative sensitivity of male and female animals to agents with known genotoxic potential (see Section 3.1.2). In addition, it was suggested that future research should explore the response of the assay to chemicals that affect the methylation status of the genome, as such changes may be more problematic for an assay using a reporter on the X-chromosome. Until these questions are resolved, it is advisable to use results from females with a certain degree of caution. However, where human exposure to chemicals may be sex-specific, as with some pharmaceutical agents, the Workgroup also recommends that the test should be performed with animals of the appropriate sex.

2.4.2. Number of doses, maximum doses, and age and number of animals

The Workgroup recommends following criteria similar to those given in OECD Test Guidelines (TGs) for conducting acute and/or repeat dose *in vivo* toxicology and genetic toxicology studies (e.g., TG407, TG408, and TG488 [48–50]). In general, young adult animals should be used (note that evidence indicates that treating animals *in utero* produces unstable responses [51]). A minimum of 3 dose levels of the test article should be employed, with the high dose being the MTD or, if toxicity does not limit the dose, 1000 mg/kg/day for subchronic (14 days or more) and 2000 mg/kg/day for acute treatments. Range-finding studies to select dose levels should be conducted in the same species, strain, and sex as the main test. The MTD is defined as the highest dose that will be tolerated without evidence of excessive toxicity such as death, pain, suffering, excessive weight loss, and/or distress. At this time, there is an inadequate understanding on how perturbations in hematopoiesis affect the performance of the *Pig-a* assay. It is important, therefore, that the results of the assay are interpreted in the context of potential assay confounders.

Group sizes of 6 are recommended but 5 analyzable animals per group may be acceptable if justified (e.g., based on power calculations, see Table 2). The use of an animal group for a positive control is not considered mandatory provided an appropriate standard that "mimics" mutants [52] is used each time flow cytometric analysis is performed. The standard serves to establish the flow cytometer