

62. Kami M, Yamamoto K, Kanda Y, Saito T, Chiba S, Iga T, Hirai H: Influence of blood sampling sites on measurement of blood cyclosporine concentration. *Transplantation* 69:1750-1752,2000.
63. Mitani K, Yamagata T, Iida C, Oda H, Maki K, Ichikawa M, Asai T, Honda H, Kurokawa M, Hirai H: Non-redundant roles of the elongation factor MEN in postimplantation development. *Biochem Biophys Res Commun* 279:563-567,2000.
64. Kanzawa F, Koizumi F, Koh Y, Nakamura T, Tatsumi Y, Fukumoto H, Saijo N, Yoshioka T, Nishio K: In vitro synergistic interactions between the cisplatin analog nedaplatin and the DNA topoisomerase I inhibitor irinotecan and the mechanism of this interaction. *Clin Cancer Res* 7:202-209,2001.
65. Fukumoto H, Tamura T, Kamiya Y, Usuda J, Suzuki T, Kanzawa F, Kuh H-J, Ohe Y, Saijo N, Nishio K: Activation-induced apoptosis of peripheral lymphocytes treated with 7-hydroxystaurosporine, UCN-01. *Invest New Drugs* 17:335-341,2000.
66. Nishio K, Nakamura T, Koh Y, Kanzawa F, Tamura T, Saijo N: Oncoprotein 18-overexpression increases the sensitivity against vindesine in the human lung cancer cells against vindesine. *Cancer*, in press, 2001
67. Ellmeier W, Jung S, Sunshine MJ, Hatam F, Xu Y, Baltimore D, Mano H, Littman DR: Severe B cell deficiency in mice lacking the tec kinase family members Tec and Btk. *J Exp Med* 192:1611-1624,2000.
68. van Dijk TB, van Den Akker E, Amelsvoort MP, Mano H, Lowenberg B, von Lindern M: Stem cell factor induces phosphatidylinositol 3'-kinase-dependent Lyn/Tec/Dok-1 complex formation in hematopoietic cells. *Blood* 96:3406-3413,2000.
69. Yoshida K, Yamashita Y, Miyazato A, Ohya K, Kitanaka A, Ikeda U, Shimada K, Yamanaka T, Ozawa K, Mano H: Mediation by the protein-tyrosine kinase Tec of signaling between the B cell antigen receptor and Dok-1. *J Biol Chem* 275:24945-24952,2000.
70. Matsuda KM, Madoiwa S, Hasumi Y, Kanazawa T, Saga Y, Kume A, Mano H, Ozawa K, Matsuda M: A novel strategy for the tumor angiogenesis-targeted gene therapy: generation of angiostatin from endogenous plasminogen by protease gene transfer. *Cancer Gene Ther* 7:589-596,2000.
71. Maeda Y, Ikeda U, Oya Ki, Shimpo M, Ueno S, Okada K, Saito T, Mano H, Ozawa K, Shimada K: Endogenously generated nitric oxide by nitric-oxide synthase gene transfer inhibits cellular proliferation. *J Pharmacol Exp Ther* 292:387-393,2000.
72. Oyama M, Yazaki T, Ohigashi T, Hoshi M, Horiguchi Y, Oya M, Asakura H, Nakashima J, Tachibana M, Uyemura K, Murai M: Application of conditionally replicating herpes vector for gene therapy treatment of urologic neoplasms. *Mol Urol* 4:83-87,2000.
73. Hoshi M, Harada A, Kawase T, Uyemura K, Yazaki T: Antitumor effects of defective herpes simplex virus-mediated transfer of the TIMP-2 gene in malignant glioma U87 in vitro: consequences for an anticancer gene therapy. *Cancer Gene Ther* 7:799-805,2000.
74. Oyama M, Ohigashi T, Hoshi M, Nakashima J, Tachibana M, Murai M, Uyemura K, Yazaki T: Intravesical and intravenous therapy of human bladder cancer by the herpes vector G207.

Hum Gene Ther 11:1683-1693,2000.

75. Oyama M, Ohigashi T, Hoshi M, Murai M, Uyemura K, Yazaki T: Oncolytic viral therapy for human prostate cancer with attenuated, replication-competent herpes simplex virus type 1 mutant G207. *Jpn J Cancer Res* 91:1-6,2000.
76. Fujimori EK, Takeuchi K, Yazaki T, Uyemura K, Nojyo Y, Tamamaki N: Expression of L1 and TAG-1 in the corticospinal, callosal, and hippocampal commissural neurons in the developing rat telencephalon as revealed by retrograde and in situ hybridization double labeling. *J Comp Neurol* 417:275-288,2000.
77. Oyama M, Ohigashi T, Yazaki T, Murai M: Experimental therapy for bladder cancers using conditionally-replicating herpes vector. *Hum Cell* 13,34,2000.
78. Oyama M, Ohigashi T, Murai M, Uyemura K, Yazaki T: Treatment of human renal cell carcinomas by G207, a replication-competent multimitated herpes simplex virus 1. *J Urol*, in press, 2001.
79. Yano L, Shimura M, Taniguchi M, Hayashi Y, Suzuki T, Hatake K, Takaku, Ishizaka Y: Improved gene transfer to neuroblastoma cells by a monoclonal antibody targeting RET, a receptor tyrosine kinase. *Hum Gene Ther* 11:995-1004,2000.
80. Shimura M, Osawa Y, Yuo A, Hatake K, Takaku F, Ishizaka Y: Oxidative stress as a necessary factor in room temperature induced apoptosis of HL-60 cells, *J Leukoc Biol* 68:87-96,2000.
81. Okui N, Sakuma R, Kobayashi N, Yoshikura H, Kitamura T, Chiba J, Kitamura Y: Packageable antiviral therapeutics against human immunodeficiency virus type 1: virion-targeted virus inactivation by incorporation of a single-chain antibody against viral integrase into progeny virions. *Hum Gene Ther* 11:537-546,2000.
82. Kobayashi S, Shirai T, Kiyokawa E, Mochizuki N, Matsuda M, Fukui Y: Membrane recruitment of DOCK180 by binding to PtdIns(3,4,5)P3. *Biochem J* 354(Pt1):73-78,2001.
83. Ohba Y, Mochizuki N, Matsuo K, Yamashita S, Nakaya M, Hashimoto Y, Hamaguchi M, Kurata T, Nagashima K, Matsuda M. Rap2 as a slowly responding molecular switch in the Rap1 signaling cascade. *Mol Cell Biol* 20:6074-6083,2000.
84. Yamashita S, Mochizuki N, Ohba Y, Tobiume M, Okada Y, Sawa H, Nagashima K, Matsuda M. CalDAG-GEFIII activation of Ras, R-ras, and Rap1. *J Biol Chem* 275:25488-25493,2000.
85. Mochizuki N, Ohba Y, Kobayashi S, Otsuka N, Graybiel AM, Tanaka S, Matsuda M. Crk activation of JNK via C3G and R-Ras. *J Biol Chem* 275:12667-12671,2000.
86. Ohba Y, Mochizuki N, Yamashita S, Chan AM, Schrader JW, Hattori S, Nagashima K, Matsuda M. Regulatory proteins of R-Ras, TC21/R-Ras2, and M-Ras/R-Ras3. *J Biol Chem* 275:20020-20026,2000.