

Hige S, Sakamoto M, Kage M, Mizokami M, Narimatsu H.	assessment in patients with viral hepatitis.				
Miyagi Y, Nomura H, Yamashita N, Tanimoto H, <u>Ito K</u> , Masaki N, Mizokami M, Shibuya T.	Estimation of two real-time RT-PCR assays for quantitation of hepatitis C virus RNA during PEG-IFN plus ribavirin therapy by HCV genotypes and IL28B genotype.	J Infect Chemother	19(1)	63-69	2013
<u>Ito K</u> , Yotsuyanagi H, Yatsuhashi H, Krino Y, Takikawa Y, Saito t, AraseY, Imazeki F, Kurosaki M, Umemura T, Ichida T, Toyoda H, <u>Yoneda M</u> , Mita E, Yamamoto K, Michitaka K, Meshiro T, Tanuma J, Tanaka Y, Sugiyama M, Murata K, Masaki N, Mizokami M; Japanese AHB Study Group.	Risk factors for long-term persistence of serum hepatitis B surface antigen following acute hepatitis B virus infection in Japanese adults.	Hepatology	59(1)	89-97	2014
Murata K, Sugiyama M, Kimura T, Yoshio S, Kanto T, Kirikae I, Saito H, Aoki Y, Hiramine S, Matsui T, <u>Ito K</u> , Korenaga M, Imamura M, Masaki N, Mizokami M.	Ex vivo induction of IFN- λ 3 by a TLR7 agonist determines response to Peg-IFN/Ribavirin therapy in chronic hepatitis C patients.	J Gastroenterol	49(1)	126-37	2014
<u>Yoneda M</u>	Is the measurement of serum cytokeratin-18 useful for diagnosis or screening nonalcoholic steatohepatitis?	Hepatology Res (in press).			
Nakade Y, <u>Yoneda M</u>	Prognostic factors for Regression from impaired glucose tolerance to Normal Glucose Regulation in	Intern Med (in press)			

	Japanese patients with nonalcoholic fatty liver disease.				
Hayashi K, Katano Y, Masuda H, Ishizu Y, Kuzuya T, Honda T, <u>Ishigami M</u> , Itoh A, Hirooka Y, Nakano I, Ishikawa T, Urano F, Yoshioka K, Toyoda H, Kumada T, Goto H.	Pegylated interferon monotherapy in patients with chronic hepatitis C with low viremia and its relationship to mutations in the NS5A region and the single nucleotide polymorphism of interleukin-28B.	Hepatol Res	43	580-588	2013
田中弘教, <u>飯島尋子</u>	造影超音波検査－厚生労働省有井班の診断アルゴリズムを含めて－.	臨牀と研究	90	188-93	2013
<u>飯島尋子</u>	超音波造影剤による肝腫瘍診断	Mebio	30	39-45	2013
田中弘教, 池田直人, 岩田一也, 高嶋智之, 楊和典, 石井昭生, 坂井良行, 會澤信弘, 岩田恵典, 榎本平之, 齋藤正紀, 今西宏安, 西口修平, <u>飯島尋子</u>	穿刺針ナビゲーションシステム VirtuTRAX™ 安全使用のための短縮クリップ法.	肝臓	54	291-3	2013
田中弘教, <u>飯島尋子</u>	肝細胞がんの Vascular Invasion-画像診断から脈管浸潤の限界にせまる-.	肝胆膵	66	935-40	2013
Mitsunori Y, Tanaka S, Nakamura N, Ban D, Irie T, Noguchi N, Kudo A, <u>Iijima H</u> , Arii S.	Contrast-enhanced intraoperative ultrasound for hepatocellular carcinoma: high sensitivity of diagnosis and therapeutic impact.	J Hepatobiliary Pancreat Sci.	20	234-42	2013

Bota S, Sporea I, Peck-Radosavljevic M, Sirli R, Tanaka H, Iijima H, Saito H, Ebinuma H, Lupsor M, Badea R, Fierbinteanu-Braticevici C, Petrisor A, Friedrich-Rust M, Sarrazin C, Takahashi H, Ono N, Piscaglia F, Marinelli S, D'Onofrio M, Gallotti A, Salzl P, Popescu A, Danila M.	The influence of aminotransferase levels on liver stiffness assessed by Acoustic Radiation Force Impulse Elastography: A retrospective multicentre study.	Dig Liver Dis.			Epub ahead of print	2013
Tamura Y, Suda T, Aarii S, Sata M, Moriyasu F, Imamura H, Kawasaki S, Izumi N, Takayama T, Kokudo N, Yamamoto M, Iijima H, Aoyagi Y.	Value of Highly Sensitive Fucosylated Fraction of Alpha-Fetoprotein for Prediction of Hepatocellular Carcinoma Recurrence After Curative Treatment.	Dig Dis Sci.	58	2406-12		2013
Enomoto H, Sakai Y, Aizawa N, Iwata Y, Tanaka H, Ikeda N, Kunihiro H, You K, Ishii A, Takashima T, Iwata K, Saito M, Imanishi H, Iijima H, Nishiguchi S.	Association of amino acid imbalance with the severity of liver fibrosis and esophageal varices.	Ann Hepatol.	12	471-8		2013
Tanaka H, Iijima H, Higashiura A, Yoh K, Ishii A, Takashima T, Sakai Y, Aizawa N, Iwata K, Ikeda N, Iwata Y, Enomoto H, Saito M, Imanishi H, Hirota S, Fujimoto J, Nishiguchi S.	New malignant grading system for hepatocellular carcinoma using the Sonazoid contrast agent for ultrasonography.	J Gastroenterol.			Epub ahead of print	2013

飯島尋子, 井倉技, 中山晴夫, 小林正宏, 熊田博光, 井廻道夫.	血清アルブミン濃度が軽度～中等度に低下した肝硬変患者の QOL に及ぼすリーバクト R 配合顆粒の影響.	Medicine and Drug Journal.	49	127-39	2013
Singh S, Eaton JE, Murad MH, Tanaka H, <u>Iijima H</u> , Talwalkar JA.	Accuracy of Spleen Stiffness Measurement in Detection of Esophageal Varices in Patients With Chronic Liver Disease: Systematic Review and Meta-analysis.	Clin Gastroenterology Hepatol.		Epub ahead of print	2013
Aizawa N, Enomoto H, Takashima T, Sakai Y, Iwata K, Ikeda N, Tanaka H, Iwata Y, Saito M, Imanishi H, <u>Iijima H</u> , Nishiguchi S.	Thrombocytopenia in pegylated interferon and ribavirin combination therapy for chronic hepatitis C.	J Gastroenterology.		Epub ahead of print	2013
Inoue T, Hyodo T, Murakami T, Takayama Y, Nishie A, Higaki A, Korenaga K, Sakamoto A, Osaki Y, Aikata H, Chayama K, Suda T, Takano T, Miyoshi K, Koda M, Numata K, Tanaka H, <u>Iijima H</u> , Ochi H, Hirooka M, Imai Y, Kudo M.	Hypovascular Hepatic Nodules Showing Hypointense on the Hepatobiliary-Phase Image of Gd-EOB-DTPA-Enhanced MRI to Develop a Hypervascular Hepatocellular Carcinoma: A Nationwide Retrospective Study on Their Natural Course and Risk Factors.	Dig Dis.	31	472-9	2013
Sakata K, Hara M, Terada T, Watanabe N, Takaya D, Yaguchi S, Matsumoto T, Matsuura T, Shirouzu M, Yokoyama S, Yamaguchi T, Miyazawa K, <u>Aizaki H</u> , Suzuki T, Wakita T, Imoto M, Kojima S.	CV NS3 protease enhances liver fibrosis via binding to and activating TGF- β type I receptor.	Sci Rep.	22	3243	2013

Nakajima S, Watashi K, Kamisuki S, Tsukuda S, Takemoto K, Matsuda M, Suzuki R, <u>Aizaki H</u> , Sugawara F, Wakita T.	Specific inhibition of hepatitis C virus entry into host hepatocytes by fungi-derived sulochrin and its derivatives.	Biochem Biophys Res Commun.	440	515-20	2013
Watashi K, Liang G, Iwamoto M, Marusawa H, Uchida N, Daito T, Kitamura K, Muramatsu M, Ohashi H, Kiyohara T, Suzuki R, Li J, Tong S, Tanaka Y, Murata K, <u>Aizaki H</u> , Wakita T.	Interleukin-1 and Tumor Necrosis Factor- α Trigger Restriction of Hepatitis B Virus Infection via a Cytidine Deaminase Activation-induced Cytidine Deaminase (AID).	J Biol Chem.	288	31715-27	2013
Suzuki R, Matsuda M, Watashi K, <u>Aizaki H</u> , Matsuura Y, Wakita T, Suzuki T.	Signal peptidase complex subunit 1 participates in the assembly of hepatitis C virus through an interaction with E2 and NS2.	PLoS Pathog.	9	e1003589	2013
Matsumoto Y, Matsuura T, Aoyagi H, Matsuda M, Hmwe SS, Date T, Watanabe N, Watashi K, Suzuki R, Ichinose S, Wake K, Suzuki T, Miyamura T, Wakita T, <u>Aizaki H</u> .	Antiviral activity of glycyrrhizin against hepatitis C virus in vitro.	PLoS One.	18	e68992	2013
Akazawa D, Moriyama M, Yokokawa H, Omi N, Watanabe N, Date T, Morikawa K, <u>Aizaki H</u> , Ishii K, Kato T, Mochizuki H, Nakamura N, Wakita T.	Neutralizing antibodies induced by cell culture-derived hepatitis C virus protect against infection in mice.	Gastroenterology.	145	447-55	2013
Suzuki R, Ishikawa T, Konishi E, Matsuda M, Watashi K, <u>Aizaki H</u> , Takasaki T, Wakita T.	Production of single-round infectious chimeric flaviviruses with DNA-based Japanese encephalitis virus replicon.	J Gen Virol.	95	60-65	2014

Iwamoto M, Watashi K, Tsukuda S, Aly1 HH, Fukasawa M, Suzuki R, <u>Aizaki H</u> , Ito T, Koiwai O, Kusuhara H, Wakita T	Evaluation and Identification of hepatitis B virus entry inhibitors using HepG2 cells overexpressing a membrane transporter NTCP.	Biochem Biophys Res Commun.	443	808-13	2014
--	--	-----------------------------------	-----	--------	------

